



100 Majestic Way, Bangor, PA 18013 / [www.biospectra.us](http://www.biospectra.us)

ANALYTICAL METHOD FOR THE DETERMINATION OF  
TRACE METALS BY INDUCTIVELY COUPLED PLASMA  
OPTICAL EMISSION SPECTROSCOPY (ICP-OES) IN  
SODIUM HYDROXIDE

## TABLE OF CONTENTS

1.	PURPOSE:.....	3
2.	SCOPE:.....	3
3.	RESPONSIBILITIES: .....	3
4.	REFERENCES:.....	3
	TABLE 1: LIMITS FOR SODIUM HYDROXIDE PRODUCTS .....	4
5.	MATERIALS AND EQUIPMENT:.....	4
	TABLE 2: REFERENCE STANDARDS .....	5
6.	PROCEDURE:.....	6
	TABLE 3: INTERMEDIATE STANDARD .....	6
	TABLE 4: 15 PPB CALIBRATION STANDARD .....	7
	TABLE 5: 45 PPB CALIBRATION STANDARD .....	8
	TABLE 6: 60 PPB CALIBRATION STANDARD .....	9
	TABLE 7: CALIBRATION BLANK .....	10
7.	INSTRUMENT PROCEDURE:.....	11
	TABLE 8: EXAMPLE SAMPLE ANALYSIS SEQUENCE .....	11
	TABLE 9: ICP-OES PARAMETERS.....	12
	TABLE 10: LINEAR RANGE AND CORRESPONDING WAVELENGTH.....	12
8.	REPORTING: .....	12
	TABLE 11: RESULT REPORTING .....	13

## 1. PURPOSE:

- 1.1. To provide a procedure for the assessment of Trace Metals in Sodium Hydroxide products via ICP-OES. This procedure was assessed as a full quantitative option-1 procedure as per validation report BSI-RPT-2113 and follows the validation parameters for quantitation procedures as outlined in USP <233>.
- 1.2. Elements under USP <232> validated for this test method are as follows:
  - 1.2.1. Class 1: Cd and Pb
  - 1.2.2. Class 2A: Co and Ni
  - 1.2.3. Class 3: Ba, Cr, Cu, and Mo
  - 1.2.4. Class 4: Ca, Fe, Mg, Mn, and Zn
  - 1.2.5. Other: Bi and Sr

## 2. SCOPE:

- 2.1. Applies to Sodium Hydroxide along with related products manufactured at BioSpectra.
- 2.2. Applies to the Perkin Elmer Avio 500 S/N 081S1905062 ICP-OES, or qualified ICP-OES, located in the Quality Control (QC) Laboratory at the BioSpectra Bangor, PA facility.

## 3. RESPONSIBILITIES:

- 3.1. The Laboratory Technology Manager, or other designated individual, is responsible for the control, implementation, training, and maintenance of this method.
- 3.2. The Laboratory Services Staff are responsible for complying with the requirements of this procedure
- 3.3. If any abnormalities are determined during routine use of the ICP-OES or during calibration, Management shall be promptly notified. If necessary, the ICP-OES will be serviced and recalibrated by Perkin Elmer before being approved for use.

## 4. REFERENCES:

- 4.1. BSI-PRL-0852, Analytical Method Validation Protocol: Trace Metals in Sodium Hydroxide
- 4.2. BSI-RPT-2113, Analytical Method Validation Report: Trace Metals in Sodium Hydroxide
- 4.3. BSI-SOP-0362, Operation and Maintenance of the Perkin Elmer Avio 500 ICP-OES
- 4.4. ICH Guideline for Elemental Impurities Q3D Current
- 4.5. NexION Operation with Syngistix Software Guide
- 4.6. USP <232>, Elemental Impurities- Limits
- 4.7. USP <233>, Elemental Impurities- Procedures
- 4.8. USP <730> Plasma Spectrochemistry
- 4.9. USP <1730> Plasma Spectrochemistry—Theory and Practice

**TABLE 1: LIMITS FOR SODIUM HYDROXIDE PRODUCTS**

Element	ICH Class	30% LOQ (µg/g) in sample	50% Target (µg/g) in sample	100% Target (µg/g) in sample	150% Target (µg/g) in sample
Cd	1	0.90	1.5	3.0	4.5
Pb	1	0.90	1.5	3.0	4.5
Co	2A	0.90	1.5	3.0	4.5
Ni	2A	0.90	1.5	3.0	4.5
Ba	3	0.90	1.5	3.0	4.5
Cr	3	0.90	1.5	3.0	4.5
Cu	3	0.90	1.5	3.0	4.5
Mo	3	0.90	1.5	3.0	4.5
Ca	4	0.90	1.5	3.0	4.5
Fe	4	0.90	1.5	3.0	4.5
Mg	4	0.90	1.5	3.0	4.5
Mn	4	0.90	1.5	3.0	4.5
Zn	4	0.90	1.5	3.0	4.5
Bi	Not Applicable	0.90	1.5	3.0	4.5
Sr	Not Applicable	0.90	1.5	3.0	4.5

**5. MATERIALS AND EQUIPMENT:**

## 5.1. Equipment

- 5.1.1. Analytical Balance
- 5.1.2. Perkin Elmer Avio 500 ICP-OES S/N 081S1905062, or qualified ICP-OES
- 5.1.3. Micropipettes, Rainin or Eppendorf

## 5.2. Reagents

- 5.2.1. Nitric Acid, Trace metals grade or equivalent
- 5.2.2. Deionized (DI) water (Type 1 Ultrapure)

## 5.3. Consumable Supplies

- 5.3.1. SCP Digtubes<sup>®</sup> 15 mL, 50 mL, and 100 mL
- 5.3.2. Pipette Tips of various sizes

## 5.4. Personnel

- 5.4.1. All personnel that executed the protocol are trained on ICP-OES or are considered Subject Matter Experts. This test method will be assigned a mark as read training to QC analysts involved with the execution.

**TABLE 2: REFERENCE STANDARDS**

<b>Identification</b>	<b>Manufacturer</b>	<b>Concentrations / Elements</b>
1,000 µg/mL Cadmium Standard	SCP Science	Cd (1,000 µg/mL)
1,000 µg/mL Lead Standard	SCP Science	Pb (1,000 µg/mL)
1,000 µg/mL Cobalt Standard	SCP Science	Co (1,000 µg/mL)
1,000 µg/mL Nickel Standard	SCP Science	Ni (1,000 µg/mL)
1,000 µg/mL Barium Standard	SCP Science	Ba (1,000 µg/mL)
1,000 µg/mL Chromium Standard	SCP Science	Cr (1,000 µg/mL)
1,000 µg/mL Copper Standard	SCP Science	Cu (1,000 µg/mL)
1,000 µg/mL Molybdenum Standard	SCP Science	Mo (1,000 µg/mL)
1,000 µg/mL Calcium Standard	SCP Science	Ca (1,000 µg/mL)
1,000 µg/mL Iron Standard	SCP Science	Fe (1,000 µg/mL)
1,000 µg/mL Magnesium Standard	SCP Science	Mg (1,000 µg/mL)
1,000 µg/mL Manganese Standard	SCP Science	Mn (1,000 µg/mL)
1,000 µg/mL Zinc Standard	SCP Science	Zn (1,000 µg/mL)
1,000 µg/mL Strontium Standard	SCP Science	Sr (1,000 µg/mL)
1,000 µg/mL Bismuth Standard	SCP Science	Bi (1,000 µg/mL)
1,000 µg/mL Scandium Standard	SCP Science	Sc (1,000 µg/mL)
1,000 µg/mL Yttrium Standard	SCP Science	Y (1,000 µg/mL)

<sup>1</sup>Additional standards/custom standards can be used as long as the concentration remains the same in final preparations.

**6. PROCEDURE:**

- 6.1. All standards will be prepared volumetrically from stock solutions purchased from certified vendors. If the vendor supplied stock standard is within 2% of the nominal value as per the certificate of analysis, then the nominal value will be used to calculate the concentration of the standard. If the stock standard certificate of analysis value is greater than or less than 2% of the nominal value, then the certificate of analysis value will be used for the stock standard concentration.
- 6.2. Internal Standard Solution
- 6.2.1. Add 0.500 mL of Sc (1,000 µg/mL) and 0.500 mL of Y (1,000 µg/mL) to a 50 mL Digitube®.
- 6.2.2. Dilute to 50 mL final volume with deionized water.
- 6.2.3. Scale proportionally as needed for use.
- 6.3. Intermediate Standard Preparation
- 6.3.1. Prepare a standard solution containing the elements listed in Table 3, using the individual single source 1,000 µg/mL stock standards. Prepare by adding stock standards to a 15 mL Digitube®. Add deionized water to approximately 8 mL and add 1.0 mL of concentrated nitric acid. Dilute to final volume using DI Water.

**TABLE 3: INTERMEDIATE STANDARD**

Identification	Element	Stock Identification	Amount Added (mL)	Nitric Acid (mL)	Final Volume (mL)	Final Conc. (µg/mL)
Intermediate Standard	Cd	1,000 µg/mL Cd Std	0.150	1.0	10	15
	Pb	1,000 µg/mL Pb Std	0.150			15
	Co	1,000 µg/mL Co Std	0.150			15
	Ni	1,000 µg/mL Ni Std	0.150			15
	Ba	1,000 µg/mL Ba Std	0.150			15
	Cr	1,000 µg/mL Cr Std	0.150			15
	Cu	1,000 µg/mL Cu Std	0.150			15
	Mo	1,000 µg/mL Mo Std	0.150			15
	Ca	1,000 µg/mL Ca Std	0.150			15
	Fe	1,000 µg/mL Fe Std	0.150			15
	Mg	1,000 µg/mL Mg Std	0.150			15
	Mn	1,000 µg/mL Mn Std	0.150			15
	Zn	1,000 µg/mL Zn Std	0.150			15
	Bi	1,000 µg/mL Bi Std	0.150			15
	Sr	1,000 µg/mL Sr Std	0.150			15

## 6.4. 15 ppb Calibration Standard Preparation

- 6.4.1. Prepare a solution containing the elements listed in Table 4 below in 5.0% HNO<sub>3</sub> acid matrix.
- 6.4.2. Add 0.050 mL of intermediate standard to a separate 50 mL Digitube® followed by the addition of approximately 35 mL of deionized water.
- 6.4.3. Add 2.50 mL of nitric acid then dilute to 45 mL using deionized water.
- 6.4.4. Add 1.0 mL of internal standard solution and dilute to volume using deionized water.
- 6.4.5. Do not allow intermediate standard to contact concentrated acids while preparing solutions. (Standards are stable for 1 day).

**TABLE 4: 15 ppb CALIBRATION STANDARD**

Identification	Element	Intermediate Standard (mL)	Nitric Acid (mL)	Internal Standard Solution (mL)	Final Volume (mL)	Final Conc. (µg/L)
<b>15 ppb Calibration Standard</b>	Cd	0.050	2.50	1.0	50	15
	Pb					15
	Co					15
	Ni					15
	Ba					15
	Cr					15
	Cu					15
	Mo					15
	Ca					15
	Fe					15
	Mg					15
	Mn					15
	Zn					15
	Bi					15
Sr	15					

## 6.5. 45 ppb Calibration Standard Preparation

- 6.5.1. Prepare a solution containing the elements listed in Table 5 below in 5.0% HNO<sub>3</sub> acid matrix.
- 6.5.2. Add 0.150 mL of intermediate standard to a separate 50 mL Digitube<sup>®</sup> followed by the addition of approximately 35 mL of deionized water.
- 6.5.3. Add 2.50 mL of nitric acid then dilute to 45 mL using deionized water.
- 6.5.4. Add 1.0 mL of internal standard solution and dilute to volume using deionized water.
- 6.5.5. Do not allow intermediate standard to contact concentrated acids while preparing solutions. (Standards are stable for 1 day).

**TABLE 5: 45 ppb CALIBRATION STANDARD**

Identification	Element	Intermediate Standard (mL)	Nitric Acid (mL)	Internal Standard Solution (mL)	Final Volume (mL)	Final Conc. (µg/L)
<b>45 ppb Calibration Standard</b>	Cd	0.150	2.50	1.0	50	45
	Pb					45
	Co					45
	Ni					45
	Ba					45
	Cr					45
	Cu					45
	Mo					45
	Ca					45
	Fe					45
	Mg					45
	Mn					45
	Zn					45
	Bi					45
Sr	45					



## 6.6. 60 ppb Calibration Standard Preparation

- 6.6.1. Prepare a solution containing the elements listed in Table 6 below in 5.0% HNO<sub>3</sub> acid matrix.
- 6.6.2. Add 0.200 mL of intermediate standard to a separate 50 mL Digitube<sup>®</sup> followed by the addition of approximately 35 mL of deionized water.
- 6.6.3. Add 2.50 mL of nitric acid then dilute to 45 mL using deionized water.
- 6.6.4. Add 1.0 mL of internal standard solution and dilute to volume using deionized water.
- 6.6.5. Do not allow intermediate standard to contact concentrated acids while preparing solutions. (Standards are stable for 1 day).

**TABLE 6: 60 ppb CALIBRATION STANDARD**

Identification	Element	Intermediate Standard (mL)	Nitric Acid (mL)	Internal Standard Solution (mL)	Final Volume (mL)	Final Conc. (µg/L)
<b>60 ppb Calibration Standard</b>	Cd	0.200	2.50	1.0	50	60
	Pb					60
	Co					60
	Ni					60
	Ba					60
	Cr					60
	Cu					60
	Mo					60
	Ca					60
	Fe					60
	Mg					60
	Mn					60
	Zn					60
	Bi					60
Sr	60					

## 6.7. Calibration Blank

- 6.7.1. Prepare a solution containing 5.0% HNO<sub>3</sub> acid matrix as described in Table 7 below.
- 6.7.2. To a separate 50 mL Digitube<sup>®</sup>, add approximately 35 mL of DI Water.
- 6.7.3. Add 2.50 mL of nitric acid then dilute to 45 mL using DI Water.
- 6.7.4. Add 1.0 mL of internal standard solution and dilute to volume using DI Water.
- 6.7.5. Do not allow Internal Standard Solution to contact concentrated acids.

**TABLE 7: CALIBRATION BLANK**

Description	Nitric Acid (mL)	Internal Standard (mL)	Final Volume (mL)
Cal Blank	2.50	1.0	50

## 6.8. Method Blank Preparation

- 6.8.1. Add approximately 35 mL of deionized water to a 50 mL Digitube<sup>®</sup>.
- 6.8.2. Add 2.50 mL of nitric acid and swirl to mix.
- 6.8.3. Add deionized water to approximately 45 mL and then transfer 1.0 mL of Internal Standard.
- 6.8.4. Dilute to a final volume of 50 mL using deionized water and mix well.

## 6.9. Sample Preparation

- 6.9.1. Samples are stable for 1 day.
- 6.9.2. Weigh approximately 500 mg of sample into a 50 mL Digitube<sup>®</sup>.
- 6.9.3. Transfer approximately 20 mL of deionized water and swirl the solution to mix thoroughly.
- 6.9.4. Add 2.50 mL of nitric acid then swirl to mix.
- 6.9.5. Add deionized water to approximately 45 mL and then transfer 1.0 mL of Internal Standard.
- 6.9.6. Dilute to a final volume of 50 mL with deionized water and mix thoroughly.

**7. INSTRUMENT PROCEDURE:**

- 7.1. Perform the ICP-OES daily performance check prior to beginning the analytical sequence. Refer to Avio 500 ICP-OES SOP, BSI-SOP-0362, for Daily Check procedures.
- 7.2. A calibration curve of no less than two standards and a blank must be used. The calibration correlation coefficient (R) must be  $\geq 0.99$ .
- 7.3. Set up the sequence similar to that shown in Table 8.
- 7.4. Confirm the calibration by analyzing the 45 ppb standard after the calibration. The calibration check must recover  $\pm 20\%$  of the calculated theoretical concentration following calibration.
- 7.5. A re-analysis of the check standard will be performed a minimum of once every 10 samples and at the end of the analytical run.
- 7.6. Bracketing standard checks must recover NMT 20% of the calculated theoretical concentration for multi-element analysis. Additionally, the drift (calculated as absolute difference) between the bracketing standard checks are to be NMT 20% for each target element.
- 7.7. The sample concentration is calculated as:

$$\text{Conc. } (\mu\text{g/g}) = \frac{\text{Solution Conc. } (\mu\text{g/L}) \times \text{Solution vol. (L)} \times \text{Dilution Factor}}{\text{Sample Mass (g)}}$$

**TABLE 8: EXAMPLE SAMPLE ANALYSIS SEQUENCE**

ID	Type	Level
Cal Blank	Cal Blank	Level 1
15 ppb Cal Std	Cal Std	Level 2
45 ppb Cal Std	Cal Std	Level 3
60 ppb Cal Std	Cal Std	Level 4
Cal Blank Check	QC Check	Not Applicable
45 ppb Check Std 1	QC Check	Not Applicable
Method Blank	Sample	Not Applicable
Sample(s) 10 or less	Sample	Not Applicable
45 ppb Check Std 2	QC Check	Not Applicable

## 7.8. Instrument Setup and Parameters

- 7.8.1. Instrument settings are only listed as guidelines. Settings may be changed in order to accommodate changes in sample matrix or hardware configurations.
- 7.8.2. The gas flows for Plasma, Auxiliary, and Nebulizer can be set at 12 mL/min, 0.20 mL/min, and 0.70 mL/min, respectively.
- 7.8.3. The instrument method is stored under the Test Methods labelled as “NaOH\_TraceMetals” for trace metal analysis.
- 7.8.4. Instrument method can be truncated in order to selectively analyze metals as long as parameters match the full methods.
- 7.8.5. All wavelengths are tied to Sc 357.253 as the internal standard and performed in axial mode for detection. Other wavelengths that met validation parameters, as described in validation report BSI-RPT-2113, can be substituted at a later date if instrument conditions change over time.

**TABLE 9: ICP-OES PARAMETERS**

ICP-OES System	Perkin Elmer Avio 500 Inductively Coupled Plasma Optical Emission Spectrometry (ICP-OES)
Points per Peak	4
Replicates	3
Viewing Distance	15.0
Nebulizer Gas	Argon
Shear Gas	Compressed Air
Sample Rinses	Rinse-1: 30 sec at 1.0 mL/min 5.0% HNO <sub>3</sub> (or as applicable to mitigate carry over)

**TABLE 10: LINEAR RANGE AND CORRESPONDING WAVELENGTH**

Element	Wavelength (nm)	Linear Range (µg/L)	Element	Wavelength (nm)	Linear Range (µg/L)
Cd	214.440	9 - 60	Ca	396.847	9 - 60
Pb	220.353	9 - 60	Fe	259.939	9 - 60
Co	230.786	9 - 60	Mg	279.553	9 - 60
Ni	227.022	9 - 60	Mn	257.610	9 - 60
Ba	493.408	9 - 60	Zn	202.548	9 - 60
Cr	267.716	9 - 60	Bi	206.170	9 - 60
Cu	324.752	9 - 60	Sr	407.771	9 - 60
Mo	204.597	9 - 60			

## 8. REPORTING:

- 8.1. Any result below the 30% LOQ concentration will be reported as less than the corresponding LOQ value listed in Table 1. Results above the LOQ concentration will be reported in µg/g (ppm) according to Table 11 below.

**TABLE 11: RESULT REPORTING**

<b>Result</b>	<b>Reporting</b>
If < LOQ	Report as < LOQ
If $\geq$ LOQ and < 1.0 ppm	Report to two (2) decimal places
If $\geq$ LOQ and $\geq$ 1.0 ppm	Report to whole number