



100 Majestic Way, Bangor, PA 18013 / [www.biospectra.us](http://www.biospectra.us)

# GUANIDINE HYDROCHLORIDE TESTING METHODS

The information contained herein is the confidential property of BioSpectra. The recipient is responsible for its safe-keeping and the prevention of unauthorized appropriation, use, disclosure and copying.

Page **1** of **18**

Controlled Copy Number: 2, Controlled Copy Location: Website, Printed By: VIRGINIA.PENA, on 20 Jan 2026

## TABLE OF CONTENTS

1.	PURPOSE:.....	3
2.	SCOPE:.....	3
3.	RESPONSIBILITIES: .....	3
4.	REFERENCES: .....	3
5.	EQUIPMENT: .....	4
6.	REAGENTS: .....	4
7.	ANALYTICAL PROCEDURE: .....	6
	FIGURE 1: RAW SPECTRA .....	14
8.	COMPENDIAL DIFFENTIATIONS:.....	18
	TABLE 1: IN-HOUSE VALIDATED METHODS.....	18
	TABLE 2: PRODUCT QUALITY DESCRIPTION .....	18
	TABLE 3: APPROVED OUTSIDE TESTING .....	18

## **1. PURPOSE:**

- 1.1. To provide the Laboratory personnel with procedures for testing Guanidine Hydrochloride, Raw Material, In-Process, Finished Goods, and Stability samples.

## **2. SCOPE:**

- 2.1. Applies to the testing of Guanidine Hydrochloride Raw Material, In-Process, Finished Goods, and Stability samples in the Laboratory at all BioSpectra locations. Methods include testing for all types of Guanidine Hydrochloride sold by BioSpectra; only the specific tests required for the requested type must be tested.

## **3. RESPONSIBILITIES:**

- 3.1. The Laboratory Manager or qualified designee is responsible for training, maintenance and implementation of this procedure.
- 3.2. The Laboratory Technicians are responsible for compliance with the terms of this procedure. This includes notifying the appropriate personnel if any analyses fail to meet their respective specifications.
- 3.3. Laboratory Technicians are responsible for referencing the applicable summary sheet or batch record for analysis specifications.

## **4. REFERENCES:**

- 4.1. BSI-ATM-0071, Method of Analysis: Elemental Impurities by ICP-MS in Guanidine HCl 6M
- 4.2. BSI-ATM-0089, Analytical Method of Analysis: Trace Metals in Finished Goods Products by ICP-OES
- 4.3. BSI-ATM-0117, Guanidine Hydrochloride Related Substances via HPLC
- 4.4. BSI-ATM-0131, Analytical Method for the Determination of Trace Metals in BioTech Products
- 4.5. BSI-FRM-0716, Guanidine Hydrochloride Analytical Procedure
- 4.6. BSI-FRM-0728, Analytical Procedure for Gel Assays
- 4.7. BSI-FRM-0745, Analytical Procedure for Protease Assay
- 4.8. BSI-RPT-1766, Analytical Method Transfer Report: Guanidine Hydrochloride Related Substances via HPLC with UV Detection
- 4.9. BSI-SOP-0019, Result Reporting
- 4.10. BSI-SOP-0090, Lambda 25 UV/Vis Operation and Calibration
- 4.11. BSI-SOP-0094, Muffle Furnace SOP and Calibration
- 4.12. BSI-SOP-0095, DNase (Endonuclease) Assay
- 4.13. BSI-SOP-0096, RNase (Ribonuclease) Assay
- 4.14. BSI-SOP-0098, Balance SOP
- 4.15. BSI-SOP-0126, Laboratory Notebooks
- 4.16. BSI-SOP-0138, DNase (Exonuclease) Assay
- 4.17. BSI-SOP-0139, Protease Assay
- 4.18. BSI-SOP-0140, Standardization of Titrants
- 4.19. BSI-SOP-0143, Metrohm Titrando 907 Auto-Titrator SOP
- 4.20. BSI-SOP-0254, Spectrum Two UATR SOP
- 4.21. BSI-SOP-0255, XL200 pH/Conductivity Meter SOP

The information contained herein is the confidential property of BioSpectra. The recipient is responsible for its safe-keeping and the prevention of unauthorized appropriation, use, disclosure and copying.

- 4.22. BSI-SOP-0256, MP50 Melting Range Operation, Verification and Calibration SOP
- 4.23. BSI-SOP-0303, NexION 350X ICP-MS SOP
- 4.24. BSI-SOP-0362, Operation and Maintenance of the Perkin Elmer Avio 500 ICP-OES
- 4.25. BSI-SOP-0422, Empower 3 General Procedure
- 4.26. BSI-SOP-0573, MP90 Melting Range Operation, Verification, and Calibration SOP
- 4.27. BSI-SOP-0668, OPI-180 OD Handheld Colorimeter SOP
- 4.28. *Current ACS Reagent Chemicals*
- 4.29. *Current USP*

## 5. EQUIPMENT:

- 5.1. Lambda 25 UV/Vis Spectrophotometer, or equivalent
- 5.2. Analytical Balance
- 5.3. Perkin Elmer NexION 350X ICP-MS
- 5.4. Perkin Elmer Avio 500 ICP-OES
- 5.5. Perkin Elmer Spectrum Two UATR
- 5.6. Blue M Oven, or equivalent
- 5.7. MP50 Melting Point Apparatus
- 5.8. MP90 Melting Range Apparatus
- 5.9. OPI-180 OD Handheld Colorimeter
- 5.10. XL200 pH/Conductivity Meter or equivalent.
- 5.11. Muffle Furnace
- 5.12. Metrohm 907 Titrando Auto-Titrator.
- 5.13. Waters Alliance HPLC

## 6. REAGENTS:

- 6.1. **Acetate Buffer (pH 3.5):** Dissolve 62.5 g of ammonium acetate in 62.5 mL of purified water and add 47.0 mL of concentrated hydrochloric acid. Adjust, if necessary, with 6N ammonium hydroxide or 6N hydrochloric acid to a pH of 3.5, dilute with purified water to make 250 mL.
- 6.2. **Ammonium Hydroxide (6N):** Pipette 41.29 mL of concentrated Ammonium Hydroxide into a 100 mL volumetric flask. Dilute to volume with purified water.
- 6.3. **Barium Chloride TS:** Dissolve 30 g of barium chloride dihydrate in water to make 250 mL.
- 6.4. **Brucine Sulfate Solution:** Carefully dissolve 600 mg of brucine sulfate in 1000 mL of 70% sulfuric acid.
- 6.5. **Chloramine T Solution:** Dissolve 1.0 g of white, water-soluble Chloramine T in 100 mL of water. Prepare fresh for each use.
- 6.6. **Crystal Violet Indicator (1%):** Dissolve 1 g of crystal violet in 100 mL of Glacial Acetic Acid. Store in amber bottle.
- 6.7. **Cyanide Solution (stock):** Prepare from Potassium Cyanide or purchase commercially.
- 6.8. **Cyanide Stock Solution (5 ppm):** Dissolve 4 mg of potassium cyanide in purified water and dilute to 1000 mL with purified water.
- 6.9. **Ethyl Alcohol:** Purchased Commercially
- 6.10. **Eosin Y Indicator:** Dissolve 50 mg of Eosin Y in 10 mL of purified water.
- 6.11. **Ferric Chloride Solution (1 g/100 mL):** Dissolve 1.0 g of ferric chloride hexahydrate in purified water and dilute to 100 mL with purified water.
- 6.12. **Ferrous Sulfate Heptahydrate:** Purchased Commercially.

The information contained herein is the confidential property of BioSpectra. The recipient is responsible for its safe-keeping and the prevention of unauthorized appropriation, use, disclosure and copying.

- 6.13. **Glacial Acetic Acid:** Purchased commercially.
- 6.14. **Glycerin Base TS:** To 200 g of glycerol, add purified water to bring the total weight to 235 g. Add 140 mL of 1N Sodium Hydroxide and 50 mL of purified water.
- 6.15. **Guanidine Hydrochloride UATR Reference Standard:** Dry a purchased reference standard for 4 hours at 105°C. Compare to a previously approved reference standard. Correlation must achieve  $\geq 0.95$  to meet requirements.
- 6.16. **Hydrochloric Acid (HCl) (3N):** Pipette 25.75 mL of concentrated hydrochloric acid and transfer to a 100-mL volumetric flask that contains a small amount of purified water. Dilute to volume with purified water.
- 6.17. **Hydrochloric Acid HCl (0.1N):** Purchased Commercially.
- 6.18. **Lead Nitrate Stock Solution:** Dissolve 0.1598 g of lead nitrate in 100 mL purified water and add 1 mL of nitric acid. Dilute with purified water to 1000 mL. Store in glass container free from soluble lead salts.
- 6.19. **Mercuric Acetate Solution:** Weigh 50.0 g of Mercuric Acetate. Transfer to a 1000-mL volumetric flask and dissolve in Glacial Acetic Acid. Q.S to volume with Glacial Acetic Acid.
- 6.20. **Methanol:** Purchased commercially.
- 6.21. **Nitrate Standard (10 ppm):** Dilute 2.5 mL of Nitrate Stock Solution (1000 ppm NO<sub>3</sub>) to 250mL with purified water.
- 6.22. **Nitric Acid (concentrated):** Purchased Commercially.
- 6.23. **Perchloric Acid:** Purchased Commercially.
- 6.24. **Phenolphthalein Indicator:** Dissolve 1.0 g of phenolphthalein in 100 mL of reagent grade alcohol.
- 6.25. **Phosphate Buffer:** Dissolve 138 grams of sodium phosphate in water. Add 70 mL of Acetic acid.
- 6.26. **PVA (0.2%):** Dissolve 2.0 g of polyvinyl alcohol in approximately 800 mL of purified water while gently heating and stirring. Once dissolved, remove the stir bar and Q.S. to 1000 mL with purified water.
- 6.27. **Pyridine-Barbituric acid Reagent:** Place 60 g of barbituric acid in a 1000-mL volumetric flask and add just enough water to wash the sides of the flask and wet the barbituric acid. Add 300 mL of pyridine and mix. Add 60 mL of HCl, mix, and cool to room temperature. Dilute to 1000 mL with water and mix. This reagent is stable for approximately 6 months if stored in a cool, dark place. Filter when prepared or just before use, if insoluble form through Whatman No. 40 filter paper or equivalent.
- 6.28. **Silver Nitrate (0.1N):** Purchased commercially.
- 6.29. **Sodium Chloride (NaCl):** Prepare a crucible at 450°C for 30 minutes. Allow to cool in a desiccator and weigh a maximum of 10.0 g of an approved lot of sodium chloride. Dry at 450°C for 24 hours. Cool in a desiccator, transfer to a previously dried vial, and store in desiccator. Stable for 3 months.
- 6.30. **Sodium Hydroxide Solution (1N):** Dilute 160 mL of 50% NaOH to 3L with purified water or dissolve 40 g of NaOH pellets in purified water and dilute to 1L.
- 6.31. **Sodium Hydroxide Solution (0.2N):** Dilute 8 grams of NaOH pellets to 1L with purified water.
- 6.32. **Sodium Hydroxide (NaOH) (0.1N):** Purchased Commercially.
- 6.33. **Sulfuric Acid:** Purchased Commercially.

6.34. **Sulfuric Acid (0.020N H<sub>2</sub>SO<sub>4</sub>):** Slowly add 20 mL of 0.1N sulfuric acid to 80 mL of purified water to make a total volume of 100 mL.

6.35. **Thioacetamide TS:** Dissolve 4.0 g of Thioacetamide in 100 mL of purified water.

## 7. ANALYTICAL PROCEDURE:

### 7.1. IN-PROCESS ML ABSORBANCE:

- 7.1.1. Prepare 10 mL of sample by pipetting 5 mL water and 5 mL specified Mother Liquor into a 50-mL beaker. Swirl to dissolve completely.
- 7.1.2. Refer to Lambda 25 UV/Vis Operation and Calibration to determine the Absorbance of the sample.
- 7.1.3. Record results at specified wavelengths in the Guanidine Hydrochloride In-Process Testing Log Book or appropriate In-Process notebook and the appropriate laboratory documentation.
- 7.1.4. Notify the appropriate personnel if any results do not meet the specifications.

### 7.2. IN-PROCESS ML ASSAY:

- 7.2.1. Standardize 0.1N AgNO<sub>3</sub> as per Standardization of Titrants.
- 7.2.2. Accurately weigh 0.35 g of sample.
- 7.2.3. Transfer to a 250-mL beaker and dissolve with 10 mL of purified water.
- 7.2.4. Add 10 mL of glacial acetic acid, 100 mL of methanol, and 10 mL of a 0.2% polyvinyl alcohol solution.
- 7.2.5. Titrate with 0.1N AgNO<sub>3</sub> to a potentiometric end-point utilizing the Metrohm Titrando 907.

$$\% GHCL = \frac{mL_{AgNO_3} \times N_{AgNO_3} \times 9.553}{Sample\ Weight\ (g)} \qquad \% Cl = \frac{mL_{AgNO_3} \times N_{AgNO_3} \times 3.545}{Sample\ Weight\ (g)}$$

#### 7.2.6. Alternate Manual Assay Method:

- 7.2.6.1. Refer to Standardization of 0.1N Silver Nitrate (AgNO<sub>3</sub>) by hand as per Standardization of Titrants
- 7.2.6.2. Accurately weigh 0.35 g of sample.
- 7.2.6.3. Transfer to a 250-mL beaker and dissolve with 10 mL of purified water.
- 7.2.6.4. Add 10 mL of glacial acetic acid, 100 mL of methanol, and 0.5 mL of Eosin Y Indicator.
- 7.2.6.5. Titrate to a pink endpoint.

$$\% GHCL = \frac{mL_{AgNO_3} \times N_{AgNO_3} \times 9.553}{Sample\ Weight\ (g)}$$

### 7.3. IN-PROCESS WET CRYSTAL ABSORBANCE:

- 7.3.1. Prepare a 6M solution of the specified sample (Solution may be scaled as required).
- 7.3.2. Accurately weigh 14.3 g of sample. Transfer accurately weighed sample to a 50-mL graduated cylinder and q.s. to 25 mL with purified water. Swirl to dissolve completely.
- 7.3.3. Refer to Lambda 25 UV/Vis Spectrophotometer Operation and calibration to determine the absorbance of the sample.

- 7.3.4. Record results at specified wavelengths in the Guanidine Hydrochloride In-Process Testing Logbook or In-Process Laboratory notebook, the appropriate laboratory documentation, and in batch record.

#### 7.4. **IN-PROCESS DRY CRYSTAL ASSAY:**

- 7.4.1. Standardize 0.1N AgNO<sub>3</sub> as per Standardization of Titrants.
- 7.4.2. Accurately weigh 0.35 g of sample that has been previously dried at 105°C for four hours.
- 7.4.3. Transfer to a 250-mL beaker and dissolve with 10 mL of purified water.
- 7.4.4. Add 10 mL of glacial acetic acid, 100 mL of methanol, and 10 mL of a 0.2% polyvinyl alcohol solution.
- 7.4.5. Titrate with 0.1N AgNO<sub>3</sub> to a potentiometric end-point utilizing the Metrohm Titrand 907.
- 7.4.6. **Alternate Manual Assay Method:**
- 7.4.6.1. Refer to Standardization of 0.1N Silver Nitrate (AgNO<sub>3</sub>) by hand as per Standardization of Titrants
- 7.4.6.2. Accurately weigh 0.35 g of sample that has been previously dried at 105°C for four hours.
- 7.4.6.3. Transfer to a 250-mL beaker and dissolve with 10 mL of purified water.
- 7.4.6.4. Add 10 mL of glacial acetic acid, 100 mL of methanol, and 0.5 mL of Eosin Y Indicator.
- 7.4.6.5. Titrate to a pink endpoint.

$$\% GHCL = \frac{mL_{AgNO_3} \times N_{AgNO_3} \times 9.553}{Sample Weight (g)} \qquad \% Cl = \frac{mL_{AgNO_3} \times N_{AgNO_3} \times 3.545}{Sample Weight (g)}$$

#### 7.5. **ABSORBANCE (IDENTIFICATION B):**

- 7.5.1. Prepare a 6M solution of the specified sample (Solution may be scaled as needed).
- 7.5.2. Accurately weigh 57.2 g of sample into a 100-mL volumetric flask and Q.S. to the mark with purified water.
- 7.5.3. Mix thoroughly.
- 7.5.4. Refer to Lambda 25 UV/Vis Operation and Calibration to determine the absorbance of the sample.
- 7.5.5. NOTE: This 6M solution preparation may be used to perform the following tests in this procedure: Chloride Identification, Clarity and Color of Solution, Melamine, pH (6M), and Solubility (6M). Solution must be covered and used within 24 hours.

#### 7.6. **ACIDITY:**

- 7.6.1. Using an analytical balance, dissolve 20 g of sample in 200 mL purified water.
- 7.6.2. Add 0.15 mL of phenolphthalein indicator and titrate to a lasting pink endpoint with 0.01N NaOH.
- 7.6.3. Not more than 5.5 mL of titrant is required to obtain the pink endpoint (0.01% max).

#### 7.7. **AMMONIUM CHLORIDE:**

- 7.7.1. Accurately weigh and dissolve 0.8 g of sample in 100 mL purified water in an Erlenmeyer flask and mix well.
- 7.7.2. Add ~3 drops of Methyl Red Solution and mix well.
- 7.7.3. Neutralize the solution to a yellow color with dropwise additions of 0.1N NaOH VS.

The information contained herein is the confidential property of BioSpectra. The recipient is responsible for its safe-keeping and the prevention of unauthorized appropriation, use, disclosure and copying.

- 7.7.3.1. Note: If solution is red; neutralize with 0.1N HCl.
- 7.7.4. Add 15mL of neutralized formaldehyde solution.
  - 7.7.4.1. Solution Preparation Reference; Neutralized Formaldehyde Solution - Transfer 100 mL of formaldehyde solution into a suitable container. Add 2.0 g of magnesium carbonate and shake for several minutes. Filter and use the clear filtrate.
  - 7.7.4.2. If any red color due to the indicator (Methyl Red Solution) is present, titrate the solution with 0.1N Sodium Hydroxide to a yellow endpoint.
- 7.7.5. Add 10 drops of phenolphthalein TS and mix well.
- 7.7.6. Titrate with 0.1N NaOH VS until a light pink endpoint which lasts 10 minutes is visible. Record volume of titrant delivered.
  - 7.7.6.1. Note: A timer and reference solution of water against a white background is recommended to help visually detect whether the visible pink color is still present after 10 minutes.
- 7.7.7. Not more than 1.50 mL of titrant is required to obtain the light pink endpoint (1.0% max).

#### 7.8. **APPEARANCE AND COLOR:**

- 7.8.1. Place 25-50 g of the sample in a clean, dry glass beaker.
- 7.8.2. In an area with sufficient lighting, view the sample from all sides.
- 7.8.3. The sample should be white in color and characteristic of crystals.
- 7.8.4. If the appearance and color result is unable to be definitively determined visually, the sample may be analyzed using the Colorimeter. Refer to BSI-SOP-0668, OPI-180 OD Handheld Colorimeter SOP.
- 7.8.5. If the sample does not conform to these specifications, notify the appropriate personnel immediately.

#### 7.9. **APPEARANCE OF SOLUTION:**

- 7.9.1. Prepare a 6M solution of the specified sample (Solution may be scaled as needed).
  - 7.9.1.1. Accurately weigh 57.2 g of sample into a 100-mL volumetric flask and Q.S. to the mark with purified water.
  - 7.9.1.2. Mix thoroughly.
  - 7.9.1.3. Solution for “Absorbance” test may be used.
- 7.9.2. Clear (2.2.1.) Turbidimetry
  - 7.9.2.1. Rinse the sample bottle with the sample solution twice.
  - 7.9.2.2. Fill sample bottle with the sample Solution S to the white line.
  - 7.9.2.3. Coat outside of bottle with a thin coat of silicon oil.
  - 7.9.2.4. Remove any air bubbles from the solution by using a syringe.
  - 7.9.2.5. Allow the sample to sit capped for 2-3 minutes.
  - 7.9.2.6. Follow the appropriate SOP as follows:
    - 7.9.2.6.1. Stroudsburg- Measure and record the turbidity of the sample according to Portable Turbidimeter Operation and Calibration.
    - 7.9.2.6.2. Bangor- Measure and record the turbidity of the sample according to Bangor Portable Turbidimeter SOP.
  - 7.9.2.7. The sample solution must be < 3 NTU.
- 7.9.3. Colorless (2.2.2, Method II)

The information contained herein is the confidential property of BioSpectra. The recipient is responsible for its safe-keeping and the prevention of unauthorized appropriation, use, disclosure and copying.



- 7.9.3.1. Add 10 mL of Solution S into a Nessler Color Comparison Tube.
- 7.9.3.2. Add 10 mL of *USP Purified Water* into a second Nessler Color Comparison Tube.
- 7.9.3.3. Compare the colors in sufficient lighting, viewing vertically against a white background.
- 7.9.3.4. In order for the sample solution to be colorless, it must have the appearance of *USP Purified Water*.

#### 7.10. **ASSAY (DRIED BASIS) GHCl | Cl:**

---

- 7.10.1. Standardize 0.1N AgNO<sub>3</sub> as per Standardization of Titrants.
- 7.10.2. Accurately weigh 0.35 g of sample that has been previously dried at 105°C for four hours.
- 7.10.3. Transfer to a 250 mL beaker and dissolve with 10 mL of purified water.
- 7.10.4. Add 10 mL of glacial acetic acid, 100 mL of methanol, and 10 mL of a 0.2% polyvinyl alcohol solution.
- 7.10.5. Titrate with 0.1N AgNO<sub>3</sub> to a potentiometric end-point utilizing the Metrohm Titrand 907.

$$\% \text{ GHCL} = \frac{mL_{\text{AgNO}_3} \times N_{\text{AgNO}_3} \times 9.553}{\text{Sample Weight (g)}} \qquad \% \text{ Cl} = \frac{mL_{\text{AgNO}_3} \times N_{\text{AgNO}_3} \times 3.545}{\text{Sample Weight (g)}}$$

#### 7.10.6. **Alternate Manual Assay Method:**

---

- 7.10.6.1. Refer to Standardization of 0.1N Silver Nitrate (AgNO<sub>3</sub>) by hand as per Standardization of Titrants.
- 7.10.6.2. Accurately weigh 0.35 g of sample that has been previously dried at 105°C for four hours.
- 7.10.6.3. Transfer to a 250-mL beaker and dissolve with 10 mL of purified water.
- 7.10.6.4. Add 10 mL of glacial acetic acid, 100 mL of methanol, and 0.5 mL of Eosin Y Indicator.
- 7.10.6.5. Titrate to a pink endpoint.

$$\% \text{ GHCL} = \frac{mL_{\text{AgNO}_3} \times N_{\text{AgNO}_3} \times 9.553}{\text{Sample Weight (g)}}$$

#### 7.11. **ASSAY (GUANIDINE BASIS):**

---

- 7.11.1. Standardize 0.1N Perchloric by hand as per Standardization of Titrants.
- 7.11.2. NOTE: Perform a blank determination utilizing all reagents below minus the product, there will be a blank for standardization and a blank for the assay.
- 7.11.3. Accurately weigh 0.35 g of Guanidine hydrochloride into an appropriate beaker.
- 7.11.4. Add stir bar.
- 7.11.5. Add 140 mL of Glacial Acetic Acid.
- 7.11.6. Add 20 mL of 5% mercuric acetate in glacial acetic acid solution.
- 7.11.7. Add 0.1 mL of 100 mg/10 mL crystal violet in glacial acetic acid solution.
- 7.11.8. Dissolve sample while stirring.
- 7.11.9. Record Temperature of the Perchloric Acid.
- 7.11.10. Titrate to a blue green end point with 0.1N Perchloric Acid.

7.11.10.1. Note: Due to the high volumetric coefficient of expansion of glacial acetic acid, a temperature correction factor is included in the analysis as  $C_f$ .

$$C_f = [1 + 1.07 \times 10^{-3} (T_s - T_i)]$$

$T_s$  = temperature at standardization.  $T_i$  = temperature at sample titration

$$\% \text{ Guanidine} = \frac{(EP_{\text{Sample}} - EP_{\text{Blank}})(N \text{ of Perchloric Acid})(9.55318)(C_f)}{[\text{Sample (g)} - \text{Sample Weight (g)}](\% \text{ LOD})}$$

#### 7.12. **CHLORIDE IDENTIFICATION (NF IDENTIFICATION C):**

---

- 7.12.1. Prepare a 6M solution of the specified sample (Solution may be scaled as needed).  
Solution previously prepared for Absorbance (6M) may be utilized.
- 7.12.2. Accurately weigh 57.2 g of sample into a 100-mL volumetric flask and Q.S. to the mark with purified water.
- 7.12.3. Mix thoroughly.
- 7.12.4. Transfer 2 mL of the sample solution to a beaker and add 0.2 mL of 0.1N Silver Nitrate.
- 7.12.5. A white, curdy precipitate that is insoluble after the addition of 1 mL of Concentrated Nitric Acid is produced.
- 7.12.6. If no precipitate is produced, notify the appropriate personnel.
- 7.12.7. Add 4 mL of 6N Ammonium Hydroxide. The precipitate should dissolve after mild agitation.

#### 7.13. **CLARITY AND COLOR OF SOLUTION:**

---

- 7.13.1. Prepare a 6M solution of the specified sample. Solution previously prepared for Absorbance (6M) may be utilized.
- 7.13.2. Accurately weigh 57.2 g of sample into a 100-mL volumetric flask and Q.S. to the mark with purified water.
- 7.13.3. Mix thoroughly.
- 7.13.4. Sample solution should be clear and colorless and comparable to the color and turbidity of a sample of water visually.

#### 7.14. **CYANIDE:**

---

##### 7.14.1. Primary Cyanide Method:

- 7.14.1.1. Weigh 5 g of sample and transfer to a 50-mL beaker. Dissolve in 10 mL of purified water. Prepare ferrous sulfate solution immediately before use: Dissolve 0.45 g of ferrous sulfate in 50 mL of 0.1N Hydrochloric acid and dilute to 100 mL with purified water.
- 7.14.1.2. Prepare a standard solution by pipetting 0.1 mL of 5 ppm Cyanide Stock Solution and adding 9.9 mL of purified water. To the sample and standard, add 0.1 mL of 0.1N Sodium Hydroxide and 0.1 mL of the ferrous sulfate solution. Warm both the sample and standard then cool to room temperature. Add 0.2 mL of the ferric chloride solution (1g/100 mL). The sample solution must not develop more red or blue than that of the standard (0.5µg/10 mL).

##### 7.14.2. Alternative Cyanide Method:

- 7.14.2.1. Equipment and Reagents:
  - 7.14.2.1.1. Distillation apparatus: Use a suitable cyanide distillation apparatus.

The information contained herein is the confidential property of BioSpectra. The recipient is responsible for its safe-keeping and the prevention of unauthorized appropriation, use, disclosure and copying.

- 7.14.2.1.2. Pyridine- Barbituric Acid Reagent: Place 60 g of barbituric acid in a 1000-mL volumetric flask and add just enough water to wash the sides of the flask and wet the barbituric acid. Add 300 mL of pyridine and mix. Add 60 mL of HCl, mix, and cool to room temperature. Dilute to 1000 mL with water and mix. This reagent is stable for approximately 6 months if stored in a cool, dark place. Filter when prepared or just before use, if insoluble form through Whatman No. 40 filter paper or equivalent.
- 7.14.2.1.3. Sodium Hydroxide Solution, 1N- Dilute 160 mL of 50% NaOH to 3L with purified water or dissolve 40 g of NaOH pellets in purified water and dilute to 1L.
- 7.14.2.1.4. Sodium Hydroxide Solution 0.2N- Dilute 8 g of NaOH pellets to 1L with purified water.
- 7.14.2.1.5. Phosphate Buffer: Dissolve 138 g of Sodium Phosphate in water. Add 70 mL of Acetic Acid.
- 7.14.2.1.6. Chloramine T Solution: Dissolve 1.0 g of white, water-soluble Chloramine T in 100 mL of water. Prepare fresh each use.
- 7.14.2.1.7. Stock Cyanide Solution: Prepare from Potassium Cyanide or purchase commercially.
- 7.14.2.2. Distillation:
  - 7.14.2.2.1. Transfer 500 mL of the test sample [250 g of crystals plus 250 g of water mixed until dissolved] to a 1L distillation flask.
  - 7.14.2.2.2. Add several glass beads.
  - 7.14.2.2.3. Turn on the laboratory vacuum.
  - 7.14.2.2.4. Add 50 mL of 1N NaOH to the absorption bottle, and enough water to obtain an adequate depth of liquid to ensure complete scrubbing of the air stream.
  - 7.14.2.2.5. Connect it to the condenser outlet tube.
  - 7.14.2.2.6. Start cooling water through the condenser.
  - 7.14.2.2.7. Apply heat.
    - 7.14.2.2.7.1. CAUTION: Do not leave the apparatus unattended during this initial heating because the vacuum flow may need to be increased to prevent the solution from backing up into and overflowing from the air inlet tube or the absorber inlet.
- 7.14.2.3. Reflux for 2 hours. The rate of reflux should be such that the vapors rise to about 3/4-7/8 of the length of the condenser.
- 7.14.2.4. After 2 hours, turn off the heat and continue to vacuum flow for 15 minutes.
- 7.14.2.5. Disconnect the absorption bottle.
- 7.14.2.6. Transfer the contents of the absorption bottle (absorption liquid) to a 250-mL volumetric flask and dilute to the mark with water washings from the absorption bottle (at least 3 washings) and connecting tubes. The final concentration of this solution is 0.2N sodium hydroxide.

7.14.2.6.1. Note: If the following section cannot be done the same day as the distillation, transfer the solution from the volumetric flask to a clean plastic bottle for storage.

7.14.2.7. Determining CN- Concentration and Preparation:

7.14.2.7.1. Prepare a spiking solution of 1 ppm by: a) serial dilution of Potassium Cyanide diluted in 0.2N NaOH. b.) use a purchased standard solution diluted in 0.2N NaOH.

7.14.2.7.2. Transfer 50.0 mL of absorption liquid to 100-mL volumetric flask.

7.14.2.7.3. Prepare a blank of 50.0 mL of 0.2N NaOH in a 100-mL volumetric flask.

7.14.2.7.4. To 50.0 mL of 0.2N NaOH in a 100-mL volumetric flask, transfer 0.005 mg of Cyanide.

7.14.2.7.5. Add 15 mL of phosphate buffer solution to each flask and mix. (The pH of the solution now should be 5.5-6.5.)

7.14.2.7.6. Add 2.0 mL of Chloramine T solution to each flask and mix.

7.14.2.7.7. After 1 to 2 minutes, add 5.0 mL of pyridine-barbituric acid reagent to each flask and mix. Dilute to the mark with water and mix again.

7.14.2.7.8. Allow 8 minutes for color development.

7.14.2.7.9. Any red color in the sample should be less than that of the standard. (0.1 ppm maximum)

**7.15. ENDOTOXINS:**

---

7.15.1. In-House Method:

7.15.1.1. Weigh 25 mg +/-1 mg and transfer to a sterile tube.

7.15.1.2. Dissolve in ~5 mL of LAL reagent water.

7.15.1.3. Dilute to 10 mL with LAL reagent water and mix thoroughly.

7.15.1.4. Follow Endosafe nexgen-PTS Reader SOP.

7.15.2. Outside Testing:

7.15.2.1. Fill out applicable ARF in NAMSA Connect Portal.

7.15.2.2. Send 5 g of sample to NAMSA (Irvine) with SDS.

**7.16. ETHYL ALCOHOL SOLUBILITY:**

---

7.16.1. Dissolve 5 g of sample in 100 mL of reagent alcohol.

7.16.2. Sample solutions should dissolve and be a clear colorless solution to pass test.

**7.17. FREE ACID:**

---

7.17.1. Refer to Acidity.

**7.18. ENZYME ACTIVITY:**

---

7.18.1. RNase, DNase, and Protease per SOPs listed in Section 4.

**7.19. HEAVY METALS:**

---

7.19.1. Refer to section 7.34: Trace Metals.

7.19.2. Alternate Wet Method:

7.19.2.1. Standard Lead Solution: On the day of use, dilute 10.0 mL of Lead Nitrate Stock Solution with purified water to 100 mL.

- 7.19.2.2. Standard Preparation: Into a Nessler comparison tube pipette 2 mL of the Standard Lead Solution and QS to 40 mL in purified water.
- 7.19.2.3. Sample Preparation: Weigh 2 g of sample, transfer to a Nessler comparison tube, and dissolve in 40 mL of purified water.
- 7.19.2.4. Procedure: Adjust pH to between 3 and 4 using acetic acid. To all solutions, add 2 mL of pH 3.5 Acetate Buffer and 1.2 mL of thioacetamide-glycerin base TS (Mix 0.2 mL of thioacetamide TS and 1 mL of glycerin base TS, and heat in a boiling purified water bath for 20 seconds. Use immediately.). QS to 50 mL with purified water. Let stand for 2 minutes. Any brown color in the sample must not exceed that in the 0.01 mg standard (10 ppm).

## 7.20. **IDENTIFICATION (UATR):**

---

- 7.20.1. Follow Spectrum Two UATR SOP.

## 7.21. **LOSS ON DRYING @ 105°C:**

---


- 7.21.1. Dry a LOD vial in an oven at  $105 \pm 2^\circ\text{C}$  for at least 30 minutes. Cool for 15 minutes in a desiccator, weigh, and record results.
- 7.21.2. Place the vial on the analytical balance and tare the dried vial. Weigh 2.0 g of sample and record results.
- 7.21.3. Dry for 4 hours at  $105^\circ\text{C}$ . Cool for 15 minutes in desiccator.
- 7.21.4. Retain Sample for Assay, dried basis.
- 7.21.5. Reweigh and calculate the % LOD.

$$\% \text{ Loss on Drying} = \frac{\text{Initial Sample Weight (g)} - \text{Final Sample Weight (g)} \times 100}{\text{Initial Sample Weight (g)}}$$

## 7.22. **MELAMINE:**

---

- 7.22.1. Prepare a 6M solution of the specified sample (Solution may be scaled as needed). Solution previously prepared for Absorbance (6M) may be utilized.
- 7.22.2. Accurately weigh 57.2 g of sample into a 100-mL volumetric flask and Q.S. to the mark with purified water.
- 7.22.3. Mix thoroughly.
- 7.22.4. Refer to Lambda 25 UV/Vis Operation and Calibration to determine the absorbance of the sample in the range of 220-260 nm.
- 7.22.5. A shoulder in this part of the curve is acceptable but there should be no observable peak or plateau in this range to report as Passes Test.  
If any point exceeds an absorbance of 3 notify the Laboratory Manager immediately.  
Below is an example of a Melamine containing sample's spectrum.



**Figure 1: Raw Spectra**

**7.23. MELTING RANGE:**

- 7.23.1. Refer to BSI-SOP-0256, MP50 Melting Range Operation and Calibration SOP, or BSI-SOP-0573, MP90 Melting Range Operation, Verification, and Calibration SOP.

**7.24. LIMIT OF NITRATE:**

**Note: Prepare Brucine Sulfate at time of use.**

7.24.1. Sample Preparation (Solution A):

- 7.24.1.1. Add 0.40 g of sample to 2.0 mL of water, dilute to 50 mL with brucine sulfate reagent solution, and mix thoroughly.

7.24.2. Control Preparation (Solution B):

- 7.24.2.1. Add 0.40 g of sample to 2.0 mL of 10 ppm standard nitrate solution and dilute to 50 mL with brucine sulfate reagent solution, and mix thoroughly.

7.24.3. Blank Preparation (Solution C):

- 7.24.3.1. Transfer 50 mL of the Brucine Sulfate reagent solution to a 50-mL volumetric flask.

7.24.4. Procedure:

- 7.24.4.1. Heat the three solutions in a preheated (boiling) water bath for 10 min. Cool rapidly in an ice bath to room temperature. Utilizing the Lambda 25, dispense blank solution C into one of the two matched 10 mm cuvettes and place into the second cell holder (closest to the front of the instrument) in order to complete a "100% / 0A Baseline (Auto Zero)". **Note: The cell holder closest to the back of the instrument should remain empty throughout the entire analysis.** Remove blank solution C and scan sample solution A and control solution B respectively at 410 nm. Calculate the nitrate content as:

$$\% NO_3 = \frac{AU_{\text{Solution A}}}{AU_{\text{Solution B}} - AU_{\text{Solution A}}} \times \% \text{ maximum allowable}$$

**7.25. pH 5% SOLUTION @25°C:**

- 7.25.1. Dissolve 5.0 g of sample in 100 mL of purified water.

The information contained herein is the confidential property of BioSpectra. The recipient is responsible for its safe-keeping and the prevention of unauthorized appropriation, use, disclosure and copying.

- 7.25.2. Ensure the pH is stable before obtaining result.
- 7.25.3. Follow the appropriate SOP for calibration and pH measurement.

**7.26. pH 0.5M SOLUTION @25°C:**

---

- 7.26.1. Dissolve 4.8 g of sample in 100 mL of purified water.
- 7.26.2. Ensure the pH is stable before obtaining result.
- 7.26.3. Follow the appropriate SOP for calibration and pH measurement.

**7.27. pH (6M) @ 25°C:**

---

- 7.27.1. Prepare a 6M solution of the specified sample (Solution may be scaled as needed).  
Solution previously prepared for Absorbance (6M) may be utilized.
- 7.27.2. Accurately weigh 57.2 g of sample into a 100-mL volumetric flask and Q.S. to the mark with purified water.
- 7.27.3. Mix thoroughly.
- 7.27.4. Follow the appropriate SOP for calibration and pH measurement.
- 7.27.5. Set a calibrated timer for 15 minutes. Record the pH and temperature.

**7.28. RELATED SUBSTANCES:**

---

- 7.28.1. Refer to BSI-ATM-0117, Guanidine Hydrochloride Related Substances via HPLC, for HPLC method parameters and sample and standard preparation.

**7.29. RESIDUE ON IGNITION:**

---

- 7.29.1. Turn on muffle furnace and allow temperature to stabilize at 600°C. Follow muffle furnace SOP and calibration procedure for operation.
- 7.29.2. Inspect a quartz crucible for cracks, chips and discoloration.
- 7.29.3. Utilize the 10-inch forceps to insert and remove a crucible into the furnace.
- 7.29.4. Ignite the quartz crucible at  $600 \pm 50$  °C for 30 minutes. Cool in a desiccator for one hour and 30 minutes and weigh.
- 7.29.5. Weigh 5.0 g sample in the previously ignited quartz crucible. Moisten the sample with 2 mL of sulfuric acid.
- 7.29.6. Volatilize the sample with a suitable heating apparatus. Ensure that the sample does not boil over and sample is not lost.
- 7.29.7. The rate of heating should be such that from ½ to 1 hour is required to volatilize the sample.
- 7.29.8. Continue to heat the sample until all excess sulfuric acid has been volatilized.
- 7.29.9. Ignite in the muffle furnace at  $600 \pm 50$  °C for 30 minutes or until all carbon has been removed.
- 7.29.10. Cool in the desiccator for a minimum of an hour and a half and reweigh.

$$\% ROI = \frac{\text{Residue Weight (g)} \times 100}{\text{Sample Weight (g)}}$$

**7.30. SOLUBILITY (6M):**

---

- 7.30.1. Prepare a 6M solution of the specified sample (Solution may be scaled as needed).  
Solution previously prepared for Absorbance (6M) may be utilized.
- 7.30.2. Accurately weigh 57.2 g of sample into a 100-mL volumetric flask and Q.S. to the mark with purified water.
- 7.30.3. Mix thoroughly.

The information contained herein is the confidential property of BioSpectra. The recipient is responsible for its safe-keeping and the prevention of unauthorized appropriation, use, disclosure and copying.

- 7.30.4. Observe from all angles. Sample solution should be clear and colorless and comparable to the color and turbidity of a sample of water visually.

**7.31. SOLUBILITY IN WATER:**

---

- 7.31.1. Accurately weigh 60.0 g of sample and transfer to a 250-mL beaker.  
 7.31.2. Accurately measure 100 mL of purified water in a graduated cylinder and transfer to the sample beaker.  
 7.31.3. Swirl to dissolve the sample. If necessary, utilize a stirrer plate and Teflon encapsulated magnetic stirring bar.  
 7.31.4. The sample completely dissolves in purified water to report as Clear/Colorless when compared to a clear and colorless reference solution.

**7.32. CHLORIDE AND SULFATE, SULFATE:**

---

7.32.1. Sample Preparation:

- 7.32.1.1. Weigh 1.0 g of sample and transfer to a 50 mL Nessler Color Comparison Tube. Dissolve in 40 mL of purified water.

7.32.2. Standard Preparation:

- 7.32.2.1. Prepare a standard solution by pipetting 0.052 mL of 0.020 N H<sub>2</sub>SO<sub>4</sub> in a 50 mL Nessler Color Comparison Tube. QS to 40 mL with purified water.

7.32.3. Procedure:

- 7.32.3.1. To both the sample and standard solutions, add 1 mL of 3N HCl, 3 mL of Barium Chloride TS and QS to the line.  
 7.32.3.2. Cover with parafilm and mix by inversion.  
 7.32.3.3. Compare turbidity 10 minutes after addition of the barium chloride to the sample and standard solutions.  
 7.32.4. Any turbidity produced in the sample solution should not exceed that produced by the standard when viewed from above against a black surface.  
 7.32.5. If turbidity of the sample solution exceeds that of the standard notify the Laboratory Manager immediately.

**7.33. TOTAL CN:**

---

- 7.33.1. Outside Testing.

**7.34. TRACE METALS:**

---

- 7.34.1. Available Methods dependent on Product Code Requirements:

- 7.34.1.1. Primary Method of Analysis: Refer to Method of Analysis: Elemental Impurities by ICP-MS in Guanidine HCl 6M, BSI-ATM-0071 for sample preparation and analysis.  
 7.34.1.2. Refer to Analytical Method of Analysis: Trace Metals in Finished Goods Products by ICP-OES, BSI-ATM-0089 for sample preparation and analysis.  
 7.34.1.3. Refer to Analytical Method for Determination of Trace Metals in BioTech Products, BSI-ATM-0131.

**7.35. WATER INSOLUBLES:**

---

- 7.35.1. Accurately weigh 50.0 g of sample and transfer to a 150-mL beaker.  
 7.35.2. Add 25 mL of room temperature purified water. If necessary, utilize a Teflon encapsulated magnetic stirring bar and electric stir plate to dissolve sample.

The information contained herein is the confidential property of BioSpectra. The recipient is responsible for its safe-keeping and the prevention of unauthorized appropriation, use, disclosure and copying.



- 7.35.3. Dry a Gooch crucible and filter paper at  $105^{\circ} \pm 2^{\circ}\text{C}$  for 1 hour. Cool in ambient air for 15 minutes and weigh.
- 7.35.4. Filter sample solution through the Gooch crucible using a suitable vacuum pump.
- 7.35.5. Rinse sample vessel and crucible filter with at least 150 mL of purified water.
- 7.35.6. Dry the crucible at  $105^{\circ} \pm 2^{\circ}\text{C}$  for 1 hour. Cool in ambient air for 15 minutes and weigh.

$$\% \text{ Insolubles} = \frac{\text{Residue Weight (g)} \times 100}{\text{Sample Weight (g)}}$$

**7.36. WATER (BY KARL FISCHER TITRATION):**

---

- 7.36.1. Perform a standardization of the titrant (Composite 5) as per Standardization of Titrants.
- 7.36.2. Grind the sample in a dry mortar into a fine powder utilizing a pestle.
- 7.36.3. Immediately weigh 2.0 g of sample into the glass weighing spoon and tare it.
- 7.36.4. Transfer the sample to the KF vessel by removing the rubber septum and adding the sample into the titration vessel.
  - 7.36.4.1. Do not leave the rubber septum open for long periods of time as this will allow moisture to enter the titration vessel.
- 7.36.5. Return the weighing spoon to the balance, making sure not to lose any sample that was left behind. Once the weight stabilizes, press the print button on the balance.
- 7.36.6. Check to make sure there is no residual sample stuck to the sides of the titration vessel.
  - 7.36.6.1. If there is any sample stuck to the side, stop the stir bead from spinning before swirling the vessel to rinse the sides.
- 7.36.7. Once the method begins, check to ensure the sample is fully dissolved before the titration begins (i.e. before the stir command completes).
- 7.36.8. The moisture content will then be determined by the Metrohm Auto Titrando 907.

$$\% \text{ Moisture} = \frac{(\text{mL of Composite 5})\left(\frac{\text{mg}}{\text{mL}} \text{ of Composite 5}\right)(0.1)}{\text{Sample Weight (g)}}$$

**8. COMPENDIAL DIFFENTIATIONS:**

- 8.1. Compendial Analyses Not Harmonized: Not Applicable.
- 8.2. Compendial Harmonized Methods: Not Applicable.
- 8.3. In-House Validated Methods in accordance with USP General Chapter:

**Table 1: In-House Validated Methods**

<b>Analysis Name</b>
Trace Metals

- 8.4. In-House Methods for Product Quality Description:

**Table 2: Product Quality Description**

<b>Analysis Name</b>
Appearance and Color

- 8.5. Outside Approved Laboratory Testing (if required):

**Table 3: Approved Outside Testing**

<b>Analysis Name</b>
Endotoxins
Total CN
Microbial