

Product	Product Description	Product Status	Product size	Product UOM	Package size	Package UOM	Package type	Liner size	Liner type	Desiccant
BTRI-3250-01	Bis-Tris, 1 KG	Active	1.00 KG		1500.00	ML	White HDPE bottle, blue lid			
BTRI-3250-05	Bis-Tris, 5 KG	Active	5.00 KG		3.50	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
BTRI-3250-10	Bis-Tris, 10 KG	Active	10.00 KG		5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
BTRI-3250-25	Bis-Tris, 25 KG	Active	25.00 KG		12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
BTRI-3250-92	Bis-Tris, 100 G	Active	100.00 G		250.00	ML	White HDPE bottle, blue lid			
BTRI-3250-93	Bis-Tris, 500 G	Active	500.00 G		1000.00	ML	White HDPE bottle, blue lid			
BTRI-4250-01	Bis-Tris, Biotech, 1 KG	Active	1.00 KG		1500.00	ML	White HDPE bottle, blue lid			
BTRI-4250-05	Bis-Tris, Biotech, 5 KG	Active	5.00 KG		3.50	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
BTRI-4250-10	Bis-Tris, Biotech, 10 KG	In development	10.00 KG		5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
BTRI-4250-25	Bis-Tris, Biotech, 25 KG	Active	25.00 KG		12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
BTRI-4250-50	Bis-Tris, Biotech, 50 KG	Active	50.00 KG		30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
BTRI-4250-92	Bis-Tris, Biotech, 100 G	Active	100.00 G		250.00	ML	White HDPE bottle, blue lid			
BTRI-4250-93	Bis-Tris, Biotech, 500 G	Active	500.00 G		1000.00	ML	White HDPE bottle, blue lid			
CSMH-3250-01	Cysteamine Hydrochloride (2-MEA), 1 KG	In development	1.00 KG		1.25	GAL	White HDPE round pail, blue lid	(2x) 24 x 36	LDPE	4 Unit
CSMH-3250-05	Cysteamine Hydrochloride (2-MEA), 5 KG	Active	5.00 KG		3.50	GAL	White HDPE round pail, blue lid	(2x) 24 x 36	LDPE	4 Unit
CSMH-3250-10	Cysteamine Hydrochloride (2-MEA), 10 KG	Active	10.00 KG		5.00	GAL	White HDPE round pail, blue lid	(2x) 24 x 36	LDPE	4 Unit
CSMH-3250-25	Cysteamine Hydrochloride (2-MEA), 25 KG	Active	25.00 KG		12.00	GAL	White HDPE round pail, blue lid	(2x) 30 x 48	LDPE	
CSMH-3250-90	Cysteamine Hydrochloride (2-MEA), 100 G	Active	100.00 G		3.50	GAL	White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	4 Unit
CSMH-3250-91	Cysteamine Hydrochloride (2-MEA), 50 G	In development	50.00 G		0.60	GAL	White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	4 Unit
CSMH-3250-92	Cysteamine Hydrochloride (2-MEA), 100 G	Active	100.00 G		0.60	GAL	White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	4 Unit
CSMH-3250-94	Cysteamine Hydrochloride (2-MEA), 1 KG	In development	1.00 KG		0.60	GAL	White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	4 Unit
CSMH-3250-96	Cysteamine Hydrochloride (2-MEA), 1 KG	Active	1.00 KG		3.50	GAL	White HDPE round pail, blue lid	(2x) 24 x 36	LDPE	4 Unit
D010-3201-01	Dextran Powder, MW 10,000, 1 KG	Active	1.00 KG		0.60	GAL	White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	4 Unit
D010-3201-30	Dextran Powder, MW 10,000, 30 KG	Active	30.00 KG		30.00	GAL	Dark blue HDPE drum	(2x) 34 x 60	LDPE	4 Unit
DXS4-3250-05	Dextran Sulfate Sodium Salt, MW 40,000, 5 KG	In development	5.00 KG		3.50	GAL	White HDPE round pail, blue lid	(2x) 24 x 36	LDPE	4 Unit
DXSE-4250-01	Dextran Sulfate Sodium Salt, MW 8,000, 1 KG	Active	1.00 KG		0.60	GAL	White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	4 Unit
DXSE-4250-05	Dextran Sulfate Sodium Salt, MW 8,000, 5 KG	Active	5.00 KG		3.50	GAL	White HDPE round pail, blue lid	(2x) 24 x 36	LDPE	4 Unit
DXSE-4250-10	Dextran Sulfate Sodium Salt, MW 8,000, 10 KG	Active	10.00 KG		5.00	GAL	White HDPE round pail, blue lid	(2x) 24 x 36	LDPE	4 Unit
DXSE-4250-15	Dextran Sulfate Sodium Salt, MW 8,000, 15 KG	Active	15.00 KG		12.00	GAL	White HDPE round pail, blue lid	(2x) 30 x 48	LDPE	4 Unit
DXSE-4250-25	Dextran Sulfate Sodium Salt, MW 8,000, 25 KG	Active	25.00 KG		12.00	GAL	White HDPE round pail, blue lid	(2x) 30 x 48	LDPE	4 Unit
DXSE-4250-92	Dextran Sulfate Sodium Salt, MW 8,000, 100 G	Active	100.00 G		0.60	GAL	White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	4 Unit
DXSE-4250-93	Dextran Sulfate Sodium Salt, MW 8,000, 20 G	Active	20.00 G		0.60	GAL	White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	4 Unit
DXSE-4250-94	Dextran Sulfate Sodium Salt, MW 8,000, 500 G	Active	500.00 G		0.60	GAL	White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	4 Unit
DXSE-4252-01	Dextran Sulfate Sodium Salt, MW 8,000, 1 KG	In development	1.00 KG		0.60	GAL	White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	4 Unit
DXSE-4252-05	Dextran Sulfate Sodium Salt, MW 8,000, 5 KG	In development	5.00 KG		3.50	GAL	White HDPE round pail, blue lid	(2x) 24 x 36	LDPE	4 Unit
DXSE-4252-10	Dextran Sulfate Sodium Salt, MW 8,000, 10 KG	In development	10.00 KG		5.00	GAL	White HDPE round pail, blue lid	(2x) 24 x 36	LDPE	4 Unit
DXSE-4252-25	Dextran Sulfate Sodium Salt, MW 8,000, 25 KG	In development	25.00 KG		12.00	GAL	White HDPE round pail, blue lid	(2x) 30 x 48	LDPE	4 Unit
DXSE-4252-50	Dextran Sulfate Sodium Salt, MW 8,000, 50 KG	In development	50.00 KG		30.00	GAL	Dark blue HDPE drum	(2x) 34 x 60	LDPE	4 Unit
DXSE-4252-92	Dextran Sulfate Sodium Salt, MW 8,000, 100 G	In development	100.00 G		0.60	GAL	White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	4 Unit
DXSE-4252-93	Dextran Sulfate Sodium Salt, MW 8,000, 500 G	In development	500.00 G		0.60	GAL	White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	4 Unit
DXSE-5201-92	Dextran Sulfate Sodium Salt, MW 8,000, 100 G	Active	100.00 G		0.60	GAL	White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	4 Unit
GALP-3224-06	D-Galactose, Plant Derived, 1 KG	Active	1.00 KG		0.60	GAL	White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	
GALP-3224-25	D-Galactose, Plant Derived, 20 KG	In development	20.00 KG		12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
GALP-3250-05	D-Galactose, Plant Derived, 5 KG	In development	5.00 KG		3.50	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
GALP-3250-06	D-Galactose, Plant Derived, 1 KG	Active	1.00 KG		0.60	GAL	White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	

Product	Product Description	Product Status	Product size	Product UOM	Package size	Package UOM	Package type	Liner size	Liner type	Desiccant
GALP-3250-10	D-Galactose, Plant Derived, 10 kg	Active	10.00 KG		5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
GALP-3250-25	D-Galactose, Plant Derived, 20 kg	Active	20.00 KG		12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
GALP-3250-27	D-Galactose, Plant Derived, 25 kg	Active	25.00 KG		30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
GALP-3250-92	D-Galactose, Plant Derived, 100 G	Active	100.00 G		250.00	ML	White HDPE bottle, blue lid			
GALP-3251-05	D-Galactose, Plant Derived, 5 KG	Active	5.00 KG		3.50	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
GALP-3251-06	D-Galactose, Plant Derived, 1 KG	Active	1.00 KG		0.60	GAL	White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	
GALP-3251-10	D-Galactose, Plant Derived, 10 KG	Active	10.00 KG		5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
GALP-3251-21	D-Galactose, Plant Derived, 20 KG	Active	20.00 KG		12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
GALP-3251-27	D-Galactose, Plant Derived, 25 KG	Active	25.00 KG		30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
GALP-3251-50	D-Galactose, Plant Derived, 50 KG	Active	50.00 KG		30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
GALP-3252-06	D-Galactose, Plant Derived, 1 KG	Active	1.00 KG		0.60	GAL	White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	
GALP-3252-27	D-Galactose, Plant Derived, 25 KG	Active	25.00 KG		30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
GALP-3252-92	D-Galactose, Plant Derived, 100 G	Active	100.00 G		250.00	ML	White HDPE bottle, blue lid			
GALP-3350-06	D-Galactose, Plant Derived, 1 KG	Active	1.00 KG		0.60	GAL	White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	
GALP-3350-25	D-Galactose, Plant Derived, 20 KG	Active	20.00 KG		12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
GALP-3350-27	D-Galactose, Plant Derived, 25 KG	Active	25.00 KG		30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
GALP-3351-05	D-Galactose, Plant Derived, 5 KG	In development	5.00 KG		3.50	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
GALP-3351-06	D-Galactose, Plant Derived, 1 KG	Active	1.00 KG		0.60	GAL	White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	
GALP-3351-27	D-Galactose, Plant Derived, 25 KG	Active	25.00 KG		30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
GALP-3351-50	D-Galactose, Plant Derived, 50 KG	Active	50.00 KG		30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
GALP-3351-92	D-Galactose, Plant Derived, 100 G	In development	100.00 G		250.00	ML	White HDPE bottle, blue lid			
GALP-4250-05	D-Galactose, Plant Derived, Biotech, 5 KG	Active	5.00 KG		3.50	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
GALP-4250-06	D-Galactose, Plant Derived, Biotech, 6 KG	Active	1.00 KG		0.60	GAL	White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	
GALP-4250-10	D-Galactose, Plant Derived, Biotech, 10 KG	Active	10.00 KG		5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
GALP-4250-21	D-Galactose, Plant Derived, Biotech, 20 KG	Active	20.00 KG		12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
GALP-4250-27	D-Galactose, Plant Derived, Biotech, 25 KG	Active	25.00 KG		30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
GALP-4250-50	D-Galactose, Plant Derived, Biotech, 50 KG	Active	50.00 KG		30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
GALP-4250-92	D-Galactose, Plant Derived, Biotech, 100 G	Active	100.00 G		250.00	ML	White HDPE bottle, blue lid			
GHCL-3101-02	Guanidine Hydrochloride 6M Solution, 1 L	Active	1.00 L		1.00	L	Natural HDPE bottle			
GHCL-3101-72	Guanidine Hydrochloride 6M Solution, 200 L	Active	200.00 L		55.00	GAL	Black HDPE drum			
GHCL-3101-78	Guanidine Hydrochloride 6M Solution, 1135 L	Active	1135.00 L		330.00	GAL	Tote			
GHCL-3101-80	Guanidine Hydrochloride 6M Solution, 200 L	Active	200.00 L		55.00	GAL	Black HDPE drum			
GHCL-3201-59	Guanidine Hydrochloride, 50 KG	Active	50.00 KG		26.00	GAL	Fiber drum	(2x) 34 x 60	LDPE	4 Unit
GHCL-3202-01	Guanidine Hydrochloride, 1 KG	Active	1.00 KG		1500.00	ML	White HDPE bottle, blue lid			
GHCL-3202-25	Guanidine Hydrochloride, 25 KG	Active	25.00 KG		12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
GHCL-3202-50	Guanidine Hydrochloride, 50 KG	Active	50.00 KG		30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	4 Unit
GHCL-3202-93	Guanidine Hydrochloride, 500 G	Active	500.00 G		1000.00	ML	White HDPE bottle, blue lid			
GHCL-3220-01	Guanidine Hydrochloride, 1 KG	Active	1.00 KG		1500.00	ML	White HDPE bottle, blue lid			
GHCL-3220-10	Guanidine Hydrochloride, 10 KG	Active	10.00 KG		5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
GHCL-3220-25	Guanidine Hydrochloride, 25 KG	Active	25.00 KG		12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
GHCL-3220-50	Guanidine Hydrochloride, 50 KG	In development	50.00 KG		30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
GHCL-3220-51	Guanidine Hydrochloride, 50 KG	Active	50.00 KG		30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	4 Unit
GHCL-3220-92	Guanidine Hydrochloride, 100 G	Active	100.00 G		250.00	ML	White HDPE bottle, blue lid			
GHCL-3220-93	Guanidine Hydrochloride, 500 G	Active	500.00 G		1000.00	ML	White HDPE bottle, blue lid			
GHCL-3220-94	Guanidine Hydrochloride, 2.5 KG	Active	2.50 KG		3.50	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	

Product	Product Description	Product Status	Product size	Product UOM	Package size	Package UOM	Package type	Liner size	Liner type	Desiccant
GHCL-3224-25	Guanidine Hydrochloride, 25 KG	Active	25.00 KG		12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
GHCL-3224-93	Guanidine Hydrochloride, 500 G	Active	500.00 G		1000.00	ML	White HDPE bottle, blue lid			
GHCL-3225-01	Guanidine Hydrochloride, 1 KG	Active	1.00 KG		1500.00	ML	White HDPE bottle, blue lid			
GHCL-3225-25	Guanidine Hydrochloride, 25 KG	Active	25.00 KG		12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
GHCL-3225-92	Guanidine Hydrochloride, 100 G	Active	100.00 G		250.00	ML	White HDPE bottle, blue lid			
GHCL-3250-01	Guanidine Hydrochloride, 1 KG	Active	1.00 KG		1500.00	ML	White HDPE bottle, blue lid			
GHCL-3251-59	Guanidine Hydrochloride, 50 KG	Active	50.00 KG		26.00	GAL	Fiber drum	(2x) 34 x 60	LDPE	4 Unit
GHCL-3251-93	Guanidine Hydrochloride, 500 G	Active	500.00 G		1000.00	ML	White HDPE bottle, blue lid			
GHCL-3253-25	Guanidine Hydrochloride, 25 KG	Active	25.00 KG		12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
GHCL-3253-93	Guanidine Hydrochloride, 500 G	Active	500.00 G		1000.00	ML	White HDPE bottle, blue lid			
GHCL-3255-25	Guanidine Hydrochloride, 25 KG	In development	25.00 KG		12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
GHCL-4120-72	Guanidine Hydrochloride 6M Solution, 200L	Active	200.00 L		55.00	GAL	Black HDPE drum			
GHCL-4223-01	Guanidine Hydrochloride, 1 KG	Active	1.00 KG		1500.00	ML	White HDPE bottle, blue lid			
GHCL-4223-05	Guanidine Hydrochloride, 5 KG	Active	5.00 KG		3.50	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
GHCL-4223-10	Guanidine Hydrochloride, 10 KG	Active	10.00 KG		5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
GHCL-4223-25	Guanidine Hydrochloride, 25 KG	Active	25.00 KG		12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
GHCL-4223-50	Guanidine Hydrochloride, 50 KG	Active	50.00 KG		30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
GHCL-4223-92	Guanidine Hydrochloride, 100 G	Active	100.00 G		250.00	ML	White HDPE bottle, blue lid			
GHCL-4223-93	Guanidine Hydrochloride, 500 G	Active	500.00 G		1000.00	ML	White HDPE bottle, blue lid			
GTHI-3220-10	Guanidine Thiocyanate, 10 KG	Active	10.00 KG		5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
GTHI-3220-14-76	Guanidine Thiocyanate, 14.76 KG	Active	14.76 KG		30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE/ 25L EZI-f	
GTHI-3220-15	Guanidine Thiocyanate, 14.76 KG	In development	14.76 KG		12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
GTHI-3220-25	Guanidine Thiocyanate, 25 KG	Active	25.00 KG		12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
GTHI-3220-92	Guanidine Thiocyanate, 100 G	Active	100.00 G		250.00	ML				
GTHI-3220-95	Guanidine Thiocyanate, 10 G	In development	10.00 G		60.00	ML				
GTHI-4221-01	Guanidine Thiocyanate, Biotech, 1 KG	Active	1.00 KG		1500.00	ML	White HDPE bottle, blue lid			
GTHI-4221-05	Guanidine Thiocyanate, Biotech, 5 KG	Active	5.00 KG		3.50	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
GTHI-4221-10	Guanidine Thiocyanate, Biotech, 10 KG	Active	10.00 KG		5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
GTHI-4221-25	Guanidine Thiocyanate, Biotech, 25 KG	Active	25.00 KG		12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
GTHI-4221-50	Guanidine Thiocyanate, Biotech, 50 KG	Active	50.00 KG		30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
GTHI-4221-92	Guanidine Thiocyanate, Biotech, 100 G	Active	100.00 G		250.00	ML	White HDPE bottle, blue lid			
GTHI-4221-93	Guanidine Thiocyanate, Biotech, 500 G	Active	500.00 G		1000.00	ML	White HDPE bottle, blue lid			
GTHI-5203-10	Guanidine Thiocyanate, 10 KG	Active	10.00 KG		5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
GTHI-5220-10	Guanidine Thiocyanate, 10 KG	Active	10.00 KG		5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
HEPE-3220-01	HEPES, 1 KG	Active	1.00 KG		1500.00	ML	White HDPE bottle, blue lid			
HEPE-3220-06	HEPES, 1 KG	Active	1.00 KG		0.60	GAL	White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	
HEPE-3220-10	HEPES, 10 KG	Active	10.00 KG		5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
HEPE-3220-12	HEPES, 12 KG	Active	12.00 KG		5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
HEPE-3220-20	HEPES, 20 KG	Active	20.00 KG		12.00	GAL	White HDPE round pail, blue lid	(2x) 30 x 48	LDPE	
HEPE-3220-22	HEPES, 2 KG	Active	2.00 KG		1.25	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
HEPE-3220-25	HEPES, 25 KG	Active	25.00 KG		12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
HEPE-3220-50	HEPES, 50 KG	Active	50.00 KG		30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
HEPE-3220-91	HEPES, 250 G	Active	250.00 G		1000.00	ML	White HDPE bottle, blue lid			
HEPE-3220-92	HEPES, 100 G	Active	100.00 G		250.00	ML	White HDPE bottle, blue lid			
HEPE-3220-93	HEPES, 500 G	Active	500.00 G		1000.00	ML	White HDPE bottle, blue lid			

Product	Product Description	Product Status	Product size	Product UOM	Package size	Package UOM	Package type	Liner size	Liner type	Desiccant
HEPE-3221-10	HEPES, 10 KG	Active	10.00 KG		5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
HEPE-3250-05	HEPES, 5 KG	Active	5.00 KG		3.50	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
HEPE-3250-10	HEPES, 10 KG	Active	10.00 KG		5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
HEPE-3251-01	HEPES, 1kg	Active	1.00 KG		1500.00	ML	White HDPE bottle, blue lid			
HEPE-3251-25	HEPES, 25 KG	Active	25.00 KG		12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
HEPE-3251-50	HEPES, 50 KG	Active	50.00 KG		30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
HEPE-3251-92	HEPES, 100 G	Active	100.00 G		250.00	ML	White HDPE bottle, blue lid			
HEPE-3254-25	HEPES, 25 KG	In development	25.00 KG		12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
HEPE-3254-50	HEPES, 50 KG	In development	50.00 KG		30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
HEPE-3320-25	HEPES, 25 KG	In development	25.00 KG		12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
HEPE-3320-92	HEPES, 100 G	In development	100.00 G		250.00	ML	White HDPE bottle, blue lid			
HEPE-3351-25	HEPES, 25 KG	In development	25.00 KG		12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
HEPE-3351-92	HEPES, 100 G	In development	100.00 G		250.00	ML	White HDPE bottle, blue lid			
HEPE-3354-25	HEPES, 25 KG	In development	25.00 KG		12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
HEPE-3354-50	HEPES, 50 KG	Active	50.00 KG		30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
HEPE-4250-01	HEPES, Biotech, 1 KG	Active	1.00 KG		1500.00	ML	White HDPE bottle, blue lid			
HEPE-4250-05	HEPES, Biotech, 5 KG	Active	5.00 KG		3.50	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
HEPE-4250-10	HEPES, Biotech, 10 KG	Active	10.00 KG		5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
HEPE-4250-25	HEPES, Biotech, 25 KG	Active	25.00 KG		12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
HEPE-4250-50	HEPES, Biotech, 50 KG	Active	50.00 KG		30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
HEPE-4250-92	HEPES, Biotech, 100 G	Active	100.00 G		250.00	ML	White HDPE bottle, blue lid			
HEPE-4250-93	HEPES, Biotech, 500 G	Active	500.00 G		1000.00	ML	White HDPE bottle, blue lid			
HEPE-5201-26	HEPES, 25 KG	Active	25.00 KG		26.00	GAL	Fiber drum	(2x) 34 x 60	LDPE	
HEPE-5251-25	HEPES, 25 KG	Active	25.00 KG		12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
HEPE-5251-26	HEPES, 25 KG	In development	25.00 KG		26.00	GAL	Fiber drum	(2x) 34 x 60	LDPE	
IHCL-4101-02	6N HCl in IPA Solution, 500 ML	Active	500.00 ML		17.00	OZ	Black HDPE UN leakproof bottle			
IHCL-4101-97	6N HCl in IPA Solution, 180 KG	Active	180.00 KG		55.00	GAL	Black HDPE drum			
KBRO-2220-01	Potassium Bromide, 1 KG	Active	1.00 KG		1500.00	ML	White HDPE bottle, blue lid			
KBRO-2220-05	Potassium Bromide, 5 KG	Active	5.00 KG		3.50	GAL	White HDPE round pail, blue lid	(2x) 24 x 36	LDPE	4 Unit
KBRO-2220-10	Potassium Bromide, 10 KG	Active	10.00 KG		5.00	GAL	White HDPE round pail, blue lid	(2x) 24 x 36	LDPE	4 Unit
KBRO-2220-25	Potassium Bromide, 25 KG	Active	25.00 KG		12.00	GAL	White HDPE round pail, blue lid	(2x) 30 x 48	LDPE	4 Unit
KBRO-2220-92	Potassium Bromide, 100 G	Active	100.00 G		0.60	GAL	White HDPE round pail, blue lid	9 x 15	LDPE	
KBRO-2301-01	Potassium Bromide, 1 KG	Active	1.00 KG		1500.00	ML	White HDPE bottle, blue lid			
KBRO-2301-25	Potassium Bromide, 25 KG	Active	25.00 KG		12.00	GAL	White HDPE round pail, blue lid	(2x) 30 x 48	LDPE	
KBRO-2302-01	Potassium Bromide, 1 KG	Active	1.00 KG		1500.00	ML	White HDPE bottle, blue lid			
KBRO-2302-25	Potassium Bromide, 25 KG	Active	25.00 KG		12.00	GAL	White HDPE round pail, blue lid	(2x) 30 x 48	LDPE	
LARH-5350-92	L-Arginine Hydrochloride, 100 G	Active	100.00 G		250.00	ML	White HDPE bottle, blue lid			
LARH-6250-01	L-Arginine Hydrochloride, Biotech, 1 KG	In development	1.00 KG		1500.00	ML	White HDPE bottle, blue lid			
LARH-6250-05	L-Arginine Hydrochloride, Biotech, 5 KG	In development	5.00 KG		3.50	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
LARH-6250-10	L-Arginine Hydrochloride, Biotech, 10 KG	In development	10.00 KG		5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
LARH-6250-25	L-Arginine Hydrochloride, Biotech, 25 KG	In development	25.00 KG		12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
LARH-6250-50	L-Arginine Hydrochloride, Biotech, 50 KG	In development	50.00 KG		30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
LARH-6250-92	L-Arginine Hydrochloride, Biotech, 100 G	In development	100.00 G		250.00	ML	White HDPE bottle, blue lid			
LARH-6250-93	L-Arginine Hydrochloride, Biotech, 500 G	In development	500.00 G		1000.00	ML	White HDPE bottle, blue lid			
LARH-6350-01	L-Arginine Hydrochloride, Biotech, 1 KG	Active	1.00 KG		1500.00	ML	White HDPE bottle, blue lid			

Product	Product Description	Product Status	Product size	Product UOM	Package size	Package UOM	Package type	Liner size	Liner type	Desiccant
LARH-6350-05	L-Arginine Hydrochloride, Biotech, 5 KG	Active	5.00 KG		3.50	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
LARH-6350-10	L-Arginine Hydrochloride, Biotech, 10 KG	Active	10.00 KG		5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
LARH-6350-25	L-Arginine Hydrochloride, Biotech, 25 KG	Active	25.00 KG		12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
LARH-6350-50	L-Arginine Hydrochloride, Biotech, 50 KG	Active	50.00 KG		30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
LARH-6350-92	L-Arginine Hydrochloride, Biotech, 100 G	Active	100.00 G		250.00	ML	White HDPE bottle, blue lid			
LARH-6350-93	L-Arginine Hydrochloride, Biotech, 500 G	Active	500.00 G		1000.00	ML	White HDPE bottle, blue lid			
LCYS-6250-01	L-Cystine Dihydrochloride, Biotech, 1 KG	Active	1.00 KG		1500.00	ML	White HDPE bottle, blue lid			
LCYS-6250-05	L-Cystine Dihydrochloride, Biotech, 5 KG	Active	5.00 KG		3.50	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
LCYS-6250-10	L-Cystine Dihydrochloride, Biotech, 10 KG	Active	10.00 KG		5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
LCYS-6250-25	L-Cystine Dihydrochloride, Biotech, 25 KG	Active	25.00 KG		12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
LCYS-6250-50	L-Cystine Dihydrochloride, Biotech, 50 KG	Active	50.00 KG		30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
LCYS-6250-92	L-Cystine Dihydrochloride, Biotech, 100 G	Active	100.00 G		250.00	ML	White HDPE bottle, blue lid			
LCYS-6250-93	L-Cystine Dihydrochloride, Biotech, 500 G	Active	500.00 G		1000.00	ML	White HDPE bottle, blue lid			
LCYS-6350-01	L-Cystine Dihydrochloride, Biotech, 1 KG	Active	1.00 KG		1500.00	ML	White HDPE bottle, blue lid			
LCYS-6350-05	L-Cystine Dihydrochloride, Biotech, 5 KG	Active	5.00 KG		3.50	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
LCYS-6350-10	L-Cystine Dihydrochloride, Biotech, 10 KG	Active	10.00 KG		5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
LCYS-6350-25	L-Cystine Dihydrochloride, Biotech, 25 KG	Active	25.00 KG		12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
LCYS-6350-50	L-Cystine Dihydrochloride, Biotech, 50 KG	Active	50.00 KG		30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
LCYS-6350-92	L-Cystine Dihydrochloride, Biotech, 100 G	Active	100.00 G		250.00	ML	White HDPE bottle, blue lid			
LCYS-6350-93	L-Cystine Dihydrochloride, Biotech, 500 G	Active	500.00 G		1000.00	ML	White HDPE bottle, blue lid			
LGLM-6201-01	L-Glutamine, Biotech, 1 KG	Active	1.00 KG		1500.00	ML	White HDPE bottle, blue lid			
LGLM-6201-05	L-Glutamine, Biotech, 5 KG	Active	5.00 KG		3.50	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
LGLM-6201-10	L-Glutamine, Biotech, 10 KG	Active	10.00 KG		5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
LGLM-6201-25	L-Glutamine, Biotech, 25 KG	Active	25.00 KG		12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
LGLM-6201-50	L-Glutamine, Biotech, 50 KG	Active	50.00 KG		30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
LGLM-6201-92	L-Glutamine, Biotech, 100 G	Active	100.00 G		250.00	ML	White HDPE bottle, blue lid			
LGLM-6201-93	L-Glutamine, Biotech, 500 G	Active	500.00 G		1000.00	ML	White HDPE bottle, blue lid			
LGLM-6301-01	L-Glutamine, Biotech, 1 KG	In development	1.00 KG		1500.00	ML	White HDPE bottle, blue lid			
LGLM-6301-05	L-Glutamine, Biotech, 5 KG	In development	5.00 KG		3.50	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
LGLM-6301-10	L-Glutamine, Biotech, 10 KG	In development	10.00 KG		5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
LGLM-6301-25	L-Glutamine, Biotech, 25 KG	In development	25.00 KG		12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
LGLM-6301-50	L-Glutamine, Biotech, 50 KG	In development	50.00 KG		30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
LGLM-6301-92	L-Glutamine, Biotech, 100 G	In development	100.00 G		250.00	ML	White HDPE bottle, blue lid			
LGLM-6301-93	L-Glutamine, Biotech, 500 G	In development	500.00 G		1000.00	ML	White HDPE bottle, blue lid			
LHMM-6250-01	L-Histidine Monohydrochloride Monohydrate, 1 KG	In development	1.00 KG		1500.00	ML	White HDPE bottle, blue lid			
LHMM-6250-05	L-Histidine Monohydrochloride Monohydrate, 5 KG	In development	5.00 KG		3.50	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
LHMM-6250-10	L-Histidine Monohydrochloride Monohydrate, 10 KG	In development	10.00 KG		5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
LHMM-6250-25	L-Histidine Monohydrochloride Monohydrate, 25 KG	In development	25.00 KG		12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
LHMM-6250-50	L-Histidine Monohydrochloride Monohydrate, 50 KG	In development	50.00 KG		30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
LHMM-6250-92	L-Histidine Monohydrochloride Monohydrate, 100 G	In development	100.00 G		250.00	ML	White HDPE bottle, blue lid			
LHMM-6250-93	L-Histidine Monohydrochloride Monohydrate, 500 G	In development	500.00 G		1000.00	ML	White HDPE bottle, blue lid			
LHMM-6350-01	L-Histidine Monohydrochloride Monohydrate, 1 KG	In development	1.00 KG		1500.00	ML	White HDPE bottle, blue lid			
LHMM-6350-05	L-Histidine Monohydrochloride Monohydrate, 5 KG	In development	5.00 KG		3.50	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
LHMM-6350-10	L-Histidine Monohydrochloride Monohydrate, 10 KG	In development	10.00 KG		5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
LHMM-6350-25	L-Histidine Monohydrochloride Monohydrate, 25 KG	In development	25.00 KG		12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	

Product	Product Description	Product Status	Product size	Product UOM	Package size	Package UOM	Package type	Liner size	Liner type	Desiccant
LHMM-6350-50	L-Histidine Monohydrochloride Monohydrate, 50 KG	In development	50.00 KG		30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
LHMM-6350-92	L-Histidine Monohydrochloride Monohydrate, 100 G	In development	100.00 G		250.00	ML	White HDPE bottle, blue lid			
LHMM-6350-93	L-Histidine Monohydrochloride Monohydrate, 500 G	In development	500.00 G		1000.00	ML	White HDPE bottle, blue lid			
MESH-3201-50	MES Hydrate, 50 KG	In development	50.00 KG		30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
MESH-3201-92	MES Hydrate, 100 G	Active	100.00 G		250.00	ML	White HDPE bottle, blue lid			
MESH-3250-01	MES Hydrate, 1 KG	Active	1.00 KG		1500.00	ML	White HDPE bottle, blue lid			
MESH-3250-25	MES Hydrate, 25 KG	Active	25.00 KG		12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
MESH-3350-01	MES Hydrate, 1 KG	In development	1.00 KG		1500.00	ML	White HDPE bottle, blue lid			
MESH-4250-01	MES Hydrate, Biotech, 1 KG	Active	1.00 KG		1500.00	ML	White HDPE bottle, blue lid			
MESH-4250-05	MES Hydrate, Biotech, 5 KG	Active	5.00 KG		3.50	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
MESH-4250-10	MES Hydrate, Biotech, 10 KG	Active	10.00 KG		5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
MESH-4250-25	MES Hydrate, Biotech, 25 KG	Active	25.00 KG		12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
MESH-4250-50	MES Hydrate, Biotech, 50 KG	Active	50.00 KG		30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
MESH-4250-92	MES Hydrate, Biotech, 100 G	Active	100.00 G		250.00	ML	White HDPE bottle, blue lid			
MESH-4250-93	MES Hydrate, Biotech, 500 G	Active	500.00 G		1000.00	ML	White HDPE bottle, blue lid			
MESH-9250-91	MES Hydrate, 50 G	In development	50.00 G		250.00	ML	White HDPE bottle, blue lid			
MESM-3220-01	MES Monohydrate, 1 KG	Active	1.00 KG		1500.00	ML	White HDPE bottle, blue lid			
MESM-3220-05	MES Monohydrate, 5 KG	Active	5.00 KG		3.50	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
MESM-3220-10	MES Monohydrate, 10 KG	Active	10.00 KG		5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
MESM-3220-12	MES Monohydrate, 12 KG	Active	12.00 KG		5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
MESM-3220-25	ME, Monohydrate, 25 KG	Active	25.00 KG		12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
MESM-3220-50	MES Monohydrate, 50 KG	Active	50.00 KG		30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
MESM-3220-91	MES Monohydrate, 50G	Active	50.00 G		250.00	ML	White HDPE bottle, blue lid			
MESM-3220-92	MES Monohydrate, 100G	Active	100.00 G		250.00	ML	White HDPE bottle, blue lid			
MESM-3220-93	MES Monohydrate, 500G	Active	500.00 G		1000.00	ML	White HDPE bottle, blue lid			
MESM-3222-01	MES Monohydrate, 1 KG	Active	1.00 KG		1500.00	ML	White HDPE bottle, blue lid			
MESM-3222-03	MES Monohydrate, 3 KG	Active	3.00 KG		3.50	GAL	White HDPE round pail, blue lid	24 x 36		
MESM-3222-25	MES Monohydrate, 25 KG	Active	25.00 KG		12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
MESM-3222-50	MES Monohydrate, 50 KG	Active	50.00 KG		30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
MESM-3222-91	MES Monohydrate, 50G	Active	50.00 G		250.00	ML	White HDPE bottle, blue lid			
MESM-3222-92	MES Monohydrate, 100G	Active	100.00 G		250.00	ML	White HDPE bottle, blue lid			
MESM-3222-93	MES Monohydrate, 500G	In development	500.00 G		1000.00	ML	White HDPE bottle, blue lid			
MESM-3250-06	MES Monohydrate, 1 KG	Active	1.00 KG		0.60	GAL	White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	
MESM-3250-10	MES Monohydrate, 10 KG	Active	10.00 KG		5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
MESM-3250-12	MES Monohydrate, 12 KG	Active	12.00 KG		5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
MESM-3250-25	MES Monohydrate, 25 KG	Active	25.00 KG		12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
MESM-3250-50	MES Monohydrate, 50 KG	Active	50.00 KG		30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
MESM-3250-91	MES Monohydrate, 250G	Active	250.00 G		1000.00	ML	White HDPE bottle, blue lid			
MESM-3250-92	MES Monohydrate, 100G	Active	100.00 G		250.00	ML	White HDPE bottle, blue lid			
MESM-3250-93	MES Monohydrate, 500G	Active	500.00 G		1000.00	ML	White HDPE bottle, blue lid			
MESM-3250-95	MES Monohydrate, 10G	In development	10.00 G		250.00	ML	White HDPE bottle, blue lid			
MESM-3251-25	MES Monohydrate, 25 KG	Active	25.00 KG		12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
MESM-3251-92	MES Monohydrate, 100G	Active	100.00 G		250.00	ML	White HDPE bottle, blue lid			
MESM-3320-25	MES Monohydrate, 25 KG	In development	25.00 KG		12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
MESM-3320-92	MES Monohydrate, 100G	In development	100.00 G		250.00	ML	White HDPE bottle, blue lid			

Product	Product Description	Product Status	Product size	Product UOM	Package size	Package UOM	Package type	Liner size	Liner type	Desiccant
MESM-3321-10	MES Monohydrate, 10 KG	In development	10.00 KG		5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
MESM-3321-93	MES Monohydrate, 500G	In development	500.00 G		1000.00	ML	White HDPE bottle, blue lid			
MESM-3350-25	MES Monohydrate, 25 KG	In development	25.00 KG		12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
MESM-3350-92	MES Monohydrate, 100G	In development	100.00 G		250.00	ML	White HDPE bottle, blue lid			
MESM-4220-01	MES Monohydrate, 1 KG	Active	1.00 KG		1500.00	ML	White HDPE bottle, blue lid			
MESM-4220-10	MES Monohydrate, 10 KG	Active	10.00 KG		5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
MESM-4220-25	MES Monohydrate, 25 KG	Active	25.00 KG		12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
MESM-4220-26	MES Monohydrate, 25 KG	Active	25.00 KG		26.00	GAL	Fiber drum	(2x) 34 x 60	LDPE	
MESM-4221-01	MES Monohydrate, Biotech, 1 KG	Active	1.00 KG		1500.00	ML	White HDPE bottle, blue lid			
MESM-4221-05	MES Monohydrate, Biotech, 5 KG	Active	5.00 KG		3.50	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
MESM-4221-10	MES Monohydrate, Biotech, 10 KG	Active	10.00 KG		5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
MESM-4221-25	MES Monohydrate, Biotech, 25 KG	Active	25.00 KG		12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
MESM-4221-50	MES Monohydrate, Biotech, 50 KG	Active	50.00 KG		30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
MESM-4221-92	MES Monohydrate, Biotech, 100 G	Active	100.00 G		250.00	ML	White HDPE bottle, blue lid			
MESM-4221-93	MES Monohydrate, Biotech, 500 G	Active	500.00 G		1000.00	ML	White HDPE bottle, blue lid			
MESM-5201-26	MES Monohydrate, 25 KG	Active	25.00 KG		26.00	GAL	Fiber drum	(2x) 34 x 60	LDPE	
MESN-5220-26	MES Sodium Salt, 25 KG	Active	25.00 KG		26.00	GAL	Fiber drum	(2x) 34 x 60	LDPE	
MOPS-3201-10	MOPS, 10 KG	Active	10.00 KG		5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
MOPS-3201-51	MOPS, 50 KG	Active	50.00 KG		30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
MOPS-3201-93	MOPS, 500 G	Active	500.00 G		1000.00	ML	White HDPE bottle, blue lid			
MOPS-3220-01	MOPS, 1 KG	Active	1.00 KG		1500.00	ML	White HDPE bottle, blue lid			
MOPS-3220-25	MOPS, 25 KG	Active	25.00 KG		12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
MOPS-3220-50	MOPS, 50 KG	Active	50.00 KG		30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
MOPS-3220-51	MOPS, 50 KG	Active	50.00 KG		30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
MOPS-3220-92	MOPS, 100 G	Active	100.00 G		250.00	ML	White HDPE bottle, blue lid			
MOPS-3220-93	MOPS, 500 G	Active	500.00 G		1000.00	ML	White HDPE bottle, blue lid			
MOPS-4220-05	MOPS, 5 KG	Active	5.00 KG		3.50	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
MOPS-4221-01	MOPS, Biotech, 1 KG	Active	1.00 KG		1500.00	ML	White HDPE bottle, blue lid			
MOPS-4221-05	MOPS, Biotech, 5 KG	Active	5.00 KG		3.50	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
MOPS-4221-10	MOPS, Biotech, 10 KG	Active	10.00 KG		5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
MOPS-4221-25	MOPS, Biotech, 25 KG	Active	25.00 KG		12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
MOPS-4221-50	MOPS, Biotech, 50 KG	Active	50.00 KG		30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
MOPS-4221-92	MOPS, Biotech, 100 G	Active	100.00 G		250.00	ML	White HDPE bottle, blue lid			
MOPS-4221-93	MOPS, Biotech, 500 G	Active	500.00 G		1000.00	ML	White HDPE bottle, blue lid			
MOPS-5201-26	MOPS, 50 KG	Active	50.00 KG		26.00	GAL	Fiber drum	(2x) 34 x 60	LDPE	
MOPS-5222-25	MOPS, 25 KG	Active	25.00 KG		12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
NACL-3150-02	Sodium Chloride 5M Solution, 1 L	Active	1.00 L		1.00	L	Tepoflex sample			
NACL-3150-69	Sodium Chloride 5M Solution, 200 L	Active	200.00 L		58.00	GAL	Blue HDPE drum with metal ring for BPC		BPC bag	
NACL-3150-76	Sodium Chloride 5M Solution, 1025 L	Active	1025.00 L		1000.00	L	BPC shell	1000L	BPC bag	
NAHY-4122-70	Sodium Hydroxide 10N Solution, 1 L	Active	1.00 L		1.00	L	Non-fluorinated bottle, black lid			
NAHY-4122-72	Sodium Hydroxide 10N Solution, 200 L	Active	200.00 L		55.00	GAL	Natural HDPE drum			
NAHY-4122-75	Sodium Hydroxide 10N Solution, 940 L	Active	940.00 L		275.00	GAL	Tote			
NAHY-4122-96	Sodium Hydroxide 10N Solution, Bulk	Active	13059.00 L		13059.00	L	Bulk tanker			
NAHY-4123-70	Sodium Hydroxide 2N Solution, 1 L	Active	1.00 L		1.00	L	Non-fluorinated bottle, black lid			
NAHY-4123-71	Sodium Hydroxide 2N Solution, 1 L X 6 combo pack	Active	6.00 L		6.00	L	HDPE combo pack			

Product	Product Description	Product Status	Product size	Product UOM	Package size	Package UOM	Package type	Liner size	Liner type	Desiccant
NAHY-4123-73	Sodium Hydroxide 2N Solution, 4 x 3.78 L	In development	15.12 L		15.12	L	4x1 Gallon combo pack			
NAHY-4124-78	Sodium Hydroxide 5N Solution, 10 L	Active	10.00 L		2.60	GAL	HDPE carboy			
NAHY-4126-71	Sodium Hydroxide 10N Solution, Biotech, 6 X 1 L	Active	6.00 L		6.00	L	HDPE combo pack			
NAHY-4126-72	Sodium Hydroxide 10N Solution, Biotech, 200 L	Active	200.00 L		55.00	GAL	Black HDPE drum			
NAHY-4126-73	Sodium Hydroxide 10N Solution, Biotech, 4 X 3.78 L	Active	15.12 L		15.12	L	4x1 Gallon combo pack			
NAHY-4126-75	Sodium Hydroxide 10N Solution, Biotech, 1000 L	Active	1000.00 L		275.00	GAL	Tote			
NAHY-4126-78	Sodium Hydroxide 10N Solution, Biotech, 10 L	Active	10.00 L		2.60	GAL	HDPE carboy			
NAHY-4150-70	Sodium Hydroxide 10N Solution, 1 L	Active	1.00 L		1.00	L	Non-fluorinated bottle, black lid			
NAHY-4150-71	Sodium Hydroxide 10N Solution, 1 L X 6	Active	6.00 L		6.00	L	HDPE combo pack			
NAHY-4150-73	Sodium Hydroxide 10N Solution, 4 x 3.78 L	Active	15.12 L		15.12	L	4x1 Gallon combo pack			
NAHY-4150-78	Sodium Hydroxide 10N Solution, 10 L	Active	10.00 L		2.60	GAL	HDPE carboy			
NAHY-4150-79	Sodium Hydroxide 10N Solution, 942.5 L	Active	942.50 L		275.00	GAL	Tote			
NAHY-4150-81	Sodium Hydroxide 10N Solution, 4 x 3.78 L	Active	15.12 L		15.12	L	4x1 Gallon combo pack			
NAHY-4150-83	Sodium Hydroxide 10N Solution, 4 x 4 L	Active	16.00 L		16.00	L	4 x 4L combo pack			
NAHY-4151-70	Sodium Hydroxide 8N (25%) Solution, 1 L	Active	1.00 L		1.00	L	Non-fluorinated bottle, black lid			
NDEC-3201-01	Sodium Decanoate, 1 KG	Active	1.00 KG		3.50	GAL	White HDPE round pail, blue lid	(2x) 24 x 36	LDPE	
NDEC-3201-03	Sodium Decanoate, 3 KG	Active	3.00 KG		3.50	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
NDEC-3201-05	Sodium Decanoate, 5 KG	Active	5.00 KG		5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
NDEC-3201-10	Sodium Decanoate, 10 KG	Active	10.00 KG		12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
NDEC-3201-25	Sodium Decanoate, 25 KG	Active	25.00 KG		30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
NDEC-3201-92	Sodium Decanoate, 100 G	Active	100.00 G		1000.00	ML	White HDPE bottle, blue lid			
NDEC-3201-94	Sodium Decanoate, 250 G	Active	250.00 G		1500.00	ML	White HDPE bottle, blue lid			
THCL-3203-01	Tris Hydrochloride, 1 KG	Active	1.00 KG		1500.00	ML	White HDPE bottle, blue lid			
THCL-3203-05	Tris Hydrochloride, 5 KG	Active	5.00 KG		3.50	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
THCL-3203-10	Tris Hydrochloride, 10 KG	Active	10.00 KG		5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
THCL-3203-25	Tris Hydrochloride, 25 KG	Active	25.00 KG		12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
THCL-3203-92	Tris Hydrochloride, 100 G	Active	100.00 G		250.00	ML	White HDPE bottle, blue lid			
THCL-3203-93	Tris Hydrochloride, 500 G	Active	500.00 G		1000.00	ML	White HDPE bottle, blue lid			
THCL-3220-01	Tris Hydrochloride, 1 KG	Active	1.00 KG		1500.00	ML	White HDPE bottle, blue lid			
THCL-3220-05	Tris Hydrochloride, 5 KG	Active	5.00 KG		3.50	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
THCL-3220-10	Tris Hydrochloride, 10 KG	Active	10.00 KG		5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
THCL-3220-12	Tris Hydrochloride, 12 KG	Active	12.00 KG		5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
THCL-3220-13	Tris Hydrochloride, 12 KG	Active	12.00 KG		5.00	GAL	White HDPE round pail, blue lid	(2x) 30 x 48	LDPE	
THCL-3220-25	Tris Hydrochloride, 25 KG	Active	25.00 KG		12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
THCL-3220-50	Tris Hydrochloride, 50 KG	Active	50.00 KG		30.00	GAL	Dark blue HDPE drum	(2x) 34 x 60	LDPE	4 Unit
THCL-3220-58	Tris Hydrochloride, 50 KG	Active	50.00 KG		30.00	GAL	Dark blue HDPE drum	(1x) 34 x 60, (1x)	LDPE/Tyvek	5 x 8 Unit
THCL-3220-92	Tris Hydrochloride, 100 G	Active	100.00 G		250.00	ML	White HDPE bottle, blue lid			
THCL-3220-93	Tris Hydrochloride, 500 G	Active	500.00 G		1000.00	ML	White HDPE bottle, blue lid			
THCL-3250-01	Tris Hydrochloride, 1 KG	Active	1.00 KG		1500.00	ML	White HDPE bottle, blue lid			
THCL-3250-05	Tris Hydrochloride, 5 KG	Active	5.00 KG		3.50	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
THCL-3250-10	Tris Hydrochloride, 10 KG	Active	10.00 KG		5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
THCL-3250-12	Tris Hydrochloride, 12 KG	Active	12.00 KG		5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
THCL-3250-13	Tris Hydrochloride, 12 KG	Active	12.00 KG		5.00	GAL	White HDPE round pail, blue lid	(2x) 30 x 48	LDPE	
THCL-3250-25	Tris Hydrochloride, 25 KG	Active	25.00 KG		12.00	GAL	White HDPE round pail, blue lid	(2x) 30 x 48	LDPE	
THCL-3250-28	Tris Hydrochloride, 25 KG	Active	25.00 KG		12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	

Product	Product Description	Product Status	Product size	Product UOM	Package size	Package UOM	Package type	Liner size	Liner type	Desiccant
THCL-3250-41	Tris Hydrochloride, 10 KG	Active	10.00 KG		5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
THCL-3250-42	Tris Hydrochloride, 12 KG	Active	12.00 KG		5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
THCL-3250-58	Tris Hydrochloride, 50 KG	Active	50.00 KG		30.00	GAL	Dark blue HDPE drum	(1x) 34 x 60, (1x)	LDPE/Tyvek	5 x 8 Unit
THCL-3250-91	Tris Hydrochloride, 50 G	Active	50.00 G		250.00	ML	White HDPE bottle, blue lid			
THCL-3250-92	Tris Hydrochloride, 100 G	Active	100.00 G		250.00	ML	White HDPE bottle, blue lid			
THCL-3250-93	Tris Hydrochloride, 500 G	Active	500.00 G		1000.00	ML	White HDPE bottle, blue lid			
THCL-3253-10	Tris Hydrochloride, 10 KG	Active	10.00 KG		5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
THCL-3253-12	Tris Hydrochloride, 12 KG	Active	12.00 KG		5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
THCL-3253-13	Tris Hydrochloride, 12 KG	Active	12.00 KG		5.00	GAL	White HDPE round pail, blue lid	(2x) 30 x 48	LDPE	
THCL-3253-92	Tris Hydrochloride, 100 G	Active	100.00 G		250.00	ML	White HDPE bottle, blue lid			
THCL-3253-93	Tris Hydrochloride, 500 G	Active	500.00 G		1000.00	ML	White HDPE bottle, blue lid			
THCL-3254-01	Tris Hydrochloride, 1 KG	Active	1.00 KG		1500.00	ML	White HDPE bottle, blue lid			
THCL-3254-12	Tris Hydrochloride, 12 KG	Active	12.00 KG		5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
THCL-3254-13	Tris Hydrochloride, 12 KG	Active	12.00 KG		5.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
THCL-3254-25	Tris Hydrochloride, 25 KG	Active	25.00 KG		12.00	GAL	White HDPE round pail, blue lid	(2x) 30 x 48	LDPE	
THCL-3254-50	Tris Hydrochloride, 50 KG	Active	50.00 KG		30.00	GAL	Dark blue HDPE drum	(2x) 34 x 60	LDPE	4 Unit
THCL-3254-92	Tris Hydrochloride, 100 G	Active	100.00 G		250.00	ML	White HDPE bottle, blue lid			
THCL-3254-93	Tris Hydrochloride, 500 G	Active	500.00 G		1000.00	ML	White HDPE bottle, blue lid			
THCL-3257-10	Tris Hydrochloride, 10 KG	Active	10.00 KG		5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
THCL-3257-93	Tris Hydrochloride, 500 G	Active	500.00 G		1000.00	ML	White HDPE bottle, blue lid			
THCL-3259-01	Tris Hydrochloride, 1 KG	Active	1.00 KG		1500.00	ML	White HDPE bottle, blue lid			
THCL-3259-05	Tris Hydrochloride, 5 KG	Active	5.00 KG		3.50	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
THCL-3259-10	Tris Hydrochloride, 10 KG	Active	10.00 KG		5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
THCL-3259-12	Tris Hydrochloride, 12 KG	Active	12.00 KG		5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
THCL-3259-25	Tris Hydrochloride, 25 KG	Active	25.00 KG		12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
THCL-3259-58	Tris Hydrochloride, 50 KG	Active	50.00 KG		30.00	GAL	Dark blue HDPE drum	(1x) 34 x 60, (1x)	LDPE/Tyvek	5 x 8 Unit
THCL-3259-91	Tris Hydrochloride, 250 G	Active	250.00 G		1000.00	ML	White HDPE bottle, blue lid			
THCL-3259-92	Tris Hydrochloride, 100 G	Active	100.00 G		250.00	ML	White HDPE bottle, blue lid			
THCL-3259-93	Tris Hydrochloride, 500 G	Active	500.00 G		1000.00	ML	White HDPE bottle, blue lid			
THCL-3260-01	Tris Hydrochloride, 1 KG	Active	1.00 KG		1500.00	ML	White HDPE bottle, blue lid			
THCL-3260-92	Tris Hydrochloride, 100 G	Active	100.00 G		250.00	ML	White HDPE bottle, blue lid			
THCL-4250-92	Tris Hydrochloride, 100 G	Active	100.00 G		250.00	ML	White HDPE bottle, blue lid			
THCL-4251-01	Tris Hydrochloride, Biotech, 1 KG	Active	1.00 KG		1500.00	ML	White HDPE bottle, blue lid			
THCL-4251-05	Tris Hydrochloride, Biotech, 5 KG	Active	5.00 KG		3.50	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
THCL-4251-10	Tris Hydrochloride, Biotech, 10 KG	Active	10.00 KG		5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
THCL-4251-25	Tris Hydrochloride, Biotech, 25 KG	Active	25.00 KG		12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
THCL-4251-50	Tris Hydrochloride, Biotech, 50 KG	In development	50.00 KG		30.00	GAL	Dark blue HDPE drum	(2x) 34 x 60	LDPE	4 Unit
THCL-4251-92	Tris Hydrochloride, Biotech, 100 G	Active	100.00 G		250.00	ML	White HDPE bottle, blue lid			
THCL-4251-93	Tris Hydrochloride, Biotech, 500 G	Active	500.00 G		1000.00	ML	White HDPE bottle, blue lid			
THCL-6201-01	Tris Hydrochloride, 1KG	Active	1.00 KG		1500.00	ML	White HDPE bottle, blue lid			
THCL-6201-92	Tris Hydrochloride, 100 G	Active	100.00 G		250.00	ML	White HDPE bottle, blue lid			
THCL-6250-92	Tris Hydrochloride, 100 G	Active	100.00 G		250.00	ML	White HDPE bottle, blue lid			
TRED-3250-03	Trehalose Dihydrate, 3 KG	Active	3.00 KG		3.50	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
TRED-3250-05	Trehalose Dihydrate, 5 KG	Active	5.00 KG		3.50	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
TRED-3250-06	Trehalose Dihydrate, 1 KG	Active	1.00 KG		0.60	GAL	White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	

Product	Product Description	Product Status	Product size	Product UOM	Package size	Package UOM	Package type	Liner size	Liner type	Desiccant
TRED-3250-10	Trehalose Dihydrate, 10 KG	Active	10.00 KG		5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
TRED-3250-25	Trehalose Dihydrate, 25 KG	Active	25.00 KG		12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
TRED-3250-50	Trehalose Dihydrate, 50 KG	Active	50.00 KG		30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
TRED-3250-92	Trehalose Dihydrate, 100 G	Active	100.00 G		250.00	ML	White HDPE bottle, blue lid			
TRED-3252-06	Trehalose Dihydrate, 1 KG	Active	1.00 KG		0.60	GAL	White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	
TRED-3252-10	Trehalose Dihydrate, 10 KG	Active	10.00 KG		5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
TRED-3252-25	Trehalose Dihydrate, 25 KG	Active	25.00 KG		12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
TRED-3252-92	Trehalose Dihydrate, 100 G	Active	100.00 G		250.00	ML	White HDPE bottle, blue lid			
TRED-4250-05	Trehalose Dihydrate, Biotech, 5 KG	Active	5.00 KG		3.50	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
TRED-4250-06	Trehalose Dihydrate, Biotech, 1 KG	Active	1.00 KG		0.60	GAL	White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	
TRED-4250-10	Trehalose Dihydrate, Biotech, 10 KG	Active	10.00 KG		5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
TRED-4250-25	Trehalose Dihydrate, Biotech, 25 KG	Active	25.00 KG		12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
TRED-4250-50	Trehalose Dihydrate, Biotech, 50 KG	Active	50.00 KG		30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
TRIS-1201-01	Tris/Tromethamine, 1 KG	Active	1.00 KG		1500.00	ML	White HDPE bottle, blue lid			
TRIS-1201-10	Tris/Tromethamine, 10 KG	Active	10.00 KG		5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
TRIS-1201-25	Tris/Tromethamine, 25 KG	Active	25.00 KG		12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
TRIS-1201-91	Tris/Tromethamine, 250 G	Active	250.00 G		1000.00	ML	White HDPE bottle, blue lid			
TRIS-1201-92	Tris/Tromethamine, 100 G	Active	100.00 G		250.00	ML	White HDPE bottle, blue lid			
TRIS-1201-93	Tris/Tromethamine, 500 G	Active	500.00 G		1000.00	ML	White HDPE bottle, blue lid			
TRIS-1220-01	Tris/Tromethamine, 1 KG	Active	1.00 KG		1500.00	ML	White HDPE bottle, blue lid			
TRIS-1220-25	Tris/Tromethamine, 25 KG	Active	25.00 KG		12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
TRIS-1220-91	Tris/Tromethamine, 250 G	Active	250.00 G		1000.00	ML	White HDPE bottle, blue lid			
TRIS-1220-92	Tris/Tromethamine, 100 G	Active	100.00 G		250.00	ML	White HDPE bottle, blue lid			
TRIS-1220-93	Tris/Tromethamine, 500 G	Active	500.00 G		1000.00	ML	White HDPE bottle, blue lid			
TRIS-2255-01	Tris/Tromethamine, 1 KG	Active	1.00 KG		1500.00	ML	White HDPE bottle, blue lid			
TRIS-2255-10	Tris/Tromethamine, 10 KG	Active	10.00 KG		5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
TRIS-2255-91	Tris/Tromethamine, 250 G	Active	250.00 G		1000.00	ML	White HDPE bottle, blue lid			
TRIS-2255-92	Tris/Tromethamine, 100 G	Active	100.00 G		250.00	ML	White HDPE bottle, blue lid			
TRIS-2255-93	Tris/Tromethamine, 500 G	Active	500.00 G		1000.00	ML	White HDPE bottle, blue lid			
TRIS-2257-01	Tris/Tromethamine, 1 KG	Active	1.00 KG		1500.00	ML	White HDPE bottle, blue lid			
TRIS-2257-25	Tris/Tromethamine, 25 KG	Active	25.00 KG		12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
TRIS-2257-91	Tris/Tromethamine, 250 G	Active	250.00 G		1000.00	ML	White HDPE bottle, blue lid			
TRIS-2257-92	Tris/Tromethamine, 100 G	Active	100.00 G		250.00	ML	White HDPE bottle, blue lid			
TRIS-2257-93	Tris/Tromethamine, 500 G	Active	500.00 G		1000.00	ML	White HDPE bottle, blue lid			
TRIS-3220-01	Tris/Tromethamine, 1 KG	Active	1.00 KG		1500.00	ML	White HDPE bottle, blue lid			
TRIS-3220-10	Tris/Tromethamine, 10 KG	Active	10.00 KG		5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
TRIS-3220-25	Tris/Tromethamine, 25 KG	Active	25.00 KG		12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
TRIS-3220-50	Tris/Tromethamine, 50 KG	Active	50.00 KG		30.00	GAL	Dark blue HDPE drum	(2x) 34 x 60	LDPE	4 Unit
TRIS-3220-92	Tris/Tromethamine, 100 G	Active	100.00 G		250.00	ML	White HDPE bottle, blue lid			
TRIS-3220-93	Tris/Tromethamine, 500 G	Active	500.00 G		1000.00	ML	White HDPE bottle, blue lid			
TRIS-3251-01	Tris/Tromethamine, 1 KG	Active	1.00 KG		1500.00	ML	White HDPE bottle, blue lid			
TRIS-3251-05	Tris/Tromethamine, 5 KG	Active	5.00 KG		3.50	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
TRIS-3251-10	Tris/Tromethamine, 10 KG	Active	10.00 KG		5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
TRIS-3251-92	Tris/Tromethamine, 100 G	Active	100.00 G		250.00	ML	White HDPE bottle, blue lid			
TRIS-3251-95	Tris/Tromethamine, 80 G	Active	80.00 G		1000.00	ML	White HDPE bottle, blue lid	9 x 15	LDPE	

Product	Product Description	Product Status	Product size	Product UOM	Package size	Package UOM	Package type	Liner size	Liner type	Desiccant
TRIS-3252-25	Tris/Tromethamine, 25 KG	Active	25.00 KG		12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
TRIS-3252-60	Tris/Tromethamine, 50 KG	Active	50.00 KG		30.00	GAL	Dark blue HDPE drum	(1x) 34 x 60, (2x)	LDPE/(2x) Tyve	(2x) 5 x 8 Unit
TRIS-3254-01	Tris/Tromethamine, 1KG	Active	1.00 KG		1500.00	ML	White HDPE bottle, blue lid			
TRIS-3254-05	Tris/Tromethamine, 5 KG	Active	5.00 KG		3.50	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
TRIS-3254-10	Tris/Tromethamine, 10 KG	Active	10.00 KG		5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
TRIS-3254-25	Tris/Tromethamine, 25 KG	Active	25.00 KG		12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
TRIS-3254-50	Tris/Tromethamine, 50 KG	Active	50.00 KG		30.00	GAL	Dark blue HDPE drum	(2x) 34 x 60	LDPE	4 Unit
TRIS-3254-58	Tris/Tromethamine, 50 KG	Active	50.00 KG		30.00	GAL	Dark blue HDPE drum	(1x) 34 x 60, (1x)	LDPE/Tyvek	5 x 8 Unit
TRIS-3254-92	Tris/Tromethamine, 100 G	Active	100.00 G		250.00	ML	White HDPE bottle, blue lid			
TRIS-3254-93	Tris/Tromethamine, 500 G	Active	500.00 G		1000.00	ML	White HDPE bottle, blue lid			
TRIS-3255-01	Tris/Tromethamine, 1 KG	Active	1.00 KG		1500.00	ML	White HDPE bottle, blue lid			
TRIS-3255-05	Tris/Tromethamine, 5 KG	Active	5.00 KG		3.50	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
TRIS-3255-07	Tris/Tromethamine, 2 KG	Active	2.00 KG		1.25	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
TRIS-3255-10	Tris/Tromethamine, 10 KG	Active	10.00 KG		5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
TRIS-3255-12	Tris/Tromethamine, 12 KG	Active	12.00 KG		5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
TRIS-3255-25	Tris/Tromethamine, 25 KG	Active	25.00 KG		12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
TRIS-3255-41	Tris/Tromethamine, 10 KG	Active	10.00 KG		5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
TRIS-3255-42	Tris/Tromethamine, 12 KG	Active	12.00 KG		5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
TRIS-3255-50	Tris/Tromethamine, 50 KG	Active	50.00 KG		30.00	GAL	Dark blue HDPE drum	(2x) 34 x 60	LDPE	4 Unit
TRIS-3255-58	Tris/Tromethamine, 50 KG	Active	50.00 KG		30.00	GAL	Dark blue HDPE drum	(1x) 34 x 60, (1x)	LDPE/Tyvek	5 x 8 Unit
TRIS-3255-89	Tris/Tromethamine, 200 G	Active	200.00 G		1000.00	ML	White HDPE bottle, blue lid			
TRIS-3255-91	Tris/Tromethamine, 50 G	Active	50.00 G		250.00	ML	White HDPE bottle, blue lid			
TRIS-3255-92	Tris/Tromethamine, 100 G	Active	100.00 G		250.00	ML	White HDPE bottle, blue lid			
TRIS-3255-93	Tris/Tromethamine, 500 G	Active	500.00 G		1000.00	ML	White HDPE bottle, blue lid			
TRIS-3255-94	Tris/Tromethamine, 250 G	Active	250.00 G		1000.00	ML	White HDPE bottle, blue lid			
TRIS-3257-10	Tris/Tromethamine, 10 KG	Active	10.00 KG		5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
TRIS-4223-01	Tris/Tromethamine, Biotech, 1 KG	Active	1.00 KG		1500.00	ML	White HDPE bottle, blue lid			
TRIS-4223-05	Tris/Tromethamine, Biotech, 5 KG	Active	5.00 KG		3.50	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
TRIS-4223-10	Tris/Tromethamine, Biotech, 10 KG	Active	10.00 KG		5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
TRIS-4223-25	Tris/Tromethamine, Biotech, 25 KG	Active	25.00 KG		12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
TRIS-4223-50	Tris/Tromethamine, Biotech, 50 KG	In development	50.00 KG		30.00	GAL	Dark blue HDPE drum	(2x) 34 x 60	LDPE	4 Unit
TRIS-4223-92	Tris/Tromethamine, Biotech, 100 G	Active	100.00 G		250.00	ML	White HDPE bottle, blue lid			
TRIS-4223-93	Tris/Tromethamine, Biotech, 500 G	Active	500.00 G		1000.00	ML	White HDPE bottle, blue lid			
TRIS-4224-92	Tris/Tromethamine, 100 G	Active	100.00 G		250.00	ML	White HDPE bottle, blue lid			
TRIS-4252-01	Tris/Tromethamine, 1 KG	Active	1.00 KG		1500.00	ML	White HDPE bottle, blue lid			
TRIS-4252-92	Tris/Tromethamine, 100 G	Active	100.00 G		250.00	ML	White HDPE bottle, blue lid			
TRIS-5201-26	Tris, 25 KG	Active	25.00 KG		26.00	GAL	Fiber drum	(2x) 34 x 60	LDPE	4 Unit
TRIS-5204-26	Tris, 25 KG	Active	25.00 KG		26.00	GAL	Fiber drum	(2x) 34 x 60	LDPE	
TRIS-5207-26	Tris, 25 KG	Active	25.00 KG		26.00	GAL	Fiber drum	(2x) 34 x 60	LDPE	4 Unit
TRIS-5222-10	Tris, 10 KG	Active	10.00 KG		5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
TRIS-5223-89	Tris, 200 G	Active	200.00 G		1000.00	ML	White HDPE bottle, blue lid			
TRIS-6201-01	Tris, 1 KG	Active	1.00 KG		1500.00	ML	White HDPE bottle, blue lid			
TRIS-6201-91	Tris, 250 G	Active	250.00 G		1000.00	ML	White HDPE bottle, blue lid			
TRIS-6201-92	Tris, 100 G	Active	100.00 G		250.00	ML	White HDPE bottle, blue lid			
TRIS-6202-92	Tris, 100 G	Active	100.00 G		250.00	ML	White HDPE bottle, blue lid			

Product	Product Description	Product Status	Product size	Product UOM	Package size	Package UOM	Package type	Liner size	Liner type	Desiccant
TRIS-6252-92	Tris, 100 G	In development	100.00 G		250.00	ML	White HDPE bottle, blue lid			
URAC-4201-05	Uracil, 5 KG	Active	5.00 KG		3.50	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
URAC-4202-01	Uracil, 1 KG	Active	1.00 KG		1500.00	ML	White HDPE bottle, blue lid			
URAC-4202-05	Uracil, 5 KG	Active	5.00 KG		3.50	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
URAC-4202-91	Uracil, 50 G	Active	50.00 G		250.00	ML	White HDPE bottle, blue lid			
URAC-4202-92	Uracil, 100 G	Active	100.00 G		250.00	ML	White HDPE bottle, blue lid			
URAC-4202-95	Uracil, 10g	Active	10.00 G		60.00	ML	White polypropylene bottle			
URAC-4250-01	Uracil, 1 KG	Active	1.00 KG		1500.00	ML	White HDPE bottle, blue lid			
URAC-4250-95	Uracil, 10g	Active	10.00 G		60.00	ML	White polypropylene bottle			
URAC-4301-05	Uracil, 5 KG	Active	5.00 KG		3.50	GAL	White HDPE round pail, blue lid	(2x) 24 x 36	LDPE	
URAC-4301-06	Uracil, 5 KG	Active	5.00 KG		3.50	GAL	White HDPE round pail, blue lid	(5x) 24 x 36	LDPE	
UREA-3120-02	Urea 6M Solution, 500 ML	Active	500.00 ML		17.00	OZ	Black HDPE UN leakproof bottle			
UREA-3120-76	Urea 6M Solution, 750 L	In development	750.00 L		1000.00	L	BPC shell	1000L	BPC bag	
UREA-3120-78	Urea 6M Solution, 1200 L	Active	1200.00 L		330.00	GAL	Tote		LDPE	
UREA-3220-10	Urea, 10 KG	Active	10.00 KG		5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
UREA-3220-25	Urea, 25 KG	Active	25.00 KG		12.00	GAL	White HDPE round pail, blue lid	34 x 60	LDPE	
UREA-3220-50	Urea, 50 KG	Active	50.00 KG		30.00	GAL	Dark blue HDPE drum	(1x) 34 x 60, (1x)	LDPE/Tyvek	5 x 8 Unit
UREA-3220-89	Urea, 200 G	Active	200.00 G		1000.00	ML	White HDPE bottle, blue lid			
UREA-3220-92	Urea, 100 G	Active	100.00 G		250.00	ML	White HDPE bottle, blue lid			
UREA-3220-93	Urea, 500 G	Active	500.00 G		1000.00	ML	White HDPE bottle, blue lid			
UREA-3223-25	Urea, 25 KG	Active	25.00 KG		12.00	GAL	White HDPE round pail, blue lid	34 x 60	LDPE	
UREA-3223-58	Urea, 50 KG	Active	50.00 KG		30.00	GAL	Dark blue HDPE drum	(1x) 34 x 60, (1x)	LDPE/Tyvek	5 x 8 Unit
UREA-3250-01	Urea, 1 KG	Active	1.00 KG		1500.00	ML	White HDPE bottle, blue lid			
UREA-3250-12	Urea, 12 KG	Active	12.00 KG		5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	4 Unit
UREA-3250-25	Urea, 25 KG	Active	25.00 KG		12.00	GAL	White HDPE round pail, blue lid	34 x 60	LDPE	
UREA-3250-92	Urea, 100 G	Active	100.00 G		250.00	ML	White HDPE bottle, blue lid			
UREA-3250-93	Urea, 500 G	Active	500.00 G		1000.00	ML	White HDPE bottle, blue lid			
UREA-4221-01	Urea, Biotech, 1 KG	Active	1.00 KG		1500.00	ML	White HDPE bottle, blue lid			
UREA-4221-05	Urea, Biotech, 5 KG	Active	5.00 KG		3.50	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
UREA-4221-10	Urea, Biotech, 10 KG	Active	10.00 KG		5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
UREA-4221-25	Urea, Biotech, 25 KG	Active	25.00 KG		12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
UREA-4221-50	Urea, Biotech, 50 KG	Active	50.00 KG		30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
UREA-4221-92	Urea, Biotech, 100 G	Active	100.00 G		250.00	ML	White HDPE bottle, blue lid			
UREA-4221-93	Urea, Biotech, 500 G	Active	500.00 G		1000.00	ML	White HDPE bottle, blue lid			
UREA-5203-05	Urea, 5 KG	Active	5.00 KG		3.50	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
URID-3250-01	Uridine, 1 KG	Active	1.00 KG		1500.00	ML	White HDPE bottle, blue lid			
URID-3250-10	Uridine, 10 KG	Active	10.00 KG		5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
URID-3250-25	Uridine, 25 KG	Active	25.00 KG		12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	
URID-3250-91	Uridine, 250 G	Active	250.00 G		1000.00	ML	White HDPE bottle, blue lid			
URID-3250-92	Uridine, 100 G	Active	100.00 G		250.00	ML	White HDPE bottle, blue lid			
URID-3250-93	Uridine, 500 G	Active	500.00 G		1000.00	ML	White HDPE bottle, blue lid			
URID-4250-01	Uridine, Biotech, 1 KG	Active	1.00 KG		1500.00	ML	White HDPE bottle, blue lid			
URID-4250-05	Uridine, Biotech, 5KG	Active	5.00 KG		3.50	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
URID-4250-10	Uridine, Biotech, 10 KG	Active	10.00 KG		5.00	GAL	White HDPE round pail, blue lid	24 x 36	LDPE	
URID-4250-25	Uridine, Biotech, 25 KG	Active	25.00 KG		12.00	GAL	White HDPE round pail, blue lid	30 x 48	LDPE	

Product	Product Description	Product Status	Product size	Product UOM	Package size	Package UOM	Package type	Liner size	Liner type	Desiccant
URID-4250-50	Uridine, Biotech, 50 KG	Active	50.00 KG		30.00	GAL	Dark blue HDPE drum	34 x 60	LDPE	
URID-4250-92	Uridine, Biotech, 100 G	Active	100.00 G		250.00	ML	White HDPE bottle, blue lid			
URID-4250-93	Uridine, Biotech, 500 G	Active	500.00 G		1000.00	ML	White HDPE bottle, blue lid			
URID-4250-95	Uridine, Biotech, 50 G	Active	100.00 G		250.00	ML	White HDPE bottle, blue lid			
WAFI-3150-68	Water for Injection, 200L	In development	200.00 L		200.00	L	BPC shell	200L	Typoflex BPC b	
WAFI-3150-74	Water for Injection, 3.78L	Active	3.78 L		3.78	L	Non-fluorinated HDPE bottle			