

HEPES

BET*, GMP

CAS #: 7365-45-9

Formula: $C_8H_{18}N_2O_4S$

F.W.: 238.30 g/mol

HEPE-4250

BIO PHARMA GRADE

ANALYSIS		SPECIFICATIONS
Absorbance (0.1M)	250 nm	≤ 0.0500 a.u.
	260 nm	≤ 0.0500 a.u.
	280 nm	≤ 0.0800 a.u.
Absorbance (0.05M)	250 nm	≤ 0.0500 a.u.
	260 nm	≤ 0.0500 a.u.
	280 nm	≤ 0.0800 a.u.
Appearance and Color		White Crystals
Assay, Dried Basis		$\geq 99.5\%$
Chloride		≤ 50 ppm
Endotoxin		≤ 5 EU/g
Enzymes	DNase	None Detected
	RNase	None Detected
	Protease	None Detected
Identification, IR		Conforms to Reference Standard
Insoluble Matter		$\leq 0.01\%$
Loss on Drying		$\leq 0.5\%$
Microbial Content	TAMC	≤ 100 CFU/g
	TYMC	≤ 100 CFU/g
pH (5%)		5.0 – 6.5
pK _a		7.45 – 7.65
Residue on Ignition		$\leq 0.1\%$
Solubility	Solubility (5%)	Passes Test
	Solubility (0.05M)	Passes Test
Sulfate		≤ 50 ppm
Trace Metals	Arsenic (As)	≤ 5 ppm
	Aluminum (Al)	≤ 5 ppm
	Bismuth (Bi)	≤ 5 ppm
	Calcium (Ca)	≤ 10 ppm
	Copper (Cu)	≤ 5 ppm
	Iron (Fe)	≤ 5 ppm
	Lead (Pb)	≤ 1 ppm
	Lithium (Li)	≤ 5 ppm
	Molybdenum (Mo)	≤ 5 ppm
	Nickel (Ni)	≤ 5 ppm
	Potassium (K)	≤ 50 ppm
Water, KF		$\leq 0.1\%$

General Product Overview

HEPES, Free Acid is a zwitterionic buffer used to maintain pH of media used in cell cultures. It is one of Good's buffers that has a pK_a value similar to its pH value, making it an ideal buffer for pH maintenance. HEPES is a Good's buffer because it has low UV absorptivity, minimal reactivity, stable pH and is soluble in water.

Industry Application

Suitable for use in biological and biotech chemical process applications from R&D through scale cGMP production.

Key Product Features

- Appears as a white crystalline product
- Manufactured in accordance with IPEC
- Manufactured in an enzyme free, hormone free and animal free environment
- Contains no known major food allergens (as defined by FDA and WHO)
- The final product and its raw materials are not derived from nor come into contact with animal parts, animal products, and/or animal byproducts or derivatives.
- Is not subject to genetic modification
- Synonyms: N-(2-Hydroxyethyl) Piperazine-N'-2-Ethanesulfonic Acid; 4-(2-Hydroxyethyl) Piperazine-1-Ethanesulfonic Acid

Storage and Shipping Conditions

Refer to SDS.

Standard Shelf-Life Policy

Unless otherwise noted on the Shelf-Life Statement and CoA, this product has a 2-year retest date supported by a 3-year ICH Q1 Stability Study (if one is completed).

Package Sizes

100g, 500g, 1kg, 5kg, 10kg, 25kg, 50kg

[Click here to view SDS, CoAs and other supporting regulatory documents on our website.](#)

*Bioburden and Endotoxin Tested

This is not considered a controlled document. We are not responsible for any errors or omissions, and the user is responsible for any decisions based on the information herein.