

Packaging - Specification
container for transport


SCHÜTZ

MX 330 UN / Nat / 6" red TP, 2" solid plug
Ball-valve, FKM gasket, 2" male camlock (1pc) , steel frame
2 plates, ticket

SCHÜTZ
200 Aspen Hill Road
North Branch, NJ 08876

Date: September 17, 2018
Page 1 of 2

Material number **MX 1208233**
 RMX 1096583



Technical data:

| | | |
|---------------------------|-----|-------|
| Rated volume: | 330 | gal |
| Overflow volume: | 334 | gal |
| Length: | 48 | in |
| Width: | 40 | in |
| Height with pallet: | 53 | in |
| Filling opening: | 6 | in |
| Discharge opening: | 2 | in |
| Fork opening - height | 3.5 | in |
| Label plate: | 2 | piece |
| Corner Protector: | 4 | piece |
| Weight: approx | 146 | lbs |

Nominal Specifications

Packaging - Specification
container for transport



MX 330 UN / Nat / 6" red TP, 2" solid plug
Ball-valve, FKM gasket, 2" male camlock (1pc), steel frame
2 plates, ticket

SCHÜTZ
200 Aspen Hill Road
North Branch, NJ 08876

Date: September 17, 2018
Page 2 of 2

Material number **MX 1208233**
 RMX 1096583

Construction:

Outer Container: Rectangular grid box made from tubular steel material, with bottom plate, closed with tie-bar, label plate with Schütz -Ticket on front side, additional label plate on the back side.

Material: Grid / Bottom Plate: Steel, galvanized against corrosion
 Corner Protector: 4 pieces made of HDPE, black

Inner Container: Rectangular blow molded tank of high density polyethylene (HDPE), with filling opening (6") in the middle of the top section. Valve opening in front section.

Material: High Density Polyethylene – natural

Filling Opening: Filling opening 6" with external thread, closed with red screw cap, 2" solid plug.

Material: Screw cap 6": HDPE, red
 O-ring gasket: TP
 G2-plug: Polypropylene
 O-ring 2" plug: Lucofin

Discharge Opening: Ball valve, FKM gaskets, 2" male camlock outlet, 1 pc dust cap.

Material: Housing ball-valve: HDPE
 Ball inside valve Polypropylene
 Flat gasket: FKM
 Ball gasket: FKM
 Screw cap: HDPE

Pallet: **Steel frame pallet** (1000 x 1200 mm), 4-way entry

Material: Galvanized steel

Colorant: Heavy-metal free

Delivery: Ready for filling, clean

UN-Marking: UN 31 HA1 / Y / USA / SCHUTZ

Customer Acceptance Signature: _____ **Date:** _____

Company name: _____