

## 19L CARBOY SPECIFICATIONS



2.6 Gallon / 10L HDPE Carboy	
Test	Specifications
Capacity	2.6 Gallons
Material	Natural HDPE Jerry Can
Appearance and Color	White, rectangular with built-in handle
Outside Length (Internal Tolerance)	9 .1 inches $\pm$ 1 inch
Outside Width (Internal Tolerance)	7.5 inches $\pm$ 1 inch
Outside Height (Internal Tolerance)	12.1 inches $\pm$ 1 inch
Closure	Continuous Thread Round Screw Top

2.6 Gallon / 10L HDPE Carboy Cap	
Test	Specifications
Material	50mm Black HDPE Cap
Appearance and Color	Black Round Tamper Evident Screw Cap

For further information, please contact [marketing@biospectra.us](mailto:marketing@biospectra.us).



**Connie Augustine | Director**  
**Business Integration & Retention**  
**BioSpectra Inc.**

*This is not considered a controlled document. We are not responsible for any errors or omissions, and the user is responsible for any decisions based on the information herein.*