

12/19/2025

Version 1.1

4x4L BOTTLE SPECIFICATIONS

BSI COMPONENT NUMBER: COMP-0305



4x4L Liter Non-Fluorinated HDPE Bottle Combo Pack	
Test	Specifications
Appearance and Color	HDPE Bottle Combo Pack Natural Round Plastic Bottle
Capacity	4 Bottles x 4 Liters Each
Outside Height	13.75 inches \pm 0.5 inches
Outside Width	6.0 inches \pm 0.5 inches
Air Insert (Bottom of bottle)	Arrow pointed at A
UN Rating	4G/Y27.6/S

For further information, please contact marketing@biospectra.us.



Connie Augustine | Director
Business Integration & Retention
BioSpectra Inc.

This is not considered a controlled document. We are not responsible for any errors or omissions, and the user is responsible for any decisions based on the information herein.