

10L CARBOY SPECIFICATIONS

BSI COMPONENT NUMBER: COMP-0042



| 2.6 Gallon / 10L HDPE Carboy | |
|-------------------------------------|---|
| Test | Specifications |
| Capacity | 2.6 Gallons |
| Material | Natural HDPE Jerry Can |
| Appearance and Color | White, rectangular with built-in handle |
| Outside Length (Internal Tolerance) | 9 .1 inches \pm 1 inch |
| Outside Width (Internal Tolerance) | 7.5 inches \pm 1 inch |
| Outside Height (Internal Tolerance) | 12.1 inches \pm 1 inch |
| Closure | Continuous Thread Round Screw Top |

| 2.6 Gallon / 10L HDPE Carboy Cap | |
|----------------------------------|--------------------------------------|
| Test | Specifications |
| Material | 50mm Black HDPE Cap |
| Appearance and Color | Black Round Tamper Evident Screw Cap |

For further information, please contact marketing@biospectra.us.

Connie Augustine | Director
Business Integration & Retention
BioSpectra Inc.

This is not considered a controlled document. We are not responsible for any errors or omissions, and the user is responsible for any decisions based on the information herein.