

TRIS HYDROCHLORIDE

ENDOTOXIN TESTED, GMP

CAS #: 1185-53-1

Formula: $\text{NH}_2\text{C}(\text{CH}_2\text{OH})_3 \cdot \text{HCl}$

F.W.: 157.60 g/mol

THCL-3253

BIO EXCIPIENT GRADE

ANALYSIS		SPECIFICATIONS
Absorbance (1M)	260 nm	≤ 0.06 a.u.
	280 nm	≤ 0.06 a.u.
	400 nm	≤ 0.01 a.u.
Appearance and Color		White to Off White Crystalline Powder or Granules
Assay, Dried Basis		$\geq 99.0\%$
Endotoxin		≤ 2.5 EU/g
Identification, IR		Conforms to Reference Standard
Loss on Drying		$\leq 0.5\%$
Melting Range		150 – 153°C
pH (0.5M)		3.5 - 5.0
Trace Metals	Arsenic (As)	≤ 10 ppm
	Calcium (Ca)	≤ 10 ppm
	Copper (Cu)	≤ 10 ppm
	Iron (Fe)	≤ 10 ppm
	Lead (Pb)	≤ 10 ppm
	Magnesium (Mg)	≤ 10 ppm

General Product Overview

Tris Hydrochloride is a stabilizing buffer in biological applications such as electrochromatography, UV analysis and HPLC. It is used to adjust and stabilize the pH ranges for gels used in electrophoresis applications. Tris Hydrochloride is extensively used as a biological buffer or a component of buffer solutions.

Industry Application

Suitable for use as a cGMP chemical in pharmaceutical manufacturing processes and products.

Key Product Features

- The manufacturing of Tris Hydrochloride, THCL-3253 is performed at BioSpectra's Stroudsburg, PA facility.
- Appears as white crystals
- Manufactured in accordance with ICH Q7
- Manufactured in an enzyme free, hormone free and animal free environment
- Contains no known major food allergens (as defined by FDA and WHO)
- The final product and its raw materials are not derived from nor come into contact with animal parts, animal products, and/or animal byproducts or derivatives.
- Is not subject to genetic modification
- Synonyms: 2-Amino-2-(Hydroxymethyl)propane-1,3-diol hydrochloride; Tri(Hydroxymethyl) Aminomethane Hydrochloride

Storage and Shipping Conditions

Refer to SDS.

Standard Shelf-Life Policy

Unless otherwise noted on the Shelf-Life Statement and CoA, this product has a 2-year retest date supported by a 3-year ICH Q1 Stability Study (if one is completed).

Package Sizes

10kg, 25kg and 50kg pails

[Click here to view SDS, CoAs and other supporting regulatory documents on our website.](#)

This is not considered a controlled document. We are not responsible for any errors or omissions, and the user is responsible for any decisions based on the information herein.

IUO: {LIL/SAM}