

# Desiccare, Inc. Product Specification Sheet

## 8 Unit Clay Unit-Pak Strip

Desiccare Part Number: 10BQ10A16-S

Product: Unit-Pak Unit Bag Desiccator w/ 264 grams of Montmorillonite (bentonite) Clay

Desiccant: Clay Particles 3 x 80 mesh

Description: Montmorillonite (bentonite) Clay in a printed Tyvek packet

Packet Specifications:

Clay	Net Weight:	1320 grams	±10%
	Packet Length:	42.500"	±0.250"
	Packet Width:	4.500"	±0.125"

Moisture Content: Less than 3.00% maximum residual moisture by weight (loss-on-drying)

Packet Substrate: Tyvek, spunbonded olefin. Manufactured by E.I. DuPont  
Exterior - Uncoated  
Interior - heat seal coating (Hot Melt Adhesive)

Print: Conforms with Specification MIL-D-3464 Type I & II  
Reactivation Time in bag 16 hours at 245F. Nomenclature:  
Desiccant, activated, bagged, packaging use and static  
dehumidification. Manufacturer's name and address.

Shelf Life: Two years from date of MFG

Packaging: The product will be packaged as follows:  
- 5 desiccant packets per strip / 60 strips per 30 gallon fiber drum with  
a Humidity Indicating Card  
- Product label contains the following:  
Manufacturer's name  
Description of desiccant  
Quantity per container  
Manufacturer's part number  
Manufacturer's lot number  
Date of manufacture

RECEIVED  
APR 30 2020

Revision Date: 03/20/13  
ISO 9001:2008