


UN Certification: 1H2/Y19/S &amp; 1H2/X16/S


**Cover Specification | Part No. 600137**

Origin of Materials	100% US Manufactured
Material	High Density Polyethylene (HDPE)
Mfg. Method	Injection Molding
Weight	1.03 lbs / 470g
Height	2.50"
Min Wall Thickness	.090"
Bottom OD	12.69"
Gasket	Neoprene / Available w/o Gasket
Standard Color	White, Red 186, Blue 300 <i>Please inquire about custom colors</i>
Case Quantity	20 per case / 240 per pallet


**Pail Specification | Part No. 500041**

Origin of Materials	100% US Manufactured
Material	High Density Polyethylene (HDPE)
Mfg. Method	Injection Molding / Plastic Handle
Weight	1.84 lbs / 836 g
Min. Wall Thickness	.090"
Exterior Height	11.50" With Cover 12.19"
Top Inside Diameter	11.25"
Bottom Outside Diameter	10.27"
Shipping Quantity	6 Pails and Covers Per Case 144 Units per Pallet 4,752 Units per Truckload
Standard Color	White <i>Please inquire about custom colors</i>
Print Area	6.25" x 33"
Heat Transfer Area	6.63" x 31.75"
Top Load Rating	375 lbs


 RECEIVED  
 APR 28 2020
