

HDPE DRUM EU COMPLIANCE

BioSpectra has evaluated the packaging component supply through the Supplier Approval Program. The drums utilized by BioSpectra for material packaging are manufactured from virgin high-density molecular weight polyethylene resin in compliance with 21 CFR 177.1520. Correlation and harmonization of EU Framework Regulation (EC) 1935/2004, GMP Regulation EC 2023/2006, and Plastics Regulation EU/10/2011 (as amended by EU/2017/725 EU/2020/1245) with US FDA regulation is generally accepted in the United States.

These drums additionally comply with the following regulations, and are UN rated as detailed in the table below.

- EU RoHS
- EU Packaging Directive 94/62

Part Number	Description	UN Rating
COMP-0070	12 Gallon Blue Drum	1H2/X65/S
COMP-0071	20 Gallon Blue Drum	1H2/Y130/S
COMP-0193	30 Gallon Light Blue Drum with Lid and Plastic Lever Lock	Not UN Rated
COMP-0194	30 Gallon Dark Blue Drum with Lid and Plastic Lever Lock	1H2/Y80/S
COMP-0195	30 Gallon Dark Blue Drum with Lid and Plastic Lever Lock	1H2/Y120/S with liner 1H2/Y80/S without liner
COMP-0090	55 Gallon Black Drum	1H1/Y.1.9/150
COMP-0091	55 Gallon Natural Drum	1H1/Y.1.9/150

For further information, please contact info@biospectra.us

Cassie Baum
Cassie Baum
Senior Compliance Specialist