

24" x 36" POLY LINER SPECIFICATIONS

BSI COMPONENT NUMBER: COMP-0014



24" x 36" Poly Liner	
Test	Specifications
Appearance and Color	Clear with no punctures
Dimensions (Nominal)	24" x 36"
Thickness	0.0027" minimum
Materials of Construction (Film)	Hexene LLDPE Nonslip
FDA Statement	21 CFR 177.1520
UATR	Passes Test; when compared to a previously approved poly liner for this part number

For further information, please contact marketing@biospectra.us.

**Connie Augustine | Director
Business Integration & Retention
BioSpectra Inc.**

This is not considered a controlled document. We are not responsible for any errors or omissions, and the user is responsible for any decisions based on the information herein.