

TRIS / TROMETHAMINE

MEETS USP, EP, ACS, LBLE

CAS #: 77-86-1

Formula: $\text{NH}_2\text{C}(\text{CH}_2\text{OH})_3$

F.W.: 121.14 g/mol

TRIS-6250

BIO QUALIFIED GRADE

ANALYSIS		SPECIFICATIONS
Absorbance (10%)	260 nm	< = 0.03 a.u.
	280 nm	< = 0.02 a.u.
	430 nm	< = 0.004 a.u.
Absorbance (40%)	280 nm	< = 0.06 a.u.
	290 nm	< = 0.20 a.u.
	400 nm	< = 0.02 a.u.
	430 nm	< = 0.03 a.u.
APHA Color, 20% Aqueous Solution		< = 20 APHA
Appearance (EP)		Clear and Colorless
Assay, Dried Basis (USP/EP)		> = 99.9 - 100.1%
Chloride (EP)		< = 100 ppm
Elemental Impurities		< = 2000
Endotoxin		< = 2.5 EU/g
Heavy Metals		< = 5 ppm
Identification A (EP, pH 5%) (USP)		Passes Test
Identification B (USP), Melting Range (EP)		168 - 172°C
Identification C (USP/EP)		Passes Test
Insoluble Matter		< = 0.005%
Loss on Drying (EP)		< = 0.5%
Related Substances		< = 1.0%
Residual Solvents (USP)		Passes Test
Residue on Ignition (USP) / Sulfated Ash (EP)		< = 0.1%
Sulfate		< = 5 ppm
Water		< = 0.2%
Microbial Content	TAMC	< = 10 CFU/g
	TYMC	< = 10 CFU/g

ANALYSIS		SPECIFICATIONS
Trace Metals	Aluminum (Al)	< = 5 ppm
	Arsenic (As)	< = 1 ppm
	Barium (Ba)	< = 5 ppm
	Cadmium (Cd)	< = 5 ppm
	Calcium (Ca)	< = 10 ppm
	Chloride (Cl)	< = 2 ppm
	Chromium (Cr)	< = 5 ppm
	Cobalt (Co)	< = 5 ppm
	Copper (Cu)	< = 5 ppm
	Iron (Fe)	< = 2 ppm
	Lead (Pb)	< = 1 ppm
	Lithium (Li)	< = 5 ppm
	Magnesium (Mg)	< = 5 ppm
	Manganese (Mn)	< = 5 ppm
	Nickel (Ni)	< = 1 ppm
	Potassium (K)	< = 50 ppm
	Sodium (Na)	< = 50 ppm
	Strontium (Sr)	< = 5 ppm
	Zinc (Zn)	< = 5 ppm

[Click here to view SDS, CoAs and other supporting regulatory documents on our website.](#)

Refer to page 2 for additional product details.

This is not considered a controlled document. We are not responsible for any errors or omissions, and the user is responsible for any decisions based on the information herein.

General Product Overview

Tris/Tromethamine is commonly used as a buffer and excipient in downstream biological applications during the manufacture of drug products. Tris is considered a Good's buffer because it has low UV absorptivity, minimal reactivity, stable pH and is soluble in water.

Industry Application

Suitable for use as a cGMP chemical in pharmaceutical manufacturing processes and biological/biotech process applications from R&D through scale cGMP production.

Key Product Features

- Appears as a white crystalline powder
- Manufactured in accordance with the ISO Certified Quality Management System
- Manufactured in an enzyme free, hormone free and animal free environment
- Contains no known major food allergens (as defined by FDA and WHO)
- The final product and its raw materials are not derived from nor come into contact with animal parts, animal products, and/or animal byproducts or derivatives.
- Is not subject to genetic modification
- Synonyms: Trometamol; THAM; Tris (Hydroxyethyl) Aminomethane; Trisaminol; Trismethylolaminomethane; 2-Amino-2 (Hydroxymethyl) propane-1,2-diol

Storage and Shipping Conditions

Refer to SDS.

Standard Shelf-Life Policy

Please inquire for more information regarding shelf-life.

Package Sizes

1kg, 5kg, 10kg, 25kg, 50kg

This is not considered a controlled document. We are not responsible for any errors or omissions, and the user is responsible for any decisions based on the information herein.