

9/17/2026

Version 1.3

TYVEK LINER SPECIFICATIONS



Tyvek Liner Component Specifications	
Test	Specifications
Appearance and Color	White with no punctures
Width	36.0" $\pm$ 0.5"
Length (Internal Tolerance)	60.0" $\pm$ 1.0"
Weight	219.2 grams

For further information, please contact marketing@biospectra.us.



Connie Augustine | Director of Commercial Operations
BioSpectra Inc.

This is not considered a controlled document. We are not responsible for any errors or omissions, and the user is responsible for any decisions based on the information herein.