

12/19/2025

Version 1.1

275 GAL / 1,000L TOTE SPECIFICATIONS

BSI COMPONENT NUMBER: COMP-0145



275 Gallon / 1,000L Tote	
Test	Specifications
Appearance and Color	Natural tote with steel cage and red fill opening screw cap; clean and free of damage
Outside Width (Nominal)	39.37"
Outside Length (Nominal)	47.24"
Outside Height (Nominal)	45.67"
Fill Opening Screw Cap (Nominal)	6"
Discharge Opening (Nominal)	2"
Capacity	275 Gallons
Pallet Type	Plastic Skid Pallet
UN Rating	31HA1/Y

For further information, please contact marketing@biospectra.us.

**Connie Augustine | Director
Business Integration & Retention
BioSpectra Inc.**

This is not considered a controlled document. We are not responsible for any errors or omissions, and the user is responsible for any decisions based on the information herein.