

12/19/2025

Version 1.1

0.6 GAL / 1KG PAIL SPECIFICATIONS

BSI COMPONENT NUMBER: COMP-0111



0.6 Gallon / 1kg HDPE Round Screw Top Pail and Lid	
Test	Specifications
Appearance and Color	White HDPE round pail with hole(s) for tamper evident seal, blue tamper evident lid with gasket, red trigger, and slot on top of lid
Inside Height	6.75" \pm 0.5"
Inside Top Diameter	6.15" \pm 0.5"
Outside Height, with Lid	7.0" \pm 0.5"
Outside Bottom Diameter	5.60" \pm 0.5"
UN Rating	1H2/Y4/S/YR

For further information, please contact marketing@biospectra.us.



**Connie Augustine | Director
Business Integration & Retention
BioSpectra Inc.**

This is not considered a controlled document. We are not responsible for any errors or omissions, and the user is responsible for any decisions based on the information herein.