

ANTIBLOCK ADDITIVE STATEMENT

Tris GMP

BioSpectra has assessed the supply of Poly Liners from BioSpectra's approved supplier. Based on this assessment, BioSpectra can state that Tris, Bio Excipient Grade utilizes poly liners with antiblock additives talc (magnesium silicate hydrate) and diatomaceous earth (flux calcinated silica) used in the production of the non-slip liners. All poly liners supplied to BioSpectra and used to package Tris, Bio Excipient Grade are manufactured from virgin polyethylene FDA approved material conforming to 21CFR1520 regulation.

Current Product Number	Historic Product Number
TRIS-3254	TR3254
TRIS-3255	TR3255

For further information, please contact info@biospectra.us



Cassie Baun
Compliance Specialist