

	A	B	C	D	E	F	G	H	I	J	K
1	Product	Description 1	Product Status	Product size	Product UOM	Package size	Package UOM	Package type	Liner size	Liner type	Desiccant
2	BTRI-3250-01	Bis Tris, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
3	BTRI-3250-05	Bis Tris, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
4	BTRI-3250-10	Bis Tris, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
5	BTRI-3250-25	Bis Tris, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
6	BTRI-3250-92	Bis Tris, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
7	BTRI-3250-93	Bis Tris, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
8	BTRI-4250-01	Bis-Tris Biotech, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
9	BTRI-4250-05	Bis-Tris Biotech, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
10	BTRI-4250-10	Bis-Tris Biotech, 10 KG	In development	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
11	BTRI-4250-25	Bis-Tris Biotech, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
12	BTRI-4250-50	Bis-Tris Biotech, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
13	BTRI-4250-92	Bis-Tris Biotech, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
14	BTRI-4250-93	Bis-Tris Biotech, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
15	CSMH-3250-01	Cysteamine HCl (2-MEA), 1 KG	In development	1.00 KG		1.25 GAL		White HDPE round pail, blue lid	(2x) 24 x 36	LDPE	4 Unit
16	CSMH-3250-05	Cysteamine HCl (2-MEA), 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	(2x) 24 x 36	LDPE	4 Unit
17	CSMH-3250-10	Cysteamine HCl (2-MEA), 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	(2x) 24 x 36	LDPE	4 Unit
18	CSMH-3250-25	Cysteamine HCl (2-MEA), 25 KG	In development	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	(2x) 30 x 48	LDPE	
19	CSMH-3250-90	Cysteamine HCl (2-MEA), 100 G	Active	100.00 G		3.50 GAL		White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	4 Unit
20	CSMH-3250-91	Cysteamine HCl (2-MEA), 50 G	In development	50.00 G		0.60 GAL		White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	4 Unit
21	CSMH-3250-92	Cysteamine HCl (2-MEA), 100 G	Active	100.00 G		0.60 GAL		White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	4 Unit
22	CSMH-3250-94	Cysteamine HCl (2-MEA), 1 KG	In development	1.00 KG		0.60 GAL		White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	4 Unit
23	CSMH-3350-10	Cysteamine HCl (2-MEA), 10 KG	In development	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	(2x) 24 x 36	LDPE	4 Unit
24	CSMH-3350-25	Cysteamine HCl (2-MEA), 25 KG	In development	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	(2x) 30 x 48	LDPE	4 Unit
25	D010-3201-01	Dextran Powder MW 10,000, 1 KG	Active	1.00 KG		0.60 GAL		White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	4 Unit
26	D010-3201-30	Dextran Powder MW 10,000, 30 KG	Active	30.00 KG		30.00 GAL		Dark blue HDPE drum	(2x) 34 x 60	LDPE	4 Unit
27	DXSE-4250-01	Dextran Sulfate, 8000, GMP, 1kg	In development	1.00 KG		0.60 GAL		White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	4 Unit
28	DXSE-4250-05	Dextran Sulfate, 8000, GMP, 5kg	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	(2x) 24 x 36	LDPE	4 Unit
29	DXSE-4250-10	Dextran Sulfate, 8000, GMP, 10kg	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	(2x) 24 x 36	LDPE	4 Unit
30	DXSE-4250-15	Dextran Sulfate, 8000, GMP, 15kg	Active	15.00 KG		12.00 GAL		White HDPE round pail, blue lid	(2x) 30 x 48	LDPE	4 Unit
31	DXSE-4250-25	Dextran Sulfate, 8000, GMP, 25kg	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	(2x) 30 x 48	LDPE	4 Unit
32	DXSE-4250-92	Dextran Sulfate, 8000, GMP, 100g	Active	100.00 G		0.60 GAL		White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	4 Unit
33	DXSE-4250-93	Dextran Sulfate, 8000, GMP, 20g	Active	20.00 G		0.60 GAL		White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	4 Unit
34	DXSE-4250-94	Dextran Sulfate, 8000, GMP, 500g	Active	500.00 G		0.60 GAL		White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	4 Unit
35	DXSE-4252-01	Dextran Sulfate Sodium Salt MW 8000 Biotech, 1 KG	Active	1.00 KG		0.60 GAL		White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	4 Unit
36	DXSE-4252-05	Dextran Sulfate Sodium Salt MW 8000 Biotech, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	(2x) 24 x 36	LDPE	4 Unit
37	DXSE-4252-10	Dextran Sulfate Sodium Salt MW 8000 Biotech, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	(2x) 24 x 36	LDPE	4 Unit
38	DXSE-4252-25	Dextran Sulfate Sodium Salt MW 8000 Biotech, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	(2x) 30 x 48	LDPE	4 Unit
39	DXSE-4252-50	Dextran Sulfate Sodium Salt MW 8000 Biotech, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	(2x) 34 x 60	LDPE	4 Unit
40	DXSE-4252-92	Dextran Sulfate Sodium Salt MW 8000 Biotech, 100 G	Active	100.00 G		0.60 GAL		White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	4 Unit
41	DXSE-4252-93	Dextran Sulfate Sodium Salt MW 8000 Biotech, 500 G	Active	500.00 G		0.60 GAL		White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	4 Unit
42	DXSE-5201-92	Dextran Sulfate, 8000, 100g	Active	100.00 G		0.60 GAL		White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	4 Unit
43	GALP-3224-06	D-Galactose, Plant Derived, 1KG	Active	1.00 KG		0.60 GAL		White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	
44	GALP-3224-25	D-Galactose, Plant Derived, 20kg	In development	20.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
45	GALP-3250-05	D-Galactose, Plant Derived, 5kg	In development	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	

	A	B	C	D	E	F	G	H	I	J	K
1	Product	Description 1	Product Status	Product size	Product UOM	Package size	Package UOM	Package type	Liner size	Liner type	Desiccant
46	GALP-3250-06	D-Galactose, Plant Derived, 1KG	Active	1.00 KG		0.60 GAL		White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	
47	GALP-3250-10	D-Galactose, Plant Derived, 10kg	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
48	GALP-3250-25	D-Galactose, Plant Derived, 20kg	Active	20.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
49	GALP-3250-27	D-Galactose, Plant Derived, 25kg	Active	25.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
50	GALP-3251-06	D-Galactose, Plant Derived, 1KG	Active	1.00 KG		0.60 GAL		White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	
51	GALP-3251-10	D-Galactose, Plant Derived, 10kg	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
52	GALP-3251-21	D-Galactose, Plant Derived, 20kg	Active	20.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
53	GALP-3251-27	D-Galactose, Plant Derived, 25kg	Active	25.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
54	GALP-3251-50	D-Galactose, Plant Derived, 50kg	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
55	GALP-3252-06	D-Galactose, Plant Derived, 1KG	Active	1.00 KG		0.60 GAL		White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	
56	GALP-3252-27	D-Galactose, Plant Derived, 25kg	Active	25.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
57	GALP-3350-06	D-Galactose, Plant Derived, 1KG	Active	1.00 KG		0.60 GAL		White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	
58	GALP-3350-25	D-Galactose, Plant Derived, 20kg	Active	20.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
59	GALP-3350-27	D-Galactose, Plant Derived, 25kg	Active	25.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
60	GALP-3351-06	D-Galactose, Plant Derived, 1KG	Active	1.00 KG		0.60 GAL		White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	
61	GALP-3351-27	D-Galactose, Plant Derived, 25kg	Active	25.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
62	GALP-3351-50	D-Galactose, Plant Derived, 50kg	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
63	GALP-4250-05	Galactose (D) Biotech, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
64	GALP-4250-06	Galactose (D) Biotech, 1 KG	Active	1.00 KG		0.60 GAL		White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	
65	GALP-4250-10	Galactose (D) Biotech, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
66	GALP-4250-21	Galactose (D) Biotech, 20 KG	Active	20.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
67	GALP-4250-27	Galactose (D) Biotech, 25 KG	Active	25.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
68	GALP-4250-50	Galactose (D) Biotech, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
69	GALP-4250-92	Galactose (D) Biotech, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
70	GHCL-3101-02	Guanidine HCl 6M Solution, 1L	Active	1.00 L		1.00 L		Natural HDPE bottle			
71	GHCL-3101-72	Guanidine HCl 6M Solution, 200 L	Active	200.00 L		55.00 GAL		Black HDPE drum			
72	GHCL-3101-78	Guanidine HCl 6M Solution, 1135 L	Active	1135.00 L		330.00 GAL		Tote			
73	GHCL-3101-80	Guanidine HCl 6M Solution, 200 L	Active	200.00 L		55.00 GAL		Black HDPE drum			
74	GHCL-3201-59	Guanidine HCl, 50 KG	Active	50.00 KG		26.00 GAL		Fiber drum	(2x) 34 x 60	LDPE	4 Unit
75	GHCL-3202-01	Guanidine HCl, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
76	GHCL-3202-50	Guanidine HCl, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	4 Unit
77	GHCL-3202-93	Guanidine HCl, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
78	GHCL-3220-01	Guanidine HCl, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
79	GHCL-3220-10	Guanidine HCl, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
80	GHCL-3220-25	Guanidine HCl, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
81	GHCL-3220-50	Guanidine HCl, 50 KG	In development	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
82	GHCL-3220-51	Guanidine HCl, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	4 Unit
83	GHCL-3220-92	Guanidine HCl, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
84	GHCL-3220-93	Guanidine HCl, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
85	GHCL-3220-94	Guanidine HCl, 2.5 KG	Active	2.50 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
86	GHCL-3221-25	Guanidine HCl, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
87	GHCL-3224-25	Guanidine HCl, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
88	GHCL-3224-93	Guanidine HCl, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
89	GHCL-3225-01	Guanidine HCl, 1KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			

	A	B	C	D	E	F	G	H	I	J	K
1	Product	Description 1	Product Status	Product size	Product UOM	Package size	Package UOM	Package type	Liner size	Liner type	Desiccant
90	GHCL-3225-25	Guanidine HCl, 25KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
91	GHCL-3225-92	Guanidine HCl, 100G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
92	GHCL-3250-01	Guanidine HCl, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
93	GHCL-3251-59	Guanidine HCl, 50 KG	Active	50.00 KG		26.00 GAL		Fiber drum	(2x) 34 x 60	LDPE	4 Unit
94	GHCL-3251-93	Guanidine HCl, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
95	GHCL-3253-25	Guanidine HCl, 25KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
96	GHCL-3253-93	Guanidine HCl, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
97	GHCL-3255-25	Guanidine HCl, 25KG	In development	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
98	GHCL-4223-01	Guanidine Hydrochloride Biotech, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
99	GHCL-4223-05	Guanidine Hydrochloride Biotech, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
100	GHCL-4223-10	Guanidine Hydrochloride Biotech, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
101	GHCL-4223-25	Guanidine Hydrochloride Biotech, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
102	GHCL-4223-50	Guanidine Hydrochloride Biotech, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
103	GHCL-4223-92	Guanidine Hydrochloride Biotech, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
104	GHCL-4223-93	Guanidine Hydrochloride Biotech, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
105	GHCL-5220-01	Guanidine HCl, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
106	GHCL-5220-25	Guanidine HCl, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
107	GTHI-3220-10	Guanidine Thiocyanate, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
108	GTHI-3220-14-76	Guanidine Thiocyanate, 14.76 KG	Active	14.76 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE/ 25L EZI-	
109	GTHI-3220-15	Guanidine Thiocyanate, 14.76 KG	In development	14.76 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
110	GTHI-3220-25	Guanidine Thiocyanate, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
111	GTHI-3220-92	Guanidine Thiocyanate, 100 G	Active	100.00 G		250.00 ML					
112	GTHI-3220-95	Guanidine Thiocyanate, 10 G	In development	10.00 G		60.00 ML					
113	GTHI-4220-05	Guanidine Thiocyanate, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
114	GTHI-4220-10	Guanidine Thiocyanate, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
115	GTHI-4220-25	Guanidine Thiocyanate, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
116	GTHI-4221-01	Guanidine Thiocyanate Biotech, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
117	GTHI-4221-05	Guanidine Thiocyanate Biotech, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
118	GTHI-4221-10	Guanidine Thiocyanate Biotech, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
119	GTHI-4221-25	Guanidine Thiocyanate Biotech, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
120	GTHI-4221-50	Guanidine Thiocyanate Biotech, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
121	GTHI-4221-92	Guanidine Thiocyanate Biotech, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
122	GTHI-4221-93	Guanidine Thiocyanate Biotech, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
123	GTHI-5203-10	Guanidine Thiocyanate, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
124	GTHI-5220-10	Guanidine Thiocyanate, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
125	GTHI-5221-05	Guanidine Thiocyanate, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
126	GTHI-5221-10	Guanidine Thiocyanate, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
127	GTHI-5221-25	Guanidine Thiocyanate, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
128	GTHI-5221-92	Guanidine Thiocyanate, 100 G	Active	100.00 G		250.00 ML					
129	HEHH-5201-25	HEPES Hemi Sodium, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
130	HEPE-3220-01	HEPES, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
131	HEPE-3220-06	HEPES, 1 KG	Active	1.00 KG		0.60 GAL		White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	
132	HEPE-3220-10	HEPES, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
133	HEPE-3220-20	HEPES, 20 KG	Active	20.00 KG		12.00 GAL		White HDPE round pail, blue lid	(2x) 30 x 48	LDPE	

	A	B	C	D	E	F	G	H	I	J	K
1	Product	Description 1	Product Status	Product size	Product UOM	Package size	Package UOM	Package type	Liner size	Liner type	Desiccant
134	HEPE-3220-22	HEPES, 2 KG	Active	2.00 KG		1.25 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
135	HEPE-3220-25	HEPES, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
136	HEPE-3220-50	HEPES, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
137	HEPE-3220-91	HEPES, 250 G	Active	250.00 G		1000.00 ML		White HDPE bottle, blue lid			
138	HEPE-3220-92	HEPES, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
139	HEPE-3220-93	HEPES, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
140	HEPE-3221-10	HEPES, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
141	HEPE-3250-05	HEPES, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
142	HEPE-3250-10	HEPES, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
143	HEPE-3251-01	HEPES, 1kg	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
144	HEPE-3251-25	HEPES, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
145	HEPE-3251-50	HEPES, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
146	HEPE-3251-92	HEPES, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
147	HEPE-3254-25	HEPES, 25 KG	In development	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
148	HEPE-3254-50	HEPES, 50 KG	In development	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
149	HEPE-3320-25	HEPES, 25 KG	In development	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
150	HEPE-3320-92	HEPES, 100 G	In development	100.00 G		250.00 ML		White HDPE bottle, blue lid			
151	HEPE-3351-25	HEPES, 25 KG	In development	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
152	HEPE-3351-92	HEPES, 100 G	In development	100.00 G		250.00 ML		White HDPE bottle, blue lid			
153	HEPE-3354-25	HEPES, 25 KG	In development	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
154	HEPE-3354-50	HEPES, 50 KG	In development	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
155	HEPE-4250-01	HEPES Biotech, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
156	HEPE-4250-05	HEPES Biotech, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
157	HEPE-4250-10	HEPES Biotech, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
158	HEPE-4250-25	HEPES Biotech, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
159	HEPE-4250-50	HEPES Biotech, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
160	HEPE-4250-92	HEPES Biotech, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
161	HEPE-4250-93	HEPES Biotech, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
162	HEPE-5201-26	HEPES, 25 KG	Active	25.00 KG		26.00 GAL		Fiber drum	(2x) 34 x 60	LDPE	
163	HEPE-5203-40	HEPES, 40 KG	Active	40.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
164	HEPE-5220-01	HEPES, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
165	HEPE-5220-10	HEPES, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
166	HEPN-5201-26	HEPES Sodium, 25 KG	Active	25.00 KG		26.00 GAL		Fiber drum	(2x) 34 x 60	LDPE	
167	HEPN-5203-50	HEPES Sodium, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
168	HEPN-5220-25	HEPES Sodium, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
169	HEPN-5251-25	HEPES Sodium, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
170	IHCL-4101-02	6N HCl in IPA Solution, 500ML	Active	500.00 ML		17.00 OZ		Black HDPE UN leakproof bottle			
171	IHCL-4101-97	6N HCl in IPA Solution, 180 KG	Active	180.00 KG		55.00 GAL		Black HDPE drum			
172	KBRO-2220-01	Potassium Bromide, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
173	KBRO-2220-05	Potassium Bromide, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	(2x) 24 x 36	LDPE	4 Unit
174	KBRO-2220-10	Potassium Bromide, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	(2x) 24 x 36	LDPE	4 Unit
175	KBRO-2220-25	Potassium Bromide, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	(2x) 30 x 48	LDPE	4 Unit
176	KBRO-2220-92	Potassium Bromide, 100 G	Active	100.00 G		0.60 GAL		White HDPE round pail, blue lid	9 x 15	LDPE	
177	KBRO-2301-01	Potassium Bromide, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			

	A	B	C	D	E	F	G	H	I	J	K
1	Product	Description 1	Product Status	Product size	Product UOM	Package size	Package UOM	Package type	Liner size	Liner type	Desiccant
178	KBRO-2301-25	Potassium Bromide, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	(2x) 30 x 48	LDPE	
179	KBRO-2302-01	Potassium Bromide, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
180	KBRO-2302-25	Potassium Bromide, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	(2x) 30 x 48	LDPE	
181	LARH-6250-01	Arginine (L) Hydrochloride Biotech, 1 KG	In development	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
182	LARH-6250-05	Arginine (L) Hydrochloride Biotech, 5 KG	In development	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
183	LARH-6250-10	Arginine (L) Hydrochloride Biotech, 10 KG	In development	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
184	LARH-6250-25	Arginine (L) Hydrochloride Biotech, 25 KG	In development	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
185	LARH-6250-50	Arginine (L) Hydrochloride Biotech, 50 KG	In development	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
186	LARH-6250-92	Arginine (L) Hydrochloride Biotech, 100 G	In development	100.00 G		250.00 ML		White HDPE bottle, blue lid			
187	LARH-6250-93	Arginine (L) Hydrochloride Biotech, 500 G	In development	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
188	LARH-6350-01	Arginine (L) Hydrochloride Biotech, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
189	LARH-6350-05	Arginine (L) Hydrochloride Biotech, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
190	LARH-6350-10	Arginine (L) Hydrochloride Biotech, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
191	LARH-6350-25	Arginine (L) Hydrochloride Biotech, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
192	LARH-6350-50	Arginine (L) Hydrochloride Biotech, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
193	LARH-6350-92	Arginine (L) Hydrochloride Biotech, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
194	LARH-6350-93	Arginine (L) Hydrochloride Biotech, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
195	LCYS-6250-01	Cystine (L) Hydrochloride Biotech, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
196	LCYS-6250-05	Cystine (L) Hydrochloride Biotech, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
197	LCYS-6250-10	Cystine (L) Hydrochloride Biotech, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
198	LCYS-6250-25	Cystine (L) Hydrochloride Biotech, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
199	LCYS-6250-50	Cystine (L) Hydrochloride Biotech, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
200	LCYS-6250-92	Cystine (L) Hydrochloride Biotech, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
201	LCYS-6250-93	Cystine (L) Hydrochloride Biotech, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
202	LCYS-6350-01	Cystine (L) Hydrochloride Biotech, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
203	LCYS-6350-05	Cystine (L) Hydrochloride Biotech, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
204	LCYS-6350-10	Cystine (L) Hydrochloride Biotech, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
205	LCYS-6350-25	Cystine (L) Hydrochloride Biotech, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
206	LCYS-6350-50	Cystine (L) Hydrochloride Biotech, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
207	LCYS-6350-92	Cystine (L) Hydrochloride Biotech, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
208	LCYS-6350-93	Cystine (L) Hydrochloride Biotech, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
209	LGLM-6201-01	Glutamine (L) Biotech, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
210	LGLM-6201-05	Glutamine (L) Biotech, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
211	LGLM-6201-10	Glutamine (L) Biotech, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
212	LGLM-6201-25	Glutamine (L) Biotech, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
213	LGLM-6201-50	Glutamine (L) Biotech, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
214	LGLM-6201-92	Glutamine (L) Biotech, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
215	LGLM-6201-93	Glutamine (L) Biotech, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
216	LGLM-6301-01	Glutamine (L) Biotech, 1 KG	In development	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
217	LGLM-6301-05	Glutamine (L) Biotech, 5 KG	In development	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
218	LGLM-6301-10	Glutamine (L) Biotech, 10 KG	In development	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
219	LGLM-6301-25	Glutamine (L) Biotech, 25 KG	In development	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
220	LGLM-6301-50	Glutamine (L) Biotech, 50 KG	In development	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
221	LGLM-6301-92	Glutamine (L) Biotech, 100 G	In development	100.00 G		250.00 ML		White HDPE bottle, blue lid			

	A	B	C	D	E	F	G	H	I	J	K
1	Product	Description 1	Product Status	Product size	Product UOM	Package size	Package UOM	Package type	Liner size	Liner type	Desiccant
222	LGLM-6301-93	Glutamine (L) Biotech, 500 G	In development	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
223	LHMM-6250-01	Histidine (L) MonoHCl Monohydrate Biotech, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
224	LHMM-6250-05	Histidine (L) MonoHCl Monohydrate Biotech, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
225	LHMM-6250-10	Histidine (L) MonoHCl Monohydrate Biotech, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
226	LHMM-6250-25	Histidine (L) MonoHCl Monohydrate Biotech, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
227	LHMM-6250-50	Histidine (L) MonoHCl Monohydrate Biotech, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
228	LHMM-6250-92	Histidine (L) MonoHCl Monohydrate Biotech, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
229	LHMM-6250-93	Histidine (L) MonoHCl Monohydrate Biotech, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
230	LHMM-6350-01	Histidine (L) MonoHCl Monohydrate Biotech, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
231	LHMM-6350-05	Histidine (L) MonoHCl Monohydrate Biotech, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
232	LHMM-6350-10	Histidine (L) MonoHCl Monohydrate Biotech, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
233	LHMM-6350-25	Histidine (L) MonoHCl Monohydrate Biotech, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
234	LHMM-6350-50	Histidine (L) MonoHCl Monohydrate Biotech, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
235	LHMM-6350-92	Histidine (L) MonoHCl Monohydrate Biotech, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
236	LHMM-6350-93	Histidine (L) MonoHCl Monohydrate Biotech, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
237	MESH-3201-50	MES Hydrate, 50KG	In development	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
238	MESH-3201-92	MES Hydrate, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
239	MESH-3250-01	MES Hydrate, 1KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
240	MESH-3250-25	MES Hydrate, 25KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
241	MESH-3350-01	MES Hydrate, 1KG	In development	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
242	MESH-4250-01	MES Hydrate Biotech, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
243	MESH-4250-05	MES Hydrate Biotech, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
244	MESH-4250-10	MES Hydrate Biotech, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
245	MESH-4250-25	MES Hydrate Biotech, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
246	MESH-4250-50	MES Hydrate Biotech, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
247	MESH-4250-92	MES Hydrate Biotech, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
248	MESH-4250-93	MES Hydrate Biotech, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
249	MESH-9250-91	MES Hydrate, 50 G	In development	50.00 G		250.00 ML		White HDPE bottle, blue lid			
250	MESM-3220-01	MES Monohydrate, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
251	MESM-3220-05	MES Monohydrate, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
252	MESM-3220-10	MES Monohydrate, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
253	MESM-3220-12	MES, Monohydrate, 12 KG	Active	12.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
254	MESM-3220-25	MES, Monohydrate, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
255	MESM-3220-50	MES, Monohydrate, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
256	MESM-3220-91	MES Monohydrate, 50G	Active	50.00 G		250.00 ML		White HDPE bottle, blue lid			
257	MESM-3220-92	MES, Monohydrate, 100G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
258	MESM-3220-93	MES, Monohydrate, 500G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
259	MESM-3221-10	MES Monohydrate, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
260	MESM-3221-93	MES, Monohydrate, 500G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
261	MESM-3222-01	MES Monohydrate, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
262	MESM-3222-03	MES Monohydrate, 3 KG	Active	3.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36		
263	MESM-3222-25	MES, Monohydrate, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
264	MESM-3222-50	MES, Monohydrate, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
265	MESM-3222-91	MES Monohydrate, 50G	Active	50.00 G		250.00 ML		White HDPE bottle, blue lid			

	A	B	C	D	E	F	G	H	I	J	K
1	Product	Description 1	Product Status	Product size	Product UOM	Package size	Package UOM	Package type	Liner size	Liner type	Desiccant
266	MESM-3222-92	MES, Monohydrate, 100G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
267	MESM-3222-93	MES, Monohydrate, 500G	In development	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
268	MESM-3223-05	MES, Monohydrate, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
269	MESM-3250-06	MES, Monohydrate, 1 KG	Active	1.00 KG		0.60 GAL		White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	
270	MESM-3250-10	MES Monohydrate, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
271	MESM-3250-12	MES, Monohydrate, 12 KG	Active	12.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
272	MESM-3250-25	MES, Monohydrate, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
273	MESM-3250-50	MES, Monohydrate, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
274	MESM-3250-91	MES, Monohydrate, 250G	Active	250.00 G		1000.00 ML		White HDPE bottle, blue lid			
275	MESM-3250-92	MES, Monohydrate, 100G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
276	MESM-3250-93	MES, Monohydrate, 500G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
277	MESM-3250-95	MES, Monohydrate, 10G	In development	10.00 G		250.00 ML		White HDPE bottle, blue lid			
278	MESM-3251-25	MES, Monohydrate, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
279	MESM-3251-92	MES, Monohydrate, 100G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
280	MESM-3320-25	MES, Monohydrate, 25 KG	In development	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
281	MESM-3320-92	MES, Monohydrate, 100G	In development	100.00 G		250.00 ML		White HDPE bottle, blue lid			
282	MESM-3321-10	MES Monohydrate, 10 KG	In development	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
283	MESM-3321-93	MES, Monohydrate, 500G	In development	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
284	MESM-3350-25	MES, Monohydrate, 25 KG	In development	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
285	MESM-3350-92	MES, Monohydrate, 100G	In development	100.00 G		250.00 ML		White HDPE bottle, blue lid			
286	MESM-4220-01	MES, Monohydrate, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
287	MESM-4220-10	MES, Monohydrate, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
288	MESM-4220-25	MES, Monohydrate, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
289	MESM-4220-26	MES, Monohydrate, 25 KG	Active	25.00 KG		26.00 GAL		Fiber drum	(2x) 34 x 60	LDPE	
290	MESM-4221-01	MES, Monohydrate Biotech, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
291	MESM-4221-05	MES, Monohydrate Biotech, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
292	MESM-4221-10	MES, Monohydrate Biotech, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
293	MESM-4221-25	MES, Monohydrate Biotech, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
294	MESM-4221-50	MES, Monohydrate Biotech, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
295	MESM-4221-92	MES, Monohydrate Biotech, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
296	MESM-4221-93	MES, Monohydrate Biotech, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
297	MESM-5201-26	MES, Monohydrate, 25 KG	Active	25.00 KG		26.00 GAL		Fiber drum	(2x) 34 x 60	LDPE	
298	MESM-5202-10	MES, Monohydrate, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
299	MESM-5220-10	MES, Monohydrate, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
300	MESM-5220-92	MES, Monohydrate, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
301	MESM-5220-93	MES, Monohydrate, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
302	MESN-5220-10	MES Sodium, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
303	MESN-5220-26	MES Sodium, 25 KG	Active	25.00 KG		26.00 GAL		Fiber drum	(2x) 34 x 60	LDPE	
304	MESN-5220-50	MES Sodium, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
305	MESN-5220-93	MES Sodium, 300 G	Active	300.00 G		1000.00 ML		White HDPE bottle, blue lid			
306	MOPS-3201-10	MOPS, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
307	MOPS-3201-51	MOPS, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
308	MOPS-3220-01	MOPS, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
309	MOPS-3220-25	MOPS, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	

	A	B	C	D	E	F	G	H	I	J	K
1	Product	Description 1	Product Status	Product size	Product UOM	Package size	Package UOM	Package type	Liner size	Liner type	Desiccant
310	MOPS-3220-50	MOPS, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
311	MOPS-3220-51	MOPS, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
312	MOPS-3220-92	MOPS, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
313	MOPS-3220-93	MOPS, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
314	MOPS-4220-05	MOPS, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
315	MOPS-4221-01	MOPS Biotech, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
316	MOPS-4221-05	MOPS Biotech, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
317	MOPS-4221-10	MOPS Biotech, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
318	MOPS-4221-25	MOPS Biotech, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
319	MOPS-4221-50	MOPS Biotech, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
320	MOPS-4221-92	MOPS Biotech, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
321	MOPS-4221-93	MOPS Biotech, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
322	MOPS-5201-26	MOPS, 50 KG	Active	50.00 KG		26.00 GAL		Fiber drum	(2x) 34 x 60	LDPE	
323	MOPS-5202-25	MOPS, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
324	MOPS-5203-10	MOPS, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
325	MOPS-5220-51	MOPS, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
326	MOPS-5220-92	MOPS, 100 G	Active	100.00 G		250.00 ML		White HDPE round pail, blue lid			
327	MOPS-5221-92	MOPS, 100 G	In development	100.00 G		250.00 ML		White HDPE round pail, blue lid			
328	NACL-3150-02	Sodium Chloride 5M Solution, 1 L	Active	1.00 L		1.00 L		Tepoflex sample			
329	NACL-3150-69	Sodium Chloride 5M Solution, 200 L	Active	200.00 L		58.00 GAL		Blue HDPE drum with metal ring for BPC			BPC bag
330	NACL-3150-76	Sodium Chloride 5M Solution, 1025 L	Active	1025.00 L		1000.00 L		BPC shell	1000L		BPC bag
331	NAHY-3150-76	Sodium Hydroxide 0.5N Solution, 1005 L	Active	1005.00 L		1000.00 L		BPC shell	1000L		BPC bag
332	NAHY-3153-75	Sodium Hydroxide 1N Solution, 1005 L	In development	1005.00 L		275.00 GAL		Tote			
333	NAHY-3153-76	Sodium Hydroxide 1N Solution, 1005 L	Active	1005.00 L		1000.00 L		BPC shell	1000L		BPC bag
334	NAHY-4122-70	Sodium Hydroxide 10N Solution, 1 L	Active	1.00 L		1.00 L		Non-fluorinated bottle, black lid			
335	NAHY-4122-72	Sodium Hydroxide 10N Solution, 200 L	Active	200.00 L		55.00 GAL		Natural HDPE drum			
336	NAHY-4122-75	Sodium Hydroxide 10N Solution, 940 L	Active	940.00 L		275.00 GAL		Tote			
337	NAHY-4122-96	Sodium Hydroxide 10N Solution, Bulk	Active	13059.00 L		13059.00 L		Bulk tanker			
338	NAHY-4123-70	Sodium Hydroxide 2N Solution, 1 L	Active	1.00 L		1.00 L		Non-fluorinated bottle, black lid			
339	NAHY-4123-71	Sodium Hydroxide 2N Solution, 1 L X 6 combo pack	Active	6.00 L		6.00 L		HDPE combo pack			
340	NAHY-4123-73	Sodium Hydroxide 2N Solution, 4 x 3.78 L	In development	15.12 L		15.12 L		4x1 Gallon combo pack			
341	NAHY-4124-78	Sodium Hydroxide 5N Solution, 10 L	Active	10.00 L		2.60 GAL		HDPE carboy			
342	NAHY-4126-71	Sodium Hydroxide 10N Solution Biotech, 6 X 1 L	Active	6.00 L		6.00 L		HDPE combo pack			
343	NAHY-4126-72	Sodium Hydroxide 10N Solution Biotech, 200 L	Active	200.00 L		55.00 GAL		Black HDPE drum			
344	NAHY-4126-73	Sodium Hydroxide 10N Solution Biotech, 4 X 3.78 L	Active	15.12 L		15.12 L		4x1 Gallon combo pack			
345	NAHY-4126-75	Sodium Hydroxide 10N Solution Biotech, 1000 L	Active	1000.00 L		275.00 GAL		Tote			
346	NAHY-4126-78	Sodium Hydroxide 10N Solution Biotech, 10 L	Active	10.00 L		2.60 GAL		HDPE carboy			
347	NAHY-4150-70	Sodium Hydroxide 10N Solution, 1 L	Active	1.00 L		1.00 L		Non-fluorinated bottle, black lid			
348	NAHY-4150-71	Sodium Hydroxide 10N Solution, 1 L X 6	Active	6.00 L		6.00 L		HDPE combo pack			
349	NAHY-4150-73	Sodium Hydroxide 10N Solution, 4 x 3.78 L	Active	15.12 L		15.12 L		4x1 Gallon combo pack			
350	NAHY-4150-78	Sodium Hydroxide 10N Solution, 10 L	Active	10.00 L		2.60 GAL		HDPE carboy			
351	NAHY-4150-79	Sodium Hydroxide 10N Solution, 942.5 L	Active	942.50 L		275.00 GAL		Tote			
352	NAHY-4150-81	Sodium Hydroxide 10N Solution, 4 x 3.78 L	Active	15.12 L		15.12 L		4x1 Gallon combo pack			
353	NAHY-4150-83	Sodium Hydroxide 10N Solution, 4 x 4 L	In development	16.00 L		16.00 L		4 x 4L combo pack			

	A	B	C	D	E	F	G	H	I	J	K
1	Product	Description 1	Product Status	Product size	Product UOM	Package size	Package UOM	Package type	Liner size	Liner type	Desiccant
354	NAHY-4151-70	Sodium Hydroxide 8N Solution, 1 L	Active	1.00 L		1.00 L		Non-fluorinated bottle, black lid			
355	NDEC-3201-01	Sodium Decanoate, 1kg	Active	1.00 KG		3.50 GAL		White HDPE round pail, blue lid	(2x) 24 x 36	LDPE	
356	NDEC-3201-03	Sodium Decanoate, 3 KG	Active	3.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
357	NDEC-3201-05	Sodium Decanoate, 5 KG	Active	5.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
358	NDEC-3201-10	Sodium Decanoate, 10 KG	Active	10.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
359	NDEC-3201-92	Sodium Decanoate, 100 gr	Active	100.00 G		1000.00 ML		White HDPE bottle, blue lid			
360	NDEC-3201-94	Sodium Decanoate, 250 gr	Active	250.00 G		1500.00 ML		White HDPE bottle, blue lid			
361	NDEC-3220-01	Sodium Decanoate, 1kg	Active	1.00 KG		3.50 GAL		White HDPE round pail, blue lid	(2x) 24 x 36	LDPE	
362	NDEC-3220-03	Sodium Decanoate, 3 KG	Active	3.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
363	NDEC-3220-05	Sodium Decanoate, 5 KG	Active	5.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
364	NDEC-3220-10	Sodium Decanoate, 10 KG	Active	10.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
365	NDEC-3220-91	Sodium Decanoate, 500 G	Active	500.00 G		0.60 GAL		White HDPE round pail, blue lid	9 x 15	LDPE	
366	NDEC-3220-92	Sodium Decanoate, 100 G	Active	100.00 G		1000.00 ML		White HDPE bottle, blue lid			
367	NDEC-3320-03	Sodium Decanoate, 3 KG	Active	3.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
368	NDEC-4201-01	Sodium Decanoate Biotech, 1 KG	Active	1.00 KG		3.50 GAL		White HDPE round pail, blue lid	(2x) 24 x 36	LDPE	
369	NDEC-4201-05	Sodium Decanoate Biotech, 5 KG	Active	5.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
370	NDEC-4201-10	Sodium Decanoate Biotech, 10 KG	Active	10.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
371	NDEC-4201-27	Sodium Decanoate Biotech, 25 KG	Active	25.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
372	NDEC-4201-92	Sodium Decanoate Biotech, 100 G	Active	100.00 G		1000.00 ML		White HDPE bottle, blue lid			
373	NDEC-4201-93	Sodium Decanoate Biotech, 500 G	Active	500.00 G		0.60 GAL		White HDPE round pail, blue lid	9 x 15	LDPE	
374	THCL-3203-01	Tris Hydrochloride, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
375	THCL-3203-05	Tris Hydrochloride, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
376	THCL-3203-10	Tris Hydrochloride, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
377	THCL-3203-25	Tris Hydrochloride, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
378	THCL-3203-92	Tris Hydrochloride, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
379	THCL-3203-93	Tris Hydrochloride, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
380	THCL-3220-01	Tris Hydrochloride, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
381	THCL-3220-05	Tris Hydrochloride, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
382	THCL-3220-10	Tris Hydrochloride, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
383	THCL-3220-12	Tris Hydrochloride, 12 KG	Active	12.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
384	THCL-3220-13	Tris Hydrochloride, 12 KG	Active	12.00 KG		5.00 GAL		White HDPE round pail, blue lid	(2x) 30 x 48	LDPE	
385	THCL-3220-25	Tris Hydrochloride, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
386	THCL-3220-50	Tris Hydrochloride, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	(2x) 34 x 60	LDPE	4 Unit
387	THCL-3220-58	Tris Hydrochloride, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	(1x) 34 x 60, (1x) 36 x 60	LDPE/Tyvek	5 x 8 Unit
388	THCL-3220-92	Tris Hydrochloride, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
389	THCL-3220-93	Tris Hydrochloride, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
390	THCL-3221-25	Tris Hydrochloride, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
391	THCL-3221-92	Tris Hydrochloride, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
392	THCL-3250-01	Tris Hydrochloride, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
393	THCL-3250-05	Tris Hydrochloride, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
394	THCL-3250-10	Tris Hydrochloride, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
395	THCL-3250-12	Tris Hydrochloride, 12 KG	Active	12.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
396	THCL-3250-13	Tris Hydrochloride, 12 KG	Active	12.00 KG		5.00 GAL		White HDPE round pail, blue lid	(2x) 30 x 48	LDPE	
397	THCL-3250-25	Tris Hydrochloride, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	(2x) 30 x 48	LDPE	

	A	B	C	D	E	F	G	H	I	J	K
1	Product	Description 1	Product Status	Product size	Product UOM	Package size	Package UOM	Package type	Liner size	Liner type	Desiccant
398	THCL-3250-28	Tris Hydrochloride, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
399	THCL-3250-41	Tris Hydrochloride, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
400	THCL-3250-42	Tris Hydrochloride, 12 KG	Active	12.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
401	THCL-3250-58	Tris Hydrochloride, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	(1x) 34 x 60, (1x) 36 x 60	LDPE/Tyvek	5 x 8 Unit
402	THCL-3250-91	Tris Hydrochloride, 50 G	Active	50.00 G		250.00 ML		White HDPE bottle, blue lid			
403	THCL-3250-92	Tris Hydrochloride, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
404	THCL-3250-93	Tris Hydrochloride, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
405	THCL-3253-10	Tris Hydrochloride, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
406	THCL-3253-12	Tris Hydrochloride, 12 KG	Active	12.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
407	THCL-3253-13	Tris Hydrochloride, 12 KG	Active	12.00 KG		5.00 GAL		White HDPE round pail, blue lid	(2x) 30 x 48	LDPE	
408	THCL-3253-92	Tris Hydrochloride, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
409	THCL-3253-93	Tris Hydrochloride, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
410	THCL-3254-01	Tris Hydrochloride, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
411	THCL-3254-12	Tris Hydrochloride, 12 KG	Active	12.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
412	THCL-3254-13	Tris Hydrochloride, 12 KG	Active	12.00 KG		5.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
413	THCL-3254-25	Tris Hydrochloride, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	(2x) 30 x 48	LDPE	
414	THCL-3254-50	Tris Hydrochloride, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	(2x) 34 x 60	LDPE	4 Unit
415	THCL-3254-92	Tris Hydrochloride, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
416	THCL-3254-93	Tris Hydrochloride, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
417	THCL-3257-10	Tris Hydrochloride, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
418	THCL-3257-93	Tris Hydrochloride, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
419	THCL-3259-01	Tris Hydrochloride, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
420	THCL-3259-05	Tris Hydrochloride, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
421	THCL-3259-10	Tris Hydrochloride, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
422	THCL-3259-12	Tris Hydrochloride, 12 KG	Active	12.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
423	THCL-3259-25	Tris Hydrochloride, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
424	THCL-3259-58	Tris Hydrochloride, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	(1x) 34 x 60, (1x) 36 x 60	LDPE/Tyvek	5 x 8 Unit
425	THCL-3259-91	Tris Hydrochloride, 250 G	Active	250.00 G		1000.00 ML		White HDPE bottle, blue lid			
426	THCL-3259-92	Tris Hydrochloride, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
427	THCL-3259-93	Tris Hydrochloride, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
428	THCL-3260-01	Tris Hydrochloride, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
429	THCL-3260-92	Tris Hydrochloride, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
430	THCL-4251-01	Tris Hydrochloride Biotech, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
431	THCL-4251-05	Tris Hydrochloride Biotech, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
432	THCL-4251-10	Tris Hydrochloride Biotech, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
433	THCL-4251-25	Tris Hydrochloride Biotech, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
434	THCL-4251-50	Tris Hydrochloride Biotech, 50 KG	In development	50.00 KG		30.00 GAL		Dark blue HDPE drum	(2x) 34 x 60	LDPE	4 Unit
435	THCL-4251-92	Tris Hydrochloride Biotech, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
436	THCL-4251-93	Tris Hydrochloride Biotech, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
437	TRED-3250-03	Trehalose Dihydrate, 3 KG	Active	3.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
438	TRED-3250-05	Trehalose Dihydrate, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
439	TRED-3250-06	Trehalose Dihydrate, 1 KG	Active	1.00 KG		0.60 GAL		White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	
440	TRED-3250-10	Trehalose Dihydrate, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
441	TRED-3250-25	Trehalose Dihydrate, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	

	A	B	C	D	E	F	G	H	I	J	K
1	Product	Description 1	Product Status	Product size	Product UOM	Package size	Package UOM	Package type	Liner size	Liner type	Desiccant
442	TRED-3250-50	Trehalose Dihydrate, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
443	TRED-3252-06	Trehalose Dihydrate, 1 KG	Active	1.00 KG		0.60 GAL		White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	
444	TRED-3252-10	Trehalose Dihydrate, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
445	TRED-3252-25	Trehalose Dihydrate, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
446	TRED-3253-06	Trehalose Dihydrate, 1 KG	In development	1.00 KG		0.60 GAL		White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	
447	TRED-3253-10	Trehalose Dihydrate, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
448	TRED-4250-05	Trehalose Dihydrate Biotech, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
449	TRED-4250-06	Trehalose Dihydrate Biotech, 1 KG	Active	1.00 KG		0.60 GAL		White HDPE round pail, blue lid	(2x) 9 x 15	LDPE	
450	TRED-4250-10	Trehalose Dihydrate Biotech, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
451	TRED-4250-25	Trehalose Dihydrate Biotech, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
452	TRED-4250-50	Trehalose Dihydrate Biotech, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
453	TRIS-1201-01	Tris/Tromethamine, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
454	TRIS-1201-10	Tris/Tromethamine, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
455	TRIS-1201-25	Tris/Tromethamine, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
456	TRIS-1201-91	Tris/Tromethamine, 250 G	Active	250.00 G		1000.00 ML		White HDPE bottle, blue lid			
457	TRIS-1201-92	Tris/Tromethamine, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
458	TRIS-1201-93	Tris/Tromethamine, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
459	TRIS-1220-01	Tris/Tromethamine, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
460	TRIS-1220-25	Tris/Tromethamine, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
461	TRIS-1220-91	Tris/Tromethamine, 250 G	Active	250.00 G		1000.00 ML		White HDPE bottle, blue lid			
462	TRIS-1220-92	Tris/Tromethamine, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
463	TRIS-1220-93	Tris/Tromethamine, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
464	TRIS-2255-01	Tris/Tromethamine, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
465	TRIS-2255-10	Tris/Tromethamine, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
466	TRIS-2255-91	Tris/Tromethamine, 250 G	Active	250.00 G		1000.00 ML		White HDPE bottle, blue lid			
467	TRIS-2255-92	Tris/Tromethamine, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
468	TRIS-2255-93	Tris/Tromethamine, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
469	TRIS-2257-01	Tris/Tromethamine, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
470	TRIS-2257-25	Tris/Tromethamine, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
471	TRIS-2257-91	Tris/Tromethamine, 250 G	Active	250.00 G		1000.00 ML		White HDPE bottle, blue lid			
472	TRIS-2257-92	Tris/Tromethamine, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
473	TRIS-2257-93	Tris/Tromethamine, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
474	TRIS-3201-01	Tris/Tromethamine, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
475	TRIS-3201-92	Tris/Tromethamine, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
476	TRIS-3201-93	Tris/Tromethamine, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
477	TRIS-3220-01	Tris/Tromethamine, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
478	TRIS-3220-10	Tris/Tromethamine, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
479	TRIS-3220-25	Tris/Tromethamine, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
480	TRIS-3220-50	Tris/Tromethamine, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	(2x) 34 x 60	LDPE	4 Unit
481	TRIS-3220-92	Tris/Tromethamine, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
482	TRIS-3220-93	Tris/Tromethamine, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
483	TRIS-3251-01	Tris/Tromethamine, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
484	TRIS-3251-05	Tris/Tromethamine, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
485	TRIS-3251-10	Tris/Tromethamine, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	

	A	B	C	D	E	F	G	H	I	J	K
1	Product	Description 1	Product Status	Product size	Product UOM	Package size	Package UOM	Package type	Liner size	Liner type	Desiccant
486	TRIS-3251-92	Tris/Tromethamine, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
487	TRIS-3251-95	Tris/Tromethamine, 80 G	Active	80.00 G		1000.00 ML		White HDPE bottle, blue lid	9 x 15	LDPE	
488	TRIS-3252-25	Tris/Tromethamine, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
489	TRIS-3252-60	Tris/Tromethamine, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	(1x) 34 x 60, (2x) 36 x 60	LDPE/(2x) Tyvek	(2x) 5 x 8 Unit
490	TRIS-3254-01	Tris/Tromethamine, 1KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
491	TRIS-3254-05	Tris/Tromethamine, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
492	TRIS-3254-10	Tris/Tromethamine, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
493	TRIS-3254-25	Tris/Tromethamine, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
494	TRIS-3254-50	Tris/Tromethamine, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	(2x) 34 x 60	LDPE	4 Unit
495	TRIS-3254-58	Tris/Tromethamine, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	(1x) 34 x 60, (1x) 36 x 60	LDPE/Tyvek	5 x 8 Unit
496	TRIS-3254-92	Tris/Tromethamine, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
497	TRIS-3254-93	Tris/Tromethamine, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
498	TRIS-3255-01	Tris/Tromethamine, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
499	TRIS-3255-05	Tris/Tromethamine, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
500	TRIS-3255-07	Tris/Tromethamine, 2 KG	Active	2.00 KG		1.25 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
501	TRIS-3255-10	Tris/Tromethamine, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
502	TRIS-3255-12	Tris/Tromethamine, 12 KG	Active	12.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
503	TRIS-3255-25	Tris/Tromethamine, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
504	TRIS-3255-41	Tris/Tromethamine, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
505	TRIS-3255-42	Tris/Tromethamine, 12 KG	Active	12.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
506	TRIS-3255-50	Tris/Tromethamine, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	(2x) 34 x 60	LDPE	4 Unit
507	TRIS-3255-58	Tris/Tromethamine, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	(1x) 34 x 60, (1x) 36 x 60	LDPE/Tyvek	5 x 8 Unit
508	TRIS-3255-89	Tris/Tromethamine, 200 G	Active	200.00 G		1000.00 ML		White HDPE bottle, blue lid			
509	TRIS-3255-91	Tris/Tromethamine, 50 G	Active	50.00 G		250.00 ML		White HDPE bottle, blue lid			
510	TRIS-3255-92	Tris/Tromethamine, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
511	TRIS-3255-93	Tris/Tromethamine, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
512	TRIS-3255-94	Tris/Tromethamine, 250 G	Active	250.00 G		1000.00 ML		White HDPE bottle, blue lid			
513	TRIS-3257-10	Tris/Tromethamine, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
514	TRIS-4223-01	Tris (Tromethamine) Biotech, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
515	TRIS-4223-05	Tris (Tromethamine) Biotech, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
516	TRIS-4223-10	Tris (Tromethamine) Biotech, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
517	TRIS-4223-25	Tris (Tromethamine) Biotech, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
518	TRIS-4223-50	Tris (Tromethamine) Biotech, 50 KG	In development	50.00 KG		30.00 GAL		Dark blue HDPE drum	(2x) 34 x 60	LDPE	4 Unit
519	TRIS-4223-92	Tris (Tromethamine) Biotech, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
520	TRIS-4223-93	Tris (Tromethamine) Biotech, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
521	TRIS-5201-26	Tris/Tromethamine, 25 KG	Active	25.00 KG		26.00 GAL		Fiber drum	(2x) 34 x 60	LDPE	4 Unit
522	TRIS-5203-59	Tris/Tromethamine, 50 KG	Active	50.00 KG		26.00 GAL		Fiber drum	(2x) 34 x 60	LDPE	4 Unit
523	TRIS-5204-26	Tris/Tromethamine, 25 KG	Active	25.00 KG		26.00 GAL		Fiber drum	(2x) 34 x 60	LDPE	
524	TRIS-5207-26	Tris/Tromethamine, 25 KG	Active	25.00 KG		26.00 GAL		Fiber drum	(2x) 34 x 60	LDPE	4 Unit
525	TRIS-5220-10	Tris/Tromethamine, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
526	TRIS-5222-10	Tris/Tromethamine, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
527	TRIS-5223-89	Tris/Tromethamine, 200 G	Active	200.00 G		1000.00 ML		White HDPE bottle, blue lid			
528	TRIS-6201-01	Tris/Tromethamine, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
529	URAC-4201-05	Uracil, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	

	A	B	C	D	E	F	G	H	I	J	K
1	Product	Description 1	Product Status	Product size	Product UOM	Package size	Package UOM	Package type	Liner size	Liner type	Desiccant
530	URAC-4202-01	Uracil, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
531	URAC-4202-05	Uracil, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
532	URAC-4202-91	Uracil, 50 G	Active	50.00 G		250.00 ML		White HDPE bottle, blue lid			
533	URAC-4202-92	Uracil, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
534	URAC-4202-95	Uracil, 10g	Active	10.00 G		60.00 ML		White polypropylene bottle			
535	URAC-4250-01	Uracil, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
536	URAC-4250-95	Uracil, 10g	Active	10.00 G		60.00 ML		White polypropylene bottle			
537	URAC-4301-05	Uracil, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	(2x) 24 x 36	LDPE	
538	URAC-4301-06	Uracil, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	(5x) 24 x 36	LDPE	
539	UREA-3120-02	Urea 6M Solution, 500 ML	Active	500.00 ML		17.00 OZ		Black HDPE UN leakproof bottle			
540	UREA-3120-76	Urea 6M Solution, 750 L	In development	750.00 L		1000.00 L		BPC shell	1000L	BPC bag	
541	UREA-3120-78	Urea 6M Solution, 1200 L	Active	1200.00 L		330.00 GAL		Tote		LDPE	
542	UREA-3220-10	Urea, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
543	UREA-3220-25	Urea, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	34 x 60	LDPE	
544	UREA-3220-50	Urea, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	(1x) 34 x 60, (1x) 36 x 60	LDPE/Tyvek	5 x 8 Unit
545	UREA-3220-89	Urea, 200g	Active	200.00 G		1000.00 ML		White HDPE bottle, blue lid			
546	UREA-3220-92	Urea, 100g	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
547	UREA-3220-93	Urea, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
548	UREA-3222-10	Urea, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
549	UREA-3222-25	Urea, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	34 x 60	LDPE	
550	UREA-3223-25	Urea, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	34 x 60	LDPE	
551	UREA-3223-58	Urea, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	(1x) 34 x 60, (1x) 36 x 60	LDPE/Tyvek	5 x 8 Unit
552	UREA-3250-01	Urea, 1KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
553	UREA-3250-12	Urea, 12 KG	Active	12.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	4 Unit
554	UREA-3250-25	Urea, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	34 x 60	LDPE	
555	UREA-3250-92	Urea, 100g	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
556	UREA-3250-93	Urea, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
557	UREA-4221-01	Urea Biotech, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
558	UREA-4221-05	Urea Biotech, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
559	UREA-4221-10	Urea Biotech, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
560	UREA-4221-25	Urea Biotech, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
561	UREA-4221-50	Urea Biotech, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
562	UREA-4221-92	Urea Biotech, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
563	UREA-4221-93	Urea Biotech, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
564	UREA-5203-05	Urea, 5 KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
565	UREA-5203-15	Urea, 15 KG	Active	15.00 KG		12.00 GAL		White HDPE round pail, blue lid	34 x 60	LDPE	
566	UREA-5220-25	Urea, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
567	URID-3250-01	Uridine, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			
568	URID-3250-10	Uridine, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
569	URID-3250-25	Uridine, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
570	URID-3250-91	Uridine, 250 G	Active	250.00 G		1000.00 ML		White HDPE bottle, blue lid			
571	URID-3250-92	Uridine, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
572	URID-3250-93	Uridine, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
573	URID-4250-01	Uridine Biotech, 1 KG	Active	1.00 KG		1500.00 ML		White HDPE bottle, blue lid			

	A	B	C	D	E	F	G	H	I	J	K
1	Product	Description 1	Product Status	Product size	Product UOM	Package size	Package UOM	Package type	Liner size	Liner type	Desiccant
574	URID-4250-05	Uridine Biotech, 5KG	Active	5.00 KG		3.50 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
575	URID-4250-10	Uridine Biotech, 10 KG	Active	10.00 KG		5.00 GAL		White HDPE round pail, blue lid	24 x 36	LDPE	
576	URID-4250-25	Uridine Biotech, 25 KG	Active	25.00 KG		12.00 GAL		White HDPE round pail, blue lid	30 x 48	LDPE	
577	URID-4250-50	Uridine Biotech, 50 KG	Active	50.00 KG		30.00 GAL		Dark blue HDPE drum	34 x 60	LDPE	
578	URID-4250-92	Uridine Biotech, 100 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
579	URID-4250-93	Uridine Biotech, 500 G	Active	500.00 G		1000.00 ML		White HDPE bottle, blue lid			
580	URID-4250-95	Uridine Biotech, 50 G	Active	100.00 G		250.00 ML		White HDPE bottle, blue lid			
581	WAFI-3150-68	Water for Injection (WFI), 200L	In development	200.00 L		200.00 L		BPC shell	200L		Typoflex BPC t
582	WAFI-3150-74	Water for Injection (WFI), 3.78L	Active	3.78 L		3.78 L		Non-fluorinated HDPE bottle			