

## 8-UNIT STRIP OF 5 CLAY DESICCANTS SPECIFICATIONS

BSI COMPONENT NUMBER: COMP-0006



8-Unit Clay in Tyvek – 5 Pack Strip of Desiccant		
Test		Specifications
Desiccant Name		4-Unit Tyvek Desiccant 5 Pack Strip – 8 Unit Clay in Tyvek Desiccant
Desiccant Type		Desi Pak® Unit Pak
Appearance and Color		Sealed White Tyvek® package with Gray print
Stamp Markings		Specification MIL-D-3464 Type I & II
Seal Integrity		Physically inspect the Seal and ensure it is fully intact with the seal firmly applied
Clay Moisture Pickup Percent	@40% RH, 25°C, 4 Liters/Min	17.3% min
	@80% RH, 25°C, 4 Liters/Min	23.0% min
Clay Free Moisture Percent Per Computrac Moisture Analyzer		3.0% max

For further information, please contact [marketing@biospectra.us](mailto:marketing@biospectra.us).

**Connie Augustine | Director**  
**Business Integration & Retention**  
**BioSpectra Inc.**

*This is not considered a controlled document. We are not responsible for any errors or omissions, and the user is responsible for any decisions based on the information herein.*