

## SUPERSACK EU COMPLIANCE

BioSpectra has evaluated the packaging component supply through the Supplier Approval Program. The supersack or flexible intermediate bulk container (FIBC) utilized by BioSpectra for material packaging is produced using LDPE materials that meet the compositional requirements of Monograph 3.1.4 “Polyethylene without additives for containers for preparations for parenteral use and for ophthalmic preparations” of the European Pharmacopeia (10<sup>th</sup> Edition).

| Product Number | Packaging Size |
|----------------|----------------|
| COMP-0166      | 600kg          |

For further information, please contact [info@biospectra.us](mailto:info@biospectra.us)

*Cassie Baun*

**Cassie Baun**

Senior Compliance Specialist