

BIOSPECTRA

100 Majestic Way, Bangor, PA 18013 / www.biospectra.us

Effective Date:	29-Sep-2021	29-Sep-2024	: Date of Next Review
Prepared By:	Jaron Hughes	16-001182 v.7.0	: Supersedes
QA/QC Approval:	Krista Rehrig	Carissa McCollan	: Management Approval
Reason for Revision:	See Revision History in ensur.		

CERTIFICATE OF ANALYSIS UREA

BIO EXCIPIENT GRADE / NEW CODE UREA-3220-93
(HISTORICAL CODE UR3220-G500)

LOT: UREA-0122-00045

NH_2CONH_2 ^ F.W. 60.06 g/mol. ^ CAS# 57-13-6

Manufacturing Date: 8/15/21 Retest Date: 8/31/23

Manufacturing Site: 1474 Rockdale Lane, Stroudsburg, PA 18360

Packaging Date: 9/17/22 Packaging Site: 100 Majestic Way, Bangor PA, 18013

Meets or Exceeds USP Specifications

ANALYSIS	SPECIFICATION	TEST RESULT	
Alcohol Insoluble Matter	0.04% max.	< 0.04%	
Appearance and Color	White / Crystals	White / Crystals	
Assay	98.0-102.0%	98.7%	
Enzymes	DNase	None Detected	
	Protease	None Detected	
	RNase	None Detected	
Heavy Metals	10 ppm max.	< 10 ppm	
Identification A(IR)	Passes Test	Passes Test	
Identification B	Passes Test	Passes Test	
Impurities	Urea RCA	≤ 0.1%	
	Total	≤ 2.0%	
	Unspecified	≤ 0.1%	
Insoluble Matter	0.010% max.	< 0.001%	
Loss on Drying	1.0% max.	0.2%	
Melting Range	132-135 °C	133 - 134°C	
Residue on Ignition	0.010% max.	< 0.005%	
	Arsenic (As)	5 ppm max.	< 5 ppm
	Copper (Cu)	5 ppm max.	< 5 ppm
	Iron (Fe)	5 ppm max.	< 5 ppm
	Lead (Pb)	5 ppm max.	< 5 ppm

COUNTRY OF ORIGIN: U.S.A.

TEST METHOD REFERENCE: DCN: 16-000495

The information contained herein is the property of BioSpectra. The recipient is responsible for its safe-keeping, and the prevention of unauthorized appropriation, use, disclosure and copying.

INTENDED USE: Material represented by this Certificate of Analysis is suitable for use as an excipient. It is manufactured in accordance with the ICH Q7 Good Manufacturing Practice Guide. The material represented by this Certificate of Analysis is not suitable to be used as an Active Pharmaceutical Ingredient, Drug Product or Household Item.

RESIDUAL SOLVENTS STATEMENT: Based on the manufacturing process and the controlled handling, storage and analysis of this product, this product complies with the requirements and specifications listed in the current USP method <467> Tables 1, 2, 3, or 4.

Prepared by: BB Date: 9/20/22 Job Title: QA Specialist

Reviewed by: Cassie Allert Date: 9/20/22 Job Title: QA Manager