


48" x 48" x 48" WOOD CRATE SPECIFICATIONS

BSI COMPONENT NUMBER: COMP-0144



48x48x48 Wood Crate with Pallet		
Test		Specifications
Appearance and Color		6-Piece plywood box with wooden pallet. Container is free of dirt, debris, and pests.
Panel Dimensions (LxWxT) (Nominal)	Side 1	48" x 43" x 1/4"
	Side 2	48" x 43" x 1/4"
	Top	48" x 43" x 1/4"
Pallet Height (Nominal)		5"
Wood Pallet Heat Treatment Stamp Presence		 <div>CA-00009 HT</div>

For further information, please contact marketing@biospectra.us.

**Connie Augustine | Director
Business Integration & Retention
BioSpectra Inc.**

This is not considered a controlled document. We are not responsible for any errors or omissions, and the user is responsible for any decisions based on the information herein.