

BIOBUFFER SOLUTIONS

100 Majestic Way, Bangor, PA 18013 / www.biospectra.us

March 7th, 2024

Revision 1

RESIDUAL SOLVENT STATEMENT

L-Arginine Hydrochloride

L-Arginine Hydrochloride, Bio Pharma Grade for BioBuffer Solutions Product line is manufactured using the solvent Ethyl Acetate, and is stated to comply with the ICH Q3C (R8) Guideline for Residual Solvents required concentration limit of 5000ppm for a Class 3 residual solvent. BioSpectra has evaluated the material supply through the Supplier Qualification Program.

Current Product Number
LARH-4220

For further information, please contact info@biospectra.us



Cassie Baun

Senior Compliance Specialist