



100 Majestic Way, Bangor, PA 18013 / www.biospectra.us

HDPE ROUND SCREW TOP PAIL AND LID COMPONENT SPECIFICATIONS

0.6 Gallon HDPE Round Screw Top Pail and Lid	
Test	Specifications
Appearance and Color	White HDPE round pail with hole(s) for tamper evident seal, blue tamper evident lid with gasket, red trigger, and slot on top of lid
Inside Height	6.75" \pm 0.5"
Inside Top Diameter	6.15" \pm 0.5"
Outside Height, with Lid	7.0" \pm 0.5"
Outside Bottom Diameter	5.60" \pm 0.5"
UN Rating	1H2/Y4/S/YR
Lid Engraving (top)	Engraving must reflect manufacture by M&M Industries. Specific M&M logo may vary.
Pail Engraving (bottom)	Engraving must reflect manufacture by M&M Industries. Specific M&M logo may vary.

M&M Industries, Inc. Part #: 530031
M&M Industries, Inc. Lid Part #: 600010
Uline Part # Pail: S-20029P
Uline Part # Lid: S-20029L-BLU
Uline Part # Kit (Pail and Lid): S-20029 BLU
Supplied by: M&M Industries, Inc., or Uline
Manufactured by: M&M Industries, Inc.

1.25 Gallon HDPE Round Screw Top Pail and Lid	
Test	Specifications
Appearance and Color	White HDPE round pail with hole(s) for tamper evident seal, blue tamper evident lid with gasket, red trigger, and slot on top of lid
Inside Height	8.25" ± 0.5"
Inside Top Diameter	7.77" ± 0.5"
Outside Height, with Lid	9.14" ± 0.5"
Outside Bottom Diameter	7.06" ± 0.5"
UN Rating	1H2/Y6/S/YR
Lid Engraving (top)	Engraving must reflect manufacture by M&M Industries. Specific M&M logo may vary.
Pail Engraving (bottom)	Engraving must reflect manufacture by M&M Industries. Specific M&M logo may vary.

Uline Part # Pail: S-18114P
 Uline Part # Lid: S-18114L-BLU
 Uline Part # Kit (Pail and Lid): S-18114BLU
 Supplied by: Uline
 Manufactured by: M&M Industries, Inc.

3.5 Gallon HDPE Round Screw Top Pharma Pail and Lid	
Test	Specifications
Appearance and Color	White HDPE round pail with holes for tamper evident seal, blue tamper evident lid with gasket, red trigger, and slot on top of lid
Inside Height	11.0" ± 0.5"
Inside Top Diameter	11.25" ± 0.5"
Outside Height, with Lid	12.5" ± 0.5"
Outside Bottom Diameter	10.5" ± 0.5"
UN Rating	1H2/Y19/S/YR
Lid Engraving (top)	Engraving must reflect manufacture by M&M Industries. Specific M&M logo may vary.
Pail Engraving (bottom)	Engraving must reflect manufacture by M&M Industries. Specific M&M logo may vary.

Part #: 530049
 Lid Part #: 630054
 Supplied by: M&M Industries, Inc
 Manufactured by: M&M Industries, Inc.

5 Gallon HDPE Round Screw Top Pharma Pail and Lid	
Test	Specifications
Appearance and Color	White HDPE round pail with holes for tamper evident seal, white tamper evident lid with gasket, red trigger, and slot on top of lid
Inside Height	15.0" ± 0.5"
Inside Top Diameter	11.25" ± 0.5"
Outside Height, with Lid	16.0" ± 0.5"
Outside Bottom Diameter	10.25" ± 0.5"
UN Rating	1H2/Y30/S/YR 1H2/X11.5/S/YR
Lid Engraving (top)	Engraving must reflect manufacture by M&M Industries. Specific M&M logo may vary.
Pail Engraving (bottom)	Engraving must reflect manufacture by M&M Industries. Specific M&M logo may vary.

Uline Part # Kit (Pail and Lid): S-13652W

Supplied by: Uline

Manufactured by: M&M Industries, Inc.

5 Gallon HDPE Round Screw Top Pharma Pail and Lid	
Test	Specifications
Appearance and Color	White HDPE round pail with holes for tamper evident seal, blue tamper evident lid with gasket, red trigger, and slot on top of lid
Inside Height	15.0" ± 0.5"
Inside Top Diameter	11.25" ± 0.5"
Outside Height, with Lid	16.0" ± 0.5"
Outside Bottom Diameter	10.25" ± 0.5"
UN Rating	1H2/Y30/S/YR 1H2/X11.5/S/YR
Lid Engraving (top)	Engraving must reflect manufacture by M&M Industries. Specific M&M logo may vary.
Pail Engraving (bottom)	Engraving must reflect manufacture by M&M Industries. Specific M&M logo may vary.

M&M Industries, Inc. Part #: 511589

M&M Industries, Inc. Lid Part #: 630054

Uline Part # Pail: S-13652P

Uline Part # Lid: S-13652L-BLU

Uline Part # Kit (Pail and Lid): S-13652BLU

Supplied by: M&M Industries, Inc. or Uline

Manufactured by: M&M Industries, Inc.

12 Gallon HDPE Round Screw Top Pail and Lid	
Test	Specifications
Appearance and Color	White HDPE round pail with tab insert, blue tamper evident lid with gasket, red trigger, and slot on side of lid
Inside Height	21.0" \pm 0.5"
Inside Top Diameter	14.0" \pm 0.5"
Outside Height, with Lid	21.5" \pm 0.5"
Outside Bottom Diameter	13.0" \pm 0.5"
UN Rating	1H2/Y40/S/YR with gasket 1H2/Y53/S/YR without gasket
Lid Engraving (top)	Engraving must reflect manufacture by M&M Industries. Specific M&M logo may vary.
Pail Engraving (bottom)	Engraving must reflect manufacture by M&M Industries. Specific M&M logo may vary.

M&M Industries, Inc. Part #: 530029

M&M Industries, Inc. Blue Lid Part Number: 630141 (Pharma) or 600112 (Non-Pharma)

Uline Part # Pail: S-23730P

Uline Part # Lid: S-23730L-BLU

Uline Part # Kit (Pail and Lid): S-23730BLU

Supplied by: M&M Industries, Inc., or Uline

Manufactured by: M&M Industries, Inc.

12 Gallon HDPE Round Screw Top Pharma Pail and Lid	
Test	Specifications
Appearance and Color	White HDPE round pail with tab insert, green tamper evident lid with gasket, report slot location (if applicable)
Inside Height	21.0" \pm 0.5"
Inside Top Diameter	14.0" \pm 0.5"
Outside Height, with Lid	21.5" \pm 0.5"
Outside Bottom Diameter	13.0" \pm 0.5"
UN Rating	1H2/Y40/S/YR with gasket 1H2/Y53/S/YR without gasket
Lid Engraving (top)	Engraving must reflect manufacture by M&M Industries. Specific M&M logo may vary.
Pail Engraving (bottom)	Engraving must reflect manufacture by M&M Industries. Specific M&M logo may vary.

Pail Part #: 530029

Green Lid Part Number: 600649

Supplied by: M&M Industries, Inc.

Manufactured by: M&M Industries, Inc.

10 Liter HDPE Round Click Pack y Pail and Lid	
Test	Specifications
Appearance and Color	White HDPE round pail with holes for tamper evident seal, blue tamper evident lid with white gasket, slot on side of lid
Inside Height	8.0" \pm 0.5"
Inside Top Diameter	11.0" \pm 0.5"
Outside Height, with Lid	8.5" \pm 0.5"
Outside Bottom Diameter	10.5" \pm 0.5"
UN Rating	1H2/X8/Y12/Z18/S/YR
Lid Engraving (top)	Engraving must reflect manufacture by CurTec
Pail Engraving (bottom)	Engraving must reflect manufacture by CurTec

Pail Part #: 4410

Lid Part #: 4402

Supplied by: CurTec

Manufactured by: CurTec

15 Liter HDPE Round Click Pack y Pail and Lid	
Test	Specifications
Appearance and Color	White HDPE round pail with holes for tamper evident seal, blue tamper evident lid with white gasket, slot on side of lid
Inside Height	10.0" \pm 0.5"
Inside Top Diameter	11.0" \pm 0.5"
Outside Height, with Lid	11.0" \pm 0.5"
Outside Bottom Diameter	10.5" \pm 0.5"
UN Rating	1H2/X14/Y21/Z30/S/YR
Lid Engraving (top)	Engraving must reflect manufacture by CurTec
Pail Engraving (bottom)	Engraving must reflect manufacture by CurTec

Pail Part #: 4415

Lid Part #: 4402

Supplied by: CurTec

Manufactured by: CurTec

20 Liter HDPE Round Click Pack x Pail and Lid	
Test	Specifications
Appearance and Color	White HDPE round pail with holes for tamper evident seal, blue tamper evident lid with white gasket, slot on side of lid
Inside Height	14.0" \pm 0.5"
Inside Top Diameter	11.5" \pm 0.5"
Outside Height, with Lid	15.0" \pm 0.5"
Outside Bottom Diameter	10.0" \pm 0.5"
UN Rating	1H2/X30/S/YR
Lid Engraving (top)	Engraving must reflect manufacture by CurTec
Pail Engraving (handgrip)	Engraving must reflect manufacture by CurTec

Pail Part #: 4520

Lid Part #: 4502

Supplied by: CurTec

Manufactured by: CurTec