

12/19/2025

Version 1.1

TYVEK LINER SPECIFICATIONS

BSI COMPONENT NUMBER: COMP-0019



Tyvek Liner Component Specifications	
Test	Specifications
Appearance and Color	White with no punctures
Width	36.0" \pm 0.5"
Length (Internal Tolerance)	60.0" \pm 1.0"

For further information, please contact marketing@biospectra.us.



**Connie Augustine | Director
Business Integration & Retention
BioSpectra Inc.**

This is not considered a controlled document. We are not responsible for any errors or omissions, and the user is responsible for any decisions based on the information herein.