

# BIOBUFFER SOLUTIONS

100 Majestic Way, Bangor, PA 18013 / [www.biospectra.us](http://www.biospectra.us)

January 10<sup>th</sup>, 2025

Revision 1

## SUPPLY CHAIN STATEMENT

### L-Histidine Monochloride Monohydrate

L-Histidine Monochloride Monohydrate, Bio Qualified Grade is manufactured and packaged by BioSpectra's approved supplier in accordance with the approved supplier's ISO 9001:2015 certified management system. Additional packaging of L-Histidine Monochloride Monohydrate, Bio Qualified Grade is performed at BioSpectra's Bangor, PA Facility. The primary storage and shipment location for L-Histidine Monochloride Monohydrate, Bio Qualified Grade is BioSpectra's Supply Chain Center: 51 North 3<sup>rd</sup> Street, Stroudsburg, PA 18360. BioSpectra may additionally utilize the Bangor, PA facility: 100 Majestic Way, Bangor, PA 18013 as a qualified secondary storage and shipment location as needed.

Current Product Number
LHMM-6250

For further information, please contact [info@biospectra.us](mailto:info@biospectra.us)



**Cassie Baun**

Senior Compliance Specialist