


12/19/2025

Version 1.1

WHITE STYROFOAM CROSS SPECIFICATIONS

BSI COMPONENT NUMBER: COMP-0057



White Styrofoam Cross	
Test	Specifications
Material	2# EPS Foam
Appearance and Color	White Cross
Length (Nominal)	24" \pm 0.25"
Width (Nominal)	24" \pm 0.25"
Thickness (Nominal)	4.0" \pm 0.25"

For further information, please contact marketing@biospectra.us.

**Connie Augustine | Director
Business Integration & Retention
BioSpectra Inc.**

This is not considered a controlled document. We are not responsible for any errors or omissions, and the user is responsible for any decisions based on the information herein.