

Bis-TRIS

BET*, GMP

CAS #:6976-37-0

Formula: $C_8H_{19}NO_5$

F.W.: 209.24 g/mol

BTRI-4250

BIO PHARMA GRADE

ANALYSIS		SPECIFICATIONS
Absorbance (0.1M)	280 nm	≤ 0.04 a.u.
	340 nm	≤ 0.02 a.u.
Ammonium		≤ 500 ppm
Appearance and Color		White Crystalline Powder
Assay, Dried Basis		99.0 – 101.0%
Microbial Content	TAMC	≤ 100 CFU/g
	TYMC	≤ 100 CFU/g
Chloride		≤ 50 ppm
Endotoxin		≤ 50 EU/g
Enzymes	DNase	None Detected
	Protease	None Detected
	RNase	None Detected
Identification, IR		Conforms to Reference Standard
Loss on Drying		$\leq 1.0\%$
Melting Point		100 – 105°C
pH (1%)		8.8 – 9.8
pKa		6.3 – 6.7
Residue on Ignition		$\leq 0.1\%$
Solubility (0.1M)		Clear and Complete
Sulfate		≤ 50 ppm
Trace Metals	Aluminum (Al)	≤ 5 ppm
	Calcium (Ca)	≤ 500 ppm
	Iron (Fe)	≤ 3 ppm
	Lead (Pb)	≤ 3 ppm
	Phosphorus (P)	≤ 500 ppm
Water, KF		$\leq 1.0\%$

*Bioburden and Endotoxin Tested

General Product Overview

Bis-Tris is a zwitterionic biological “Goods” buffer, often used in biological pharmaceutical manufacturing and biochemistry for extraction, separation, and purification, as well as an intermediate in synthetic reactions for both proteins and nucleic acid. Bis-Tris can also be used as a process aid in peptide synthesis, or as a scavenger agent for metal ions such as Cu and Pb.

Industry Application

Suitable for use in biological and biotech chemical process applications from R&D through scale cGMP production.

Key Product Features

- Appears as a white crystalline powder
- Manufactured in accordance with IPEC
- Manufactured in an enzyme free, hormone free and animal free environment
- Contains no known major food allergens (as defined by FDA and WHO)
- The final product and its raw materials are not derived from nor come into contact with animal parts, animal products, and/or animal byproducts or derivatives.
- Is not subject to genetic modification
- Synonyms: 2,2-Bis(hydroxymethyl)-2,2',2''- nitrilotriethanol; Bis(2-hydroxyethyl) amino tris (hydroxymethyl)methane

Storage and Shipping Conditions

Refer to SDS.

Standard Shelf-Life Policy

Unless otherwise noted on the Shelf-Life Statement and CoA, this product has a 2-year retest date supported by a 3-year ICH Q1 Stability Study (if one is completed).

Package Sizes

100g, 500g, 1kg, 5kg, 10kg, 25kg, 50kg

[Click here to view SDS, CoAs and other supporting regulatory documents on our website.](#)

This is not considered a controlled document. We are not responsible for any errors or omissions, and the user is responsible for any decisions based on the information herein.