

12/19/2025

Version 1.1

## 4-UNIT DESICCANT SPECIFICATIONS

BSI COMPONENT NUMBER: COMP-0004



4-Unit Tyvek Desiccant	
Test	Specifications
Desiccant Name	4-Unit Tyvek Desiccant
Desiccant Type	Desi Pak®
Appearance and Color	Sealed White Tyvek® package with Black or Blue print
Stamp Markings	Specification MIL-D-3464 Type I & II
Seal Integrity	Physically inspect the Seal and ensure it is fully intact with the seal firmly applied

For further information, please contact [marketing@biospectra.us](mailto:marketing@biospectra.us).

**Connie Augustine | Director  
Business Integration & Retention  
BioSpectra Inc.**

*This is not considered a controlled document. We are not responsible for any errors or omissions, and the user is responsible for any decisions based on the information herein.*