

BIOSPECTRA

100 Majestic Way, Bangor, PA 18013 / www.biospectra.us

Effective Date:	22-May-2020	22-May-2023	: Date of Next Review
Prepared By:	Kyle Snyder	16-001173 v.3.0	: Supersedes
QA/QC Approval:	Carissa McCollian	Hannah Bernier	: Management Approval
Reason for Revision:	See Revision History in ensur.		

CERTIFICATE OF ANALYSIS

MES MONOHYDRATE

BIO EXCIPIENT GRADE / ME3220-K012

LOT: ME3220-031-0520

C₆H₁₃NO₄S·H₂O ▲ F.W. 213.25 g/mol. ▲ CAS# 145224-94-8

Manufacturing Date: 5/25/2020 Retest Date: 5/31/2022

Manufacturing Site: 100 Majestic Way, Bangor PA, 18013

Packaging Date: 5/28/2020 Packaging Site: 100 Majestic Way, Bangor PA, 18013

ANALYSIS	SPECIFICATION	TEST RESULT
Absorbance (1M)	260 nm 0.1000 a.u. max. 280 nm 0.1000 a.u. max.	0.0057 a.u. 0.0045 a.u.
Appearance and Color	White / Crystals	White / Crystals
Assay	99.0% min.	100.0%
Chloride	0.005% max.	<0.005%
Color (1M, Alkaline)	Colorless	Colorless
Enzymes	DNase	None Detected
	RNase	None Detected
	Protease	None Detected
Heavy Metals (as Pb)	2 ppm max.	< 2 ppm
Identification (IR)	Passes Test	Passes Test
Loss on Drying @ 130°C	7 – 10%	9%
pH (5% Soln.)	3.1 – 3.5	3.3
pH (0.5M)	2.5 – 4.5	3.1
pK _a	5.9 – 6.3	6.2
Residue on Ignition	0.05% max.	<0.02%
Solubility (5%)	Passes Test	Passes Test
Sulfate	0.005% max.	<0.005%
Trace Elements	Arsenic (As)	2 ppm max. < 2 ppm
	Copper (Cu)	2 ppm max. < 2 ppm
	Iron (Fe)	2 ppm max. < 2 ppm
	Lead (Pb)	2 ppm max. < 2 ppm

COUNTRY OF ORIGIN: U.S.A.


TEST METHOD REFERENCE: DCN: 16-001016

The information contained herein is the property of BioSpectra. The recipient is responsible for its safe-keeping, and the prevention of unauthorized appropriation, use, disclosure and copying.

INTENDED USE: Material represented by this Certificate of Analysis is suitable for use as an excipient. It is manufactured in accordance with the ICH Q7 Good Manufacturing Practice Guide. The material represented by this Certificate of Analysis is not suitable to be used as an Active Pharmaceutical Ingredient, Drug Product or Household Item.

OVI STATEMENT: Based on the manufacturing process and the controlled handling, storage and analysis of this product, this product complies with the requirements and specifications listed in the current USP method <467> Tables 1, 2, 3, or 4.

Prepared by:  Date: 5/29/20 Job Title: QA Supervisor

Reviewed by:  Date: 05/29/20 Job Title: QA Manager