

55 GAL/200L NATURAL POLY DRUM SPECIFICATIONS

BSI COMPONENT NUMBER: COMP-0091



55-Gallon Natural Poly Drum	
Test	Specifications
Appearance and Color	Natural, closed head drum with white click lock caps
Outside Height (Nominal)	34.8"
Outside Diameter, Top (Nominal)	23.3"
Weight (Nominal)	9.52 kg. (including caps)
Capacity	55.0 U.S. Gallons (USG)
UN Rating	1H1/Y1.9/150/YR

For further information, please contact marketing@biospectra.us.



Connie Augustine | Director
Business Integration & Retention
BioSpectra Inc.

This is not considered a controlled document. We are not responsible for any errors or omissions, and the user is responsible for any decisions based on the information herein.