

12/19/2025

Version 1.1

HEAVY-DUTY HEAT TREATED PALLET SPECIFICATIONS

BSI COMPONENT NUMBER: COMP-0128



48" x 42" Heavy-Duty Wooden Pallets, Heat Treated	
Test	Specifications
Appearance and Color	Light Brown
Planks	Brown wood planks, 4 stringers
Length (Nominal)	48.0"
Width (Nominal)	42.0"
Height (Nominal)	4.75"
Stamp	Heat Treated

For further information, please contact marketing@biospectra.us.



Connie Augustine | Director
Business Integration & Retention
BioSpectra Inc.

This is not considered a controlled document. We are not responsible for any errors or omissions, and the user is responsible for any decisions based on the information herein.