



## 4x1 GALLON INDUSTRIAL SHIPPER EU COMPLIANCE

BioSpectra has evaluated the packaging component supply through the Supplier Approval Program. The 1 Gallon Non-Fluorinated HDPE Natural Bottle utilized by BioSpectra for material packaging in the 4x1 Shipping configuration meets the requirements of the Food and Drug Administration regulation 21 CFR 177.1520, which allows use of this olefin polymer in articles or components of articles intended for use in contact with food. When packaged per the Supplier's recommendations as a pre-assembled shipping kit, these industrial shippers are UN Rated for Packing Groups II and III with UN rating as indicated in the table below.

The 4x1 Gallon Industrial Shippers additionally meet the following requirements:

- ASTM D5118
- ISO 535
- 49 CFR Parts 100-185

Product Number	Size	UN Rating
COMP-0005	4x1 Gallon	4G/Y26/S

For further information, please contact [info@biospectra.us](mailto:info@biospectra.us)

**Cassie Baun**  
Senior Compliance Specialist