

# BIOSPECTRA

100 Majestic Way, Bangor, PA 18013 / [www.biospectra.us](http://www.biospectra.us)

|                      |                                |                 |                       |
|----------------------|--------------------------------|-----------------|-----------------------|
| Effective Date:      | 01-Mar-2021                    | 01-Mar-2024     | : Date of Next Review |
| Prepared By:         | Jared L Lobb                   | 16-001173 v.4.1 | : Supersedes          |
| QA/QC Approval:      | Jaron Hughes                   | Wendy Santay    | : Management Approval |
| Reason for Revision: | See Revision History in ensur. |                 |                       |

**CERTIFICATE OF ANALYSIS**  
**MES MONOHYDRATE**  
**BIO EXCIPIENT GRADE / NEW CODE MESM-3220-25**  
**(HISTORICAL CODE ME3220-K025)**  
**LOT: MESM-0121-00136**

$C_6H_{13}NO_4S \cdot H_2O$  ^ F.W. 213.3 g/mol. ^ CAS# 145224-94-8

Manufacturing Date: 7/7/21      Retest Date: 7/31/23

Manufacturing Site: 100 Majestic Way, Bangor PA, 18013

Packaging Date: 7/8/21 Packaging Site: 100 Majestic Way, Bangor PA, 18013

| ANALYSIS               | SPECIFICATION    | TEST RESULT                          |                            |
|------------------------|------------------|--------------------------------------|----------------------------|
| Absorbance (1M)        | 260 nm<br>280 nm | 0.1000 a.u. max.<br>0.1000 a.u. max. | 0.0024 a.u.<br>0.0015 a.u. |
| Appearance and Color   | White / Crystals | White / Crystals                     |                            |
| Assay                  | 99.0% min.       | 99.8%                                |                            |
| Chloride               | 0.005% max.      | <0.005%                              |                            |
| Color (1M, Alkaline)   | Colorless        | Colorless                            |                            |
| Enzymes                | DNase            | None Detected                        | None Detected              |
|                        | RNase            | None Detected                        | None Detected              |
|                        | Protease         | None Detected                        | None Detected              |
| Heavy Metals (as Pb)   | 2 ppm max.       | < 2 ppm                              |                            |
| Identification (IR)    | Passes Test      | Passes Test                          |                            |
| Loss on Drying @ 130°C | 7 – 10%          | 9%                                   |                            |
| pH (5% Soln.)          | 3.1 – 3.5        | 3.4                                  |                            |
| pH (0.5M)              | 2.5 – 4.5        | 3.2                                  |                            |
| pK <sub>a</sub>        | 5.9 – 6.3        | 6.1                                  |                            |
| Residue on Ignition    | 0.05% max.       | <0.01%                               |                            |
| Solubility (5%)        | Passes Test      | Passes Test                          |                            |
| Sulfate                | 0.005% max.      | <0.005%                              |                            |
| Trace Elements         | Arsenic (As)     | 2 ppm max.                           | < 2 ppm                    |
|                        | Copper (Cu)      | 2 ppm max.                           | < 2 ppm                    |
|                        | Iron (Fe)        | 2 ppm max.                           | < 2 ppm                    |
|                        | Lead (Pb)        | 2 ppm max.                           | < 2 ppm                    |

**COUNTRY OF ORIGIN:** U.S.A.

*The information contained herein is the property of BioSpectra. The recipient is responsible for its safe-keeping, and the prevention of unauthorized appropriation, use, disclosure and copying.*

TEST METHOD REFERENCE: DCN: 16-001016

INTENDED USE: Material represented by this Certificate of Analysis is suitable for use as an excipient. It is manufactured in accordance with the ICH Q7 Good Manufacturing Practice Guide. The material represented by this Certificate of Analysis is not suitable to be used as an Active Pharmaceutical Ingredient, Drug Product or Household Item.

OVI STATEMENT: Based on the manufacturing process and the controlled handling, storage and analysis of this product, this product complies with the requirements and specifications listed in the current USP method <467> Tables 1, 2, 3, or 4.

Prepared by: *Janet Hughes* Date: 7/27/21 Job Title: QA Specialist

Reviewed by: *C. [Signature]* Date: 7/27/21 Job Title: QA Manager