

9/17/2026

Version 1.3

30" x 48" POLY LINER SPECIFICATIONS



30" x 48" Poly Liner	
Test	Specifications
Appearance and Color	Clear with no punctures
Dimensions (Nominal)	30" x 48"
Thickness	0.0027" minimum
Weight	130 grams
Materials of Construction (Film)	Hexene LLDPE Nonslip
FDA Statement	21 CFR 177.1520
UATR	Passes Test; when compared to a previously approved poly liner for this part number

For further information, please contact marketing@biospectra.us.

A handwritten signature in black ink, appearing to read "Connie Augustine".

Connie Augustine | Director of Commercial Operations
BioSpectra Inc.

This is not considered a controlled document. We are not responsible for any errors or omissions, and the user is responsible for any decisions based on the information herein.