

4x1 GALLON (3.78L) BOTTLE SPECIFICATIONS

BSI COMPONENT NUMBER: COMP-0005



4x1 Gallon Industrial Shippers	
Test	Specifications
Appearance and Color	4x1 Gallon Non-Fluorinated HDPE Natural Round Bottle (Includes carton, corrugated pads, poly bags, locking ties and divider)
Lid	Polypropylene Cap with Polyethylene Foam Liner
Capacity	4 Bottles x 1 Gallon Each
Outside Height	11.64 inches \pm 0.5 inches
Outside Width	6.10 inches \pm 0.5 inches
UN Rating	4G/Y26/S

For further information, please contact marketing@biospectra.us.

**Connie Augustine | Director
Business Integration & Retention
BioSpectra Inc.**

This is not considered a controlled document. We are not responsible for any errors or omissions, and the user is responsible for any decisions based on the information herein.