

RESIDUAL SOLVENTS STATEMENT

L-Cystine DiHydrochloride

L-Cystine DiHydrochloride, Bio Qualified Grade is manufactured using the solvent Isopropyl Alcohol, and is stated to comply with the ICH Q3C required concentration limit of 5000ppm for a Class 3 residual solvent. BioSpectra has evaluated the material supply through the Supplier Qualification Program.

Current Product Number
LCYS-6250

For further information, please contact info@biospectra.us



Cassie Baun

Senior Compliance Specialist